



Laboratory Furniture

Co/Struc® System and Casework

General Services Administration, Federal Supply Service, Authorized Federal Supply Schedule Price List

Online access to contract ordering information, terms and conditions, up-to-date pricing, and the option to create an electronic delivery order is available through GSA Advantage!®, a menu-driven database system. The Internet address for GSA Advantage! is www.gsaadvantage.gov. For more information on ordering from Federal Supply Schedules, click on the FSS Schedules button at www.fss.gsa.gov.

FSC Group 66, Part II, Section N
FSC Classes 6640
Contract Number: GS-25F-6069D
Contract Period: February 1, 1996, through January 31, 2005

Herman Miller, Inc.
855 East Main Avenue
Zeeland, Michigan 49464
(616) 654 3000 Phone
(616) 654 8278 Fax
www.HermanMiller.com/government

Business Size: o

Price Book

Prices effective May 1, 2004
Modification No. 8
Published June 2004

Contract Information

1a. Special Item Numbers (SIN).

566-1: Modular Laboratory Furniture Systems, all product except non-modular Laboratory Tables, see 566-2 below.

566-2: Individual non-modular Laboratory Tables (C4115. and C4215.).

566-5: Workstation Design/Installation option. Including:

Workstation Design/Layout Option: Includes the design/redesign of lab furniture. Design service will be provided by Herman Miller, Inc., or its participating dealer based on a fee of \$65 per hour.

Workstation Furniture and Equipment Installation: Includes delivery coordination of furniture and related equipment, unpacking, inspection, inventory and assembly of furniture components. Pricing shall be negotiated on a project-by-project basis.

Workstation Disassembly and Reinstallation Option: Includes the taking down of installed furniture and related hardware for the purpose of packing and moving or storage. Pricing shall be negotiated on a project-by-project basis.

Workstation Repair: Includes the repair and/or replacement of damaged furniture and related hardware. Pricing shall be negotiated on a project-by-project basis.

1b. Lowest Price Model:

566-1 = C0215.27 \$14.00

566-2 = C4115.24 \$1145.00

1c. Not applicable.

2. Maximum Order. \$200,000 net product value per SIN.

3. Minimum Order. \$50 net.

4. Geographic Coverage. 48 contiguous states and District of Columbia. Other geographic locations see box below.

5. Production Points. Zeeland, Michigan, County of Ottawa.

6. Pricing. Prices shown are list. To calculate contract prices, \$0-\$25,000 list, deduct 55.2%, \$25,001 to MO, deduct 58.2%.

7. Quantity Discounts. Reference item #6.

8. Payment Terms. Net 30 days.

9a. Government Purchase Cards. ARE accepted below the micropurchase threshold.

9b. Government Purchase Cards. ARE accepted above the micropurchase threshold.

10. Foreign Items. Not applicable.

11a. Time of Delivery. 56 days ARO.

11b. Expedited Delivery. Select items may be available for expedited delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.

11c. Overnight and 2-Day Delivery. Select items may be available for overnight or 2-day delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.

11d. Urgent Requirements. Select items may be available for urgent delivery. Contact Herman Miller, Inc., or local authorized dealer for availability.

12. F.O.B. Point. Destination within the 48 contiguous states and Dist. of Columbia.

13. Order Address. Herman Miller, Inc., Government Customer Care 0161, 855 East Main Avenue, Zeeland, MI 49464.

14. Payment Address. Herman Miller, Inc., 22764 Network Place, Chicago, IL 60673-1227.

15. Warranty Provisions. Herman Miller, Inc., commercial warranty applies.

16. Export Packing Charges. Prices supplied on request.

17. Terms and Conditions of Government Commercial Credit Card Acceptance. None.

18. Terms and Conditions of Rental, Maintenance, and Repair. Not applicable.

19. Terms and Conditions of Installation. Installation services shall be negotiated between ordering agency and Herman Miller local authorized dealer.

20. Terms and Conditions of Repair Parts. Not applicable.

20a. Terms and Conditions for Any Other Services. Services shall be negotiated between ordering agency and Herman Miller local authorized dealer.

21. List of Service and Distribution Points. Access Herman Miller web site for current list of participating dealers: www.HermanMiller.com.

22. List of Participating Dealers. Access Herman Miller web site for current list of participating dealers: www.HermanMiller.com.

23. Preventive Maintenance. Not applicable.

24. Special Attributes. Reconfigurable, renewable.

25. DUNS Number. 928858430

26. Registered in Central Contractor Registration (CCR) Database. "Herman Miller, Inc."

Herman Miller Miscellaneous Information:

Federal ID Number: 38-0837640

Commercial and Government Entity Code (CAGE): OBAV4

Government facilities located in geographic areas other than those identified in item 4 above, "Geographic Coverage", write or call:
Herman Miller, Inc.


855 East Main Avenue Zeeland, MI 49464
(616) 654-3000


For answers to any contract questions, call toll free

1-800-221-7158

Introduction	4
Co/Struc® System	5
Structure	7
Work Surfaces	98
Storage	120
Transport	236
Casework	241
Storage	243
Work Surfaces	289
Appendices	295
Packaging Information	295
Keyed-Alike Information	295
Laminate Information	
Laminate Open-Line Program	297
Customer's Own Laminate Information	297
Casework Laminates	299
20-Day or Less and Assigned Lead-Time Order Information	307
10-Day or Less Order Information	309
10-Day or Less Product List	311
Indices	313
By Name	313
By Number	317

20-Day or Less/Assigned Lead-Time Shipments

All products and options in this price book not designated by an Assigned Program icon  will ship in **20 business days or less** after being acknowledged by Herman Miller.

Products and options designated by an  are on the assigned lead-time program and may ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

Introduction

A Note on the Organization of This Book

There are a few signals that will help you understand the organization of this book. Once you know them, you should be able to find your way around easily.

Like a newspaper, this book is formatted with columns of text and illustrations that run vertically.

Information wraps from one column to the next and continues for as many pages as it needs to.

Black bars are clues.

A black bar at the top of a page signals the start of information about a product.

This information is divided into two sections, each signaled by a black bar. Product Information includes a description, additional notes, and dimensional drawings. Specification Information, which always begins at the top of a column, includes a series of steps. Each step represents a decision you need to make in order to complete your product number. Steps must be followed sequentially. The number of pages needed to complete information on a product varies; some will be complete in one page; others may fill three or four pages. Just continue going through the steps until you run into the next product, signaled by the black bar at the top of the page.

Black tabs running along the outside edge of the page help you locate what section you are in or what section you want to go to. You can thumb forwards or backwards to find a section; tabs and text appear on both sides of the page.

The additional pages in the back of the book include the appendices and indices. They are signaled by a long black bar running along the outside edge of the page.

The index is presented two ways.

The first index is in alphabetical order by product name. It can be helpful if you are new to the products and numbering scheme or if you are looking for a type of product, such as work chair or round table. The second index is by product number. It lists, in alphanumeric order, the base product number, typically a five-digit number.

General Information

This book is effective May 1, 2004, subject to change without notice. Products may be purchased from authorized Herman Miller dealers who will quote prices upon request. For more information about our products and services or to see a list of dealers, visit www.HermanMiller.com.

All prices are list prices. Additional services, such as planning services, design, storage, and installation, are not included and must be added to these prices when the additional services are requested.

Illustrations, specifications, and prices are based on the latest information at the time of publication. The right is reserved to make changes at any time without notice in prices, fabrics, finishes, materials, specifications, and models, and also to discontinue models and finishes.

Dimensions are approximate. Product is shipped set up or knocked down. Products that include “Shipped knocked down” in the description require some assembly. If this statement does not appear in the description, the product is shipped assembled.

Service problems are normally handled by a Herman Miller for Healthcare authorized dealer. Where this cannot be accomplished on a local level, service problems should be referred to the Customer Care Representative for your region, Herman Miller, Inc., Zeeland, Michigan 49464.

Structure


Work Surfaces

Storage

Transport

Standard Rail

CO282



Product Information

Description

This assembly attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 25 stand-offs, and 5 spacers.


Notes

To attach assembly to wall, order rail hardware (CO28341, CO28441, or CO28541) separately.

To finish end, order end cap for standard or seismic rail (CO338) separately.

When seismic components are required, order seismic wall rail (CO479) instead of standard rail (CO282).

Dimensions



Specification Information

Step 1.

CO282FF \$191

Step 2. Surface Finish


HF inner tone light +\$0

LT light tone +\$0

Co/Struc® Structure

Rail Hardware

C0283
C0284
C0285




Co/Struc® Structure

Product Information
Description
This package contains enough fasteners to install 12' of standard rail or stand-off plate for roller rail. Package contains 300 fastener sets.
The 3" sheet metal screws are #14, flat-head, Phillips-drive screws. The 3" sheet metal screws with plastic anchors are #14, flat-head, Phillips-drive screws with 1 3/8"-long plastic anchors. The 4" toggle bolts with wall bushings are 1/4"-20, flat-head, Phillips-drive toggle bolts.

Specification Information
Step 1.
C028
Step 2. Hardware
341 3" sheet metal screws
441 3" sheet metal screws with plastic anchors
541 4" toggle bolts with wall bushings
Prices for Steps 1-2.
C028341 \$130
C028441 \$190
C028541 \$369

Rail for Roller Rail

CO345



Product Information


Description
This 72"-wide rail attaches to a wall to create a roller rail assembly. It allows lateral movement of lockers. Package contains 2.

Notes
Order rail components separately:

- End cap for roller rail (CO348)
- Roller for roller rail (CO34731)
- Stand-off plate for roller rail (CO346)

Order rail hardware (CO28341, CO28441, or CO28541) separately.

Dimensions




Specification Information

Step 1.		
CO345FF		\$273
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Structure

Roller for Roller Rail


C0347



Co/Struc® Structure


Product Information
Description
This wheel bolts into the rail for roller rail.
Notes
48 rollers are required for 12' rail; rollers are spaced 3" on center along rail.
Dimensions

Specification Information
Step 1.
C034731
\$15



End Cap for Roller Rail

CO348




Product Information

Description
This end cap attaches to each end of a roller rail and provides a positive stop. Package contains 2. Attachment hardware is included.
Dimensions


Specification Information

Step 1.		
CO348FF		\$139
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0


Co/Struc® Structure



TOP




SIDE



FRONT


Stand-Off Plate for Roller Rail

CO346




Co/Struc® Structure

Product Information
Description
This plate attaches to a wall at the bottom of a locker. It prevents friction when sliding the locker on a roller rail assembly. The plate includes 2 end caps.
Notes
Order rail hardware (CO28341, CO28441, or CO28541) separately.
Dimensions




Specification Information		
Step 1.		
CO346FF		\$108
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Product Information
Description
This end cap fastens to the end of a stand-off plate for roller rail. It provides a positive stop and a finished end. Package contains 2.
Dimensions


Specification Information		
Step 1.		
CO437FF		\$20
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Structure




Adapter Rail

C1610.




Co/Struc® Structure

Product Information
Description
This rail attaches to a panel, a module, or wall strips to support Co/Struc hanging components or a storage case in a straight-line configuration. It does not hang directly on a wall. 2 rails are required to hang a locker from an open panel frame or glazed panel.
Dimensions



Specification Information
Step 1.
C1610.
Step 2. Width
24 24" wide
30 30" wide
36 36" wide
42 42" wide
48 48" wide
60 60" wide
Prices for Steps 1-2.
C1610. 24 \$99
30 \$110
36 \$121
42 \$132
48 \$143
60 \$162
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0



Product Information

Description

This 48"-wide, aluminum rail attaches to an architectural wall. It holds work tools or accessories and can be field cut to the appropriate width. The rail includes a soft white laminate strip and 2 pair of soft white, plastic end caps to accommodate 2 field-cut rail sections.

Notes

Customer must supply attachment hardware for wall conditions.

Order rail-mounted accessories separately:


- Chart holder (CT730.)
- Chest tube rack (CT737.)
- Glove box holder (CT726.)
- Sharps container holder (CT734.)
- Supply basket (CT720. or CT722.)
- Suture basket (CT723.)
- Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)
- Utility hook (CT776.)
- Wastebasket and holder (CT728.)

Order work tools separately:

- Arch (G7150.)
- Diagonal tray (G7130.)
- Form tray (G7111.)
- Message holder (G7411.)
- Mini-shelf (G7310.)
- Organizer tray (G7210.)
- Paper tray (G7110.)
- Rail divider (G7320.)
- Tape dispenser (G7414.)
- Vertical tray (G7120.)

To mount customer-supplied accessories, order rail clamp (CT729.) separately.

Dimensions



Equipment Rail, Wall Mounted


continued

Co/Struc® Structure

Specification Information	
Step 1.	
C1170.48	\$120

Equipment Rail, Systems Mounted

C1172.



Product Information

Description

This 24"-wide aluminum rail attaches to a panel, freestanding module, or wall strips. It holds work tools or accessories. The rail includes a soft white laminate strip and 1 pair of soft white, steel end caps. Attachment hardware is included.

Notes

Order rail-mounted accessories separately:


- Chart holder (CT730.)
- Chest tube rack (CT737.)
- Glove box holder (CT726.)
- Sharps container holder (CT734.)
- Supply basket (CT720. or CT722.)
- Suture basket (CT723.)
- Tilt-out storage bin (CT9303., CT9304., CT9305., CT9306.)
- Utility hook (CT776.)
- Wastebasket and holder (CT728.)

Order work tools separately:

- Arch (G7150.)
- Diagonal tray (G7130.)
- Form tray (G7111.)
- Message holder (G7411.)
- Mini-shelf (G7310.)
- Organizer tray (G7210.)
- Paper tray (G7110.)
- Rail divider (G7320.)
- Tape dispenser (G7414.)
- Vertical tray (G7120.)

To mount customer-supplied accessories, order rail clamp (CT729.) separately.

Dimensions




Co/Struc® Structure

Equipment Rail, Systems
Mounted *continued*

Co/Struc® Structure	Specification Information	
	Step 1.	
	C1172.24	\$70


End Cap for Standard or Seismic RailCO338Rail



Product Information
Description
This end cap provides a finished end for a standard or seismic wall rail and a positive stop for a locker.
Notes
For areas requiring seismic components, order 1 end cap for each end of rail, and order seismic wall rail (CO479) separately.
For areas not requiring seismic components, order 1 end cap for each end of rail, and order standard rail (CO282) separately.
Dimensions


Specification Information		
Step 1.		
CO338FF		\$28
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Structure



Wall Rail, Seismic


CO479



Co/Struc® Structure

Product Information
Description
This aluminum rail attaches to a wall and supports Co/Struc hanging components. It includes a 12' rail, a 12' cover strip, a 12' back plate, 26 stand-offs, and 5 spacers.
Notes
Order rail for areas requiring seismic components, when attaching to studs 16" on center, or when attaching to plaster in accordance with R-0117 in California.
Order rail hardware (CO28341, CO28441, or CO28541) separately.
Order additional seismic components separately:
<ul style="list-style-type: none">• End cap for standard or seismic rail (CO338)• Toggle wall anchor for wall rail (CO483)
Field drilling is required for proper installation.
Dimensions

Specification Information		
Step 1.		
CO479FF		\$228
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0




Wall Strip

C0558

C0559

C0560



Product Information

Description

This wall strip has an aluminum channel and a steel insert. It attaches with wall fasteners to a structural wall to support hanging components. 2 wall strips are required to hang a component.

The 60"-high wall strip requires 10 fasteners; the 72"-high wall strip requires 12 fasteners; and the 84"-high wall strip requires 14 fasteners.

Notes

Order wall fasteners (AO291. or CO482) separately.

Dimensions

Specification Information

Step 1.

C05

Step 2. Height

58 60" high

59 72" high

60 84" high

Prices for Steps 1-2.


C0558 \$47

C0559 \$53

C0560 \$54




TOP



FRONT

Toggle Wall Anchor


C0482
C0483



Co/Struc® Structure


Product Information
Description
This anchor provides extra support for a wall hanger strip or wall rail. Package contains 100.
Notes
Order anchor for areas requiring seismic components.
Dimensions

Specification Information
Step 1.
C048
Step 2. Usage
3 for wall rail
2 for wall strips
Prices for Steps 1-2.
C0483
C0482



Draw Rod

CO215.



Product Information

Description

This draw rod connects Action Office® Series 1 or 2 panels or Co/Struc® C-series panels or lab modules of equal heights in a straight line. It has a threaded hole in the draw block for attachment to a seismic floor anchor bracket.


Notes

For Action Office Series 1 panels or Co/Struc C-series panels or lab modules, specify height of draw rod to match height of panel.

For Action Office Series 2 panels, specify draw rod as follows:

- 32"-high panel, CO215.27
- 39"-high panel, CO215.34
- 47"-high panel, CO215.42
- 53"-high panel, CO215.48
- 62"-high panel, CO215.57
- 67"-high panel, CO215.62
- 85"-high panel, CO215.80

Dimensions



Specification Information

Step 1.

CO215.

Step 2. Height


27	27" high
34	34" high
42	42" high
48	48" high
57	57" high
62	62" high
80	80" high

Prices for Steps 1-2.

CO215. 27	\$14
34	\$15
42	\$15
48	\$16
57	\$16
62	\$16
80	\$17

Co/Struc® Structure


Floor Anchor Bracket, Co/Struc® C0568
Lab Module



Co/Struc® Structure


Product Information
Description
This bracket fastens both sides of a module frame to the floor. Each anchor location requires 1 bracket. Finish is light tone. Package contains 10. Attachment hardware is included.
Notes
Order bracket for areas requiring seismic components.
Dimensions

Specification Information
Step 1.
C0568FF
\$542



Floor Anchor Bracket, Action
Office® Series 1

CO481



Product Information

Description

This bracket fastens both sides of an Action Office® Series 1 panel to the floor. Each anchor location requires 2 brackets. Finish is light tone. Package contains 20.

Notes

Order bracket for areas requiring seismic components.
Customer must supply required bolts and concrete anchors.

Dimensions


Specification Information

Step 1.

CO481


\$69

Co/Struc® Structure



Floor Anchor Bracket, Action Office® Series 2


CO485



Co/Struc® Structure


Product Information
Description This bracket fastens both sides of an Action Office® Series 2 panel to the floor. Package contains 10.
Notes Order bracket for areas requiring seismic components. Customer must supply required 5/8" concrete anchors.
Dimensions

Specification Information
Step 1. CO485
\$151



Floor Anchor, Ethospace® Frame

CO486




Product Information
Description
This anchor fastens an Ethospace® frame to the floor. Each frame requires 2 anchors. Package contains 10 bolts and concrete anchors.
Notes
Order anchor for areas requiring seismic components.
Order floor anchor setting tool (C1277.) separately for field installation.
Dimensions


Specification Information
Step 1.
CO486

\$179

Co/Struc® Structure




Floor Anchor Adapter KitC0552



Co/Struc® Structure


Product Information
Description
This kit replaces the bottom block of an Action Office® Series 1 or 2 spacer or 2-way, 3-way, or 4-way connector. It has a threaded hole in the draw block for attachment to a seismic floor anchor bracket. Package contains 10.
Notes
Order kit for areas requiring seismic components.
Dimensions

Specification Information
Step 1.
C0552
\$80




Panel Stiffener, Action Office®
Series 1 or Series 2 Top Cap

C0553



Product Information
Description
This top cap fits over a standard top cap and spans multiple Action Office® Series 1 or 2 panels. It provides extra rigidity along the panel run. The 12'-long top cap can be field cut to the appropriate length. Attachment hardware is included.
Notes
Order top cap for areas requiring seismic components.
Dimensions




Specification Information		
Step 1.		
C0553FF		\$428
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
MT	medium tone	+\$0

Co/Struc® Structure

Frame Stiffener, Ethospace® Top Cap


C0554



Co/Struc® Structure


Product Information
Description
This top cap replaces a standard top cap and spans multiple Ethospace® frames. It provides extra rigidity along the frame run. The 12'-long top cap can be field cut to the appropriate length. Attachment hardware is included.
Notes
Order top cap for areas requiring seismic components.
Dimensions

Specification Information
Step 1.
C0554FF
\$601
Step 2. Surface Finish
BU black umber
HF inner tone light
LU soft white



Open Support Panel


C1150.



Product Information

Description
This open panel has hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Dimensions




Specification Information

Step 1.				
C1150.				
Step 2. Height				
48	48" high			
80	80" high			
Step 3. Width				
24	24" wide			
30	30" wide			
48	48" wide			
Prices for Steps 1-3.				
		24	30	48
C1150. 48		\$325	349	423
80		\$457	494	584
Step 4. Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0

Co/Struc® Structure


Support Panel

C1151.



Co/Struc® Structure


Product Information
Description
This panel has a 48"-high recessed filler and hanger frames at 1" intervals to support components. It can be used as a return panel to stabilize a panel run. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.
Dimensions



Specification Information			
Step 1.			
C1151.			
Step 2. Height			
48	48" high		
80	80" high		
Step 3. Width			
24	24" wide		
30	30" wide		
48	48" wide		
Prices for Steps 1-3.			
	24	30	48
C1151. 48	\$437	522	668
80	\$577	668	832
Step 4. Surface Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0

Floor Entry Support Panel

C1155.



Product Information

Description

This panel has a vertical raceway for routing and concealing electrical wires and communication cables between the floor and a wire chase rail. It has $2\frac{7}{8}$ "-diameter knockouts at the floor level, a 48"-high recessed filler, and hanger frames at 1" intervals to support components. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

Order connector separately:


- 2-way 90° corner connector (C1221.)
- 3-way 90° corner connector (C1231.)

To carry 2 different types of electrical power (e.g. regular and emergency power) or to carry electrical power and low-voltage cables, vertical raceway can be divided into 2 separate internal raceways; contact local contractor.

To supply power from both sides of panel, order wire chase rail (C1310.) separately and mount 48" above floor.

For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.

Dimensions




Co/Struc® Structure

Co/Struc® Structure

Specification Information			
Step 1.			
C1155.			
Step 2. Height			
48	48" high		
80	80" high		
Step 3. Width			
24	24" wide		
30	30" wide		
48	48" wide		
Prices for Steps 1-3.			
		24	30 48
C1155. 48		\$665	761 925
80		\$822	908 1088
Step 4. Surface Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0

Ceiling Entry Support Panel

C1156.



Product Information

Description

This 80"-high panel has a vertical raceway for routing and concealing electrical wires and communication cables between the ceiling and a wire chase rail. It has 2⁷/₈"-diameter knockouts at the floor level, a 48"-high recessed filler, and hanger frames at 1" intervals to support components. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes. It cannot accept the 3- or 4-circuit electrical system and cannot be used with an Action Office® Series 2 corner connector.

Notes

Order connector separately:


- 2-way 90° corner connector (C1221.)
- 3-way 90° corner connector (C1231.)

To carry 2 different types of electrical power (e.g. regular and emergency power) or to carry electrical power and low-voltage cables, vertical raceway can be divided into 2 separate internal raceways; contact local contractor.

To supply power from both sides of panel, order wire chase rail (C1310., C1311., or C1312.) separately and mount 48" above floor.

For power feed to wire chase rail mounted at heights other than 48", field cut removable cover on 1 side of vertical raceway.

Dimensions




Ceiling Entry Support Panel *continued*

Co/Struc® Structure

Specification Information		
Step 1.		
C1156.80		
Step 2. Width		
24	24" wide	
30	30" wide	
48	48" wide	
Prices for Steps 1-2.		
C1156.80	24	\$895
	30	\$1000
	48	\$1182
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Terminal Panel

C1161.



Product Information

Description

The center point of this 60"-wide panel attaches at a 90° angle to other panels and provides support to a panel run. Each panel has adjustable glides, finished ends, and a 48"-high recessed filler and hanger frames at 1" intervals to support components. The panel has an area between the filler and the rear edge of an attached work surface for passage of wires or tubes.


The panel is available in 2 heights. The 80"-high panel supports a locker; 48"- and 80"-high panels support work surfaces and other 60"-wide components.

Notes


To support locker, order 2 adapter rails (C1610.) separately.

To connect to panel at 90° angle, order draw rod (AO215.) separately.


Dimensions



TOP



48" HIGH
FRONT




80" HIGH
FRONT

Co/Struc® Structure

Specification Information		
Step 1.		
C1161.		
Step 2. Height		
4860	48" high	
8060	80" high	
Prices for Steps 1-2.		
C1161. 4860		\$853
8060		\$1110
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Panel Wall Start


C1211.




Product Information
Description
This connector attaches a support panel to a wall or column for a peninsula configuration.
Notes
Order draw rod (AO215.) separately.
Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.
Dimensions

Specification Information
Step 1.
C1211.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1211. 48 \$151
80 \$214
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Structure




2-Way 90° Corner ConnectorC1221.




Co/Struc® Structure

Product Information
Description
This connector joins 2 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.
Notes
Order 2 draw rods (AO215,) separately.
When connecting panels of unequal heights, specify height of connector to match height of higher panel.
Dimensions



TOP



SIDE



FRONT


Specification Information
Step 1.
C1221.
Step 2. Height
4848" high
8080" high
Prices for Steps 1-2.
C1221. 48
80
Step 3. Surface Finish
HFinner tone light
LTlight tone

\$120
\$148

+\$0
+\$0

3-Way 90° Corner Connector

C1231.



Product Information

Description



This connector joins 3 support panels of equal or unequal heights at a 90° angle. It cannot be used with a panel run that supports Action Office® Series 2 work surfaces.

Notes

Order 3 draw rods (AO215.) separately.

When connecting panels of unequal heights, specify height of connector to match height of higher panel.

Dimensions



SIDE



FRONT

Specification Information

Step 1.

C1231.

Step 2. Height

48 48" high

80 80" high


Prices for Steps 1-2.

C1231. 48	\$178
80	\$220

Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0



Interface ConnectorC1261.



Co/Struc® Structure


Product Information
Description
This 12"-wide connector joins a module to a support panel. It has a removable cover strip that conceals predrilled holes used for attaching connectors.
Notes
To complete shelf surface at end of module run, order shelf end cap (C1510.E or C1511.E) separately.
Dimensions

Specification Information
Step 1.
C1261.
Step 2. Height
4812 48" high
8012 80" high
Prices for Steps 1-2.
C1261. 4812 \$417
8012 \$470
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0




Panel Connector

C1281.




Product Information
Description
This connector attaches a support panel to an interface connector.
Notes
Specify connector height to match panel height.
Order draw rod (AO215.) separately to match panel height.
Dimensions



Specification Information
Step 1.
C1281.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1281. 48 \$51
80 \$67
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Structure


Mid-Run Support, Single SidedC1272.




Co/Struc® Structure

Product Information
Description
This brace provides additional support to a straight panel run that has components hanging on 1 side. It includes a leveling glide.
Notes
To support 24"-deep work surface, specify 21"-wide mid-run support; to support 30"-deep work surface, specify 27"-wide mid-run support. Order 1 draw rod (AO215.) for each support to panel connection.
Dimensions

Specification Information		
Step 1.		
C1272.		
Step 2. Height		
48	48" high	
80	80" high	
Step 3. Width		
21	21" wide	
27	27" wide	
Prices for Steps 1-3.		
	21	27
C1272. 48	\$341	376
80	\$471	508
Step 4. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Mid-Run Support, Double Sided C1273.



Product Information

Description


This brace provides additional support to a straight panel run that has components hanging on both sides. It includes leveling glides.

Notes


To support 24"-deep work surface, specify 41"-wide mid-run support; to support 30"-deep work surface, specify 53"-wide mid-run support. To support 24"-deep work surface on 1 side and 30"-deep work surface on other side, specify 53"-wide mid-run support.

Order 1 draw rod (AO215,) for each support to panel connection.


Dimensions



TOP



FRONT



Specification Information

Step 1.

C1273.

Step 2. Height

48 48" high

80 80" high

Step 3. Width

41 41" wide


53 53" wide

Prices for Steps 1-3.

	41	53
C1273. 48	\$471	547
80	\$642	717


Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0



Support Panel Top Cap


C1274.



Co/Struc® Structure


Product Information
Description
This 20"-wide aluminum top cap attaches to 2 support panels in a straight connection. It provides additional support for a panel run.
Dimensions

Specification Information		
Step 1.		
C1274.20		\$50
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Terminal Panel Bracket

C3920.
C3921.



Product Information

Description

This bracket provides support for 2 30"-wide, 3- or 5-shelf units attached to a terminal panel or module terminal panel.


Notes

Specify 24"-high bracket for 3-shelf unit (SM303.0930 or SM303.1630); specify 34"-high bracket for 5-shelf unit (SM305.0930 or SM305.1630).

Order 1 of the following support components separately:

- Low module terminal panel (C1122.4860)
- Pass-through module terminal panel (C1122.8060)
- Terminal panel (C1161.)

Dimensions



Specification Information

Step 1.

C392

Step 2. Usage

- 0. for module terminal panel
- 1. for terminal panel

Step 3. Height

- 24 24" high
- 34 34" high


Prices for Steps 1-3.

C3920.	24	\$140
	34	\$150

C3921.	24	\$140
	34	\$150

Step 4. Surface Finish


HF	inner tone light	+\$0
LT	light tone	+\$0



Co/Struc® Structure


Product Information
Description
This 48"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services. The module includes a top shelf and 4 adjustable glides. It has hanger frames at 1" intervals to support components.
Notes
Order access panels separately.
For L, T, or X configuration, order module corner connector separately:
<ul style="list-style-type: none">• 2-way connector (C1220.)• 3-way connector (C1230.)• 4-way connector (C1240.)
Dimensions

Specification Information
Step 1.
C1112.48
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-2.
C1112.48 24 \$934
30 \$1100
48 \$1237
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0



Pass-Through Module

C1114.



Product Information

Description

This 80"-high module attaches to a wall or connects to other freestanding modules. It encloses and organizes mechanical services and provides pass-through access. The module includes an interior shelf that is field adjustable to 30", 36", or 48" high. It has hanger frames at 1" intervals to support components and 4 adjustable glides.

Notes

Order access panels separately.

For L, T, or X configuration, order module corner connector separately:

- 2-way connector (C1220.)
- 3-way connector (C1230.)
- 4-way connector (C1240.)

Dimensions

Specification Information

Step 1.

C1114.80

Step 2. Width


24	24" wide
30	30" wide
48	48" wide


Prices for Steps 1-2.

C1114.80	24	\$1163
	30	\$1300
	48	\$1507

Step 3. Surface Finish

KH	light tone	+\$0
KR	inner tone light	+\$0





Co/Struc® Structure

Product Information

Description


This 48"-high, hard-surfaced end panel attaches to a freestanding module. It supports a module run or divides space. The panel includes a top shelf end cap, adjustable glides, and finished ends.

The 60"-wide panel supports hanging components.

Notes


To attach terminal panel to module, order 2 draw rods (AO215.) separately.

Dimensions




Low Module Terminal Panel *continued*

Specification Information		
Step 1.		
C1122.48		
Step 2. Width/Position		
36L	36" wide and left-hand panel	
36R	36" wide and right-hand panel	
60	60" wide	
Prices for Steps 1-2.		
C1122.48	36L	\$868
	36R	\$868
	60	\$1036
Step 3. Surface Finish		
KH	light tone	+\$0
KR	inner tone light	+\$0



Co/Struc® Structure


Product Information
Description
This 32"-high module attaches to the top of a low module and provides overhead storage capability. Hanger frame slots are at 1" increments to support hanging components. The module has finished ends that can be removed when 2 or more add-on lab modules are used adjacent to one another. Attachment hardware is included.
Notes
Specify width of add-on module to match width of low module. Order low module (C1112.) separately.
Dimensions



Specification Information
Step 1.
C1119.32
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-2.
C1119.32 24 \$400
30 \$500
48 \$600
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Pass-Through Module Terminal
Panel

C1122.
C1124.



Product Information


Description

This 80"-high, hard-surfaced end panel attaches to a freestanding module. It supports a module run or divides space. The panel is available in 2 widths; the 60"-wide, full-height panel supports hanging components. Each panel has an interior shelf end cap, adjustable glides, and finished ends.

Notes

To attach terminal panel to module, order 2 draw rods (AO215.) separately.

Dimensions




Pass-Through Module Terminal
Panel *continued*

Co/Struc® Structure

Specification Information		
Step 1.		
C112		
Step 2. Height		
2.80	full-height	
4.80	half-height	
Step 3. Width/Position		
<i>For full-height (2.80)</i>		
60	60" wide	
<i>For half-height (4.80)</i>		
36L	36" wide and left-hand panel	
36R	36" wide and right-hand panel	
60	60" wide	
Prices for Steps 1-3.		
C1122.80	60	\$1350
C1124.80	36L	\$1035
	36R	\$1035
	60	\$1155
Step 4. Surface Finish		
KH	light tone	+\$0
KR	inner tone light	+\$0

Diplomat Panel

C1146.




Product Information

Description
This 48"-high, 48"-wide laminate panel can be field cut to fit a non-modular space in a run of freestanding modules. It includes 3 top shelf brackets, 2 floor base brackets, and 2 wall-attached wall start connectors. The panel cannot support components.

Notes
To cover both sides of a freestanding module run, specify 2 diplomat panels.
For diplomat panel used with low or pass-through module, order top shelf (C1510.) separately; shelf must be field modified.
For diplomat panel used with wall-attached module, order top shelf for wall-attached module (C1514.) separately.

Dimensions




Specification Information

Step 1.		
C1146.		\$489
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Structure

Vertical Chase Cover


C1460.




Co/Struc® Structure

Product Information
Description
This 72"-high tower encloses services between the ceiling and the top shelf of a low module, the interior shelf of a pass-through module, or a 2- or 3-way module corner connector. It includes a ceiling trim plate. The cover cannot be used with a 4-way module corner connector.
Notes
Order module or connector separately:
<ul style="list-style-type: none">• 2-way module corner connector (C1220.)• 3-way module corner connector (C1230.)• Low module (C1112.)• Pass-through module (C1114.)
Dimensions


Specification Information
Step 1.
C1460. \$445
Step 2. Surface Finish
HF inner tone light +\$0
LT light tone +\$0



TOP




SIDE



FRONT

Module Wall Start


C1212.




Product Information
Description
This connector starts a peninsula configuration by attaching a low or pass-through module to a wall or column. It includes a shelf end cap.
Notes
For low module, specify 48"-high wall start; for pass-through module, specify 80"-high wall start.
To connect wall start to module, order 2 draw rods (AO215.) separately.
Package does not contain attachment hardware; designer must select fasteners appropriate for wall conditions and installer must field install.
Dimensions

Specification Information
Step 1.
C1212.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1212. 48 \$305
80 \$402
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0

Co/Struc® Structure




2-Way Module Corner Connector C1220.




Co/Struc® Structure

Product Information
Description
This connector joins 2 modules of equal or unequal heights at a 90° angle. It includes a top shelf and filler panels to enclose the utility chase.
Notes
To attach module of equal height, order 2 draw rods (AO215.) separately.
To attach module of unequal height, order 2 change-of-height draw rods (AO219.) and 2 change-of-height panel/panel finished ends (AO251.) separately.
Dimensions

Specification Information
Step 1.
C1220.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1220. 48 \$1000
80 \$1200
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0




3-Way Module Corner Connector C1230.




Product Information
Description
This connector joins 3 modules of equal or unequal heights at 90° angles. It includes a top shelf and filler panels to enclose the utility chase.
Notes
To attach module of equal height, order 2 draw rods (AO215.) separately.
To attach module of unequal height, order 2 change-of-height draw rods (AO219.) and 2 change-of-height panel/panel finished ends (AO251.) separately.
Dimensions

Specification Information
Step 1.
C1230.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1230. 48 \$1000
80 \$1200
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0

Co/Struc® Structure




4-Way Module Corner Connector C1240.



Co/Struc® Structure


Product Information
Description
This connector joins 4 modules of equal or unequal heights at 90° angles. It includes a top shelf.
Notes
To attach module of equal height, order 2 draw rods (AO215.) separately.
To attach module of unequal height, order 2 change-of-height draw rods (AO219.) and 2 change-of-height panel/panel finished ends (AO251.) separately.
Dimensions

Specification Information
Step 1.
C1240.
Step 2. Height
48 48" high
80 80" high
Prices for Steps 1-2.
C1240. 48 \$1000
80 \$1200
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0



Module Support Leg

C1275.




Product Information


Description
This leg provides additional support for normal loading on a low or pass-through module. It has a leveling glide.
Dimensions

Specification Information

Step 1.		
C1275.		\$167
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0


Co/Struc® Structure





Anchor Kit


C1276.




Co/Struc® Structure

Product Information
Description
This kit fastens a module to the floor. It includes a floor anchor, bolt, spacer, and plate. Package contains 10.
Dimensions


Specification Information
Step 1.
C1276. \$155





Floor Anchor Setting Tool

C1277.




Product Information

Description
This steel tool is used to fasten anchors in a flush position without damaging the anchor or module.

Specification Information


Step 1.
C1277. \$28

Co/Struc® Structure



Interior Lab Bracket


C1191.



Co/Struc® Structure


Product Information
Description
This bracket supports gas, water, or vacuum lines and electrical conduit in the service chase of a module. Package contains 12.
Dimensions

Specification Information
Step 1.
C1191. \$191




Terminal Panel Bracket

C3920.
C3921.




Product Information
Description This bracket provides support for 2 30"-wide, 3- or 5-shelf units attached to a terminal panel or module terminal panel.
Notes Specify 24"-high bracket for 3-shelf unit (SM303.0930 or SM303.1630); specify 34"-high bracket for 5-shelf unit (SM305.0930 or SM305.1630). Order 1 of the following support components separately: <ul style="list-style-type: none"> • Low module terminal panel (C1122.4860) • Pass-through module terminal panel (C1122.8060) • Terminal panel (C1161.)
Dimensions



Specification Information
Step 1. C392
Step 2. Usage 0. for module terminal panel 1. for terminal panel
Step 3. Height 24 24" high 34 34" high
Prices for Steps 1-3.
C3920. 24 \$140 34 \$150
C3921. 24 \$140 34 \$150
Step 4. Surface Finish
HF inner tone light +\$0 LT light tone +\$0


Top Shelf

C1510.



Co/Struc® Structure


Product Information
Description
This 11"-deep shelf encloses the top of a low module service chase and provides a storage shelf. It has a Chemsurf® top and edge. The shelf can be drilled for fixtures or piping. It cannot attach to a pass-through module.
Dimensions




Specification Information
Step 1.
C1510.12
Step 2. Width
12 12" wide
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-2.
C1510.12 12 \$168
24 \$226
30 \$295
48 \$319
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0

Interior Shelf

C1511.



Product Information
Description
This 11"-deep shelf fits inside a pass-through module. It encloses the top of a 48"-high service chase and provides a storage shelf. The shelf has a Chemsurf® top and edge. It can be drilled for fixtures or piping. Attachment hardware is included.
Notes
To install shelf higher than 48", order interior shelf and support kit (C1515.) separately.
Dimensions




Specification Information
Step 1.
C1511.12
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-2.
C1511.12 24 \$280
30 \$327
48 \$369
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0

Co/Struc® Structure


Shelf End Cap

C1510.
C1511.




Co/Struc® Structure

Product Information
Description
This cap finishes the end of an interior or top shelf at the connection of a low or pass-through module and an interface connector. It has a Chemsurf® top and sides. Attachment hardware is included.
Dimensions



SIDE




TOP

Specification Information
Step 1.
C151
Step 2. Depth/Usage
1.E 7" deep for interior shelf
0.E 12" deep for top shelf
Prices for Steps 1-2.
C1511.E \$57
C1510.E \$60
Step 3. Surface Finish
KH light tone +\$0
KR inner tone light +\$0

Interior Shelf and Support Kit


C1515.



Product Information

Description
This 11"-deep interior shelf is installed within a pass-through module. Attachment hardware is included.

Dimensions




Specification Information

Step 1.		
C1515.12		
Step 2. Width		
24	24" wide	
30	30" wide	
48	48" wide	
Prices for Steps 1-2.		
C1515.12	24	\$450
	30	\$500
	48	\$556
Step 3. Surface Finish		
KH	light tone	+\$0
KR	inner tone light	+\$0

Co/Struc® Structure


Shelf Backsplash

C1520.



Co/Struc® Structure


Product Information
Description
This plastic barrier attaches to the back of a top shelf and creates a backsplash. It is used with a module that sits against a wall. The barrier can also compensate for variations in a wall surface. Package contains 4.
Notes
For use with interior shelf, field modify shelf backsplash.
Dimensions



Specification Information
Step 1.
C1520.
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-2.
C1520. 24 \$58
30 \$75
48 \$86
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Wall-Attached Module


C1141.



Product Information
Description This 38"-high module attaches parallel to a permanent wall. It encloses and organizes mechanical services and supports hanging components. The module includes a top shelf, 2 floor base brackets, and adjustable glides.
Notes Order access panels separately. Order backsplash and attachment hardware separately: <ul style="list-style-type: none"> • Shelf backsplash (C1520.) for each module • Wall-attachment channel for connecting wall-attached module to wall; must be customer supplied and field installed • Wall-attachment channel hardware (C1292.) for every 400' of module run • Wall-attached connector (C1293.) for each module connection
Dimensions


Specification Information
Step 1. C1141.38
Step 2. Width 24 24" wide 30 30" wide 48 48" wide
Prices for Steps 1-2. C1141.38 24 \$655 30 \$750 48 \$887
Step 3. Surface Finish KH light tone +\$0 KR inner tone light +\$0

Co/Struc® Structure



Wall Attachment Channel
Hardware

C1292.




Co/Struc® Structure

Product Information
Description
This package has fasteners for attaching 400' of wall attachment channel to a wall. Package includes 300 wall bushings, washers, screws, and toggle anchors.

Specification Information
Step 1.
C1292.
\$377

Wall-Attached Connector


C1293.



Product Information
Description This connector fastens a wall-attached module to a wall attachment channel or end closure for wall-attached module. It includes a floor base bracket. The connector fastened to a wall-attached module creates a service chase with an inside clearance of 7 ¹ / ₄ ".
Notes When attaching to wall-attached module, specify 1 connector for each module. When enclosing end of wall-attached run with end closure, specify 1 connector.
Dimensions


Specification Information
Step 1. C1293.
\$89

Co/Struc® Structure



Wall-Attached Wall Start
Connector


C1295.



Co/Struc® Structure


Product Information
Description
This connector allows a run of wall-attached modules to start directly adjacent to an architectural wall corner.
Dimensions

Specification Information
Step 1.
C1295.
\$74



Top Shelf for Wall-Attached Module

C1514.




Product Information

Description

This shelf encloses the top of a wall-attached module service chase and provides a surface for reagent storage. It has a Chemsurf® top and edge. The shelf can be drilled for fixtures or piping. It cannot be installed on a low or pass-through module.

Dimensions



Specification Information

Step 1.

C1514.10

Step 2. Width

24	24" wide
30	30" wide
48	48" wide

Prices for Steps 1-2.

C1514.10	24	\$208
	30	\$255
	48	\$293


Step 3. Surface Finish

KH	light tone	+\$0
KR	inner tone light	+\$0

Co/Struc® Structure


Filler Extrusion

C1530.
C1531.




Co/Struc® Structure

Product Information
Description
This component fills the space between a work surface and the top shelf for a wall-attached module or the interior shelf for a pass-through module. It creates a continuous work area. Attachment hardware is included.
Dimensions




Specification Information
Step 1.
C153
Step 2. Usage
0. for top shelf
1. for interior shelf
Step 3. Width
24 24" wide
30 30" wide
48 48" wide
Prices for Steps 1-3.
C1530. 24 \$44
30 \$55
48 \$75
C1531. 24 \$45
30 \$65
48 \$87
Step 4. Surface Finish
HF inner tone light +\$0
LT light tone +\$0




Floor Track

C1192.




Product Information
Description
This 8' aluminum track provides the bottom closure for a vinyl base on a module or panel. It can be field cut to the appropriate length. Finish is black. Package contains 12 tracks (96' total).
Notes
Order vinyl base (C1193.) separately.
Dimensions




Specification Information
Step 1.
C1192.96
\$717

Co/Struc® Structure



Vinyl Base


C1193.




Co/Struc® Structure


Product Information
Description
This 4"-high, 12'-long vinyl base attaches to the floor track under a module or panel. It can be field cut to the appropriate length. Package contains 8 bases (96' total).
Notes
Order floor track (C1192.) separately.
Dimensions

Specification Information		
Step 1.		
C1193.144		\$241
Step 2. Surface Finish		
MT	medium tone	+\$0






Vinyl Base, 90° Outside Corner C1194.




Product Information
Description
This 4"-high vinyl base trims the 90° corner of a module or panel.
Package contains 10.
Dimensions

Specification Information		
Step 1.		
C1194.90		\$333
Step 2. Surface Finish		
MT	medium tone	+\$0


Co/Struc® Structure





Vinyl Base, End Cap


C1195.



Co/Struc® Structure


Product Information
Description
This 4"-high, 31"-long vinyl base trims a finished end, support panel, or terminal panel. It can be field cut to the appropriate length. Package contains 10.
Dimensions

Specification Information		
Step 1.		
C1195.		\$562
Step 2. Surface Finish		
MT	medium tone	+\$0



Access Panel Assembly

C1412.



Product Information

Description

This assembly of 2 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down.

The assembly has 2 30"- or 2 42"-high access panels.


Notes

Order low module (C1112.) or pass-through module (C1114.) separately.

For wall-attached module or interior shelf mounted at 36" on pass-through module, specify 30"-high assembly.

For interior shelf mounted at 48" on pass-through module, specify 42"-high assembly.

Dimensions




Co/Struc® Structure

Specification Information			
Step 1.			
C1412.			
Step 2. Height			
30	30" high		
42	42" high		
Step 3. Width			
19	19 ¹ / ₄ " wide		
25	25 ¹ / ₄ " wide		
43	43 ¹ / ₄ " wide		
Prices for Steps 1-3.			
		19	25 43
C1412. 30		\$376	439 474
42		\$423	494 596
Step 4. Surface Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0

Combination Access Panel Assembly

C1422.



Product Information

Description

This assembly of 4 access panels attaches to a freestanding module and encloses both sides of the service chase. The panels mount in a flush or recessed position and are removable for service access. Each panel has 2 knockouts (4 per module side) for standard duplex electrical outlets. Brackets and mounting hardware are included. Shipped knocked down. The 24"-high assembly has 4 12"-high access panels; the 42"-high assembly has 2 12"- and 2 30"-high access panels.


Notes

Order low module (C1112.) or pass-through module (C1114.) separately.

For interior shelf mounted 30" high on pass-through module, specify 24"-high assembly.

For interior shelf mounted 48" high on pass-through module, specify 42"-high assembly.

Dimensions




Combination Access Panel
Assembly *continued*

Co/Struc® Structure

Specification Information				
Step 1.				
C1422.				
Step 2. Height				
24	24" high			
42	42" high			
Step 3. Width				
19	19 ¹ / ₄ " wide			
25	25 ¹ / ₄ " wide			
43	43 ¹ / ₄ " wide			
Prices for Steps 1-3.				
		19	25	43
C1422. 24		\$654	686	744
42		\$592	659	738
Step 4. Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0

Low Cup Sink Access Panel Assembly

C1432.



Product Information

Description

This assembly attaches to a pass-through module and creates a low cup sink assembly in the lower service chase. It encloses both sides of the service chase and has removable panels for service access. The assembly includes 4 12"-high access panels, 1 shelf support frame, 1 interior shelf, 1 bottom chase cover, and attachment hardware. Each 12"-high panel has 2 knockouts for standard duplex electrical outlets. The assembly has a 15³/₄" space for a cup sink; this divides the service chase into 2 12" segments. It cannot attach to a wall-attached module. Shipped knocked down.


Notes

Order pass-through module (C1114.) separately.

To enclose service chase of adjacent module, order center closure panel (C1455.1607) separately.

Customer must order sink and fixtures separately and field cut interior shelf to hold sink.

Dimensions




Low Cup Sink Access Panel
Assembly *continued*

Co/Struc® Structure

Specification Information		
Step 1.		
C1432.42		
Step 2. Width		
19	19 ¹ / ₄ " wide	
25	25 ¹ / ₄ " wide	
43	43 ¹ / ₄ " wide	
Prices for Steps 1-2.		
C1432.42	19	\$1090
	25	\$1200
	43	\$1340
Step 3. Surface Finish		
KH	light tone	+\$0
KR	inner tone light	+\$0

End Closure Panel


C1450.



Product Information
Description
This 32"-high, 11"-wide panel encloses the top section of a pass-through module.
Notes
Order pass-through module (C1114.) separately.
Dimensions


Specification Information		
Step 1.		
C1450.3211		\$100
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Structure



End Closure Assembly

C1451.
C1452.



Co/Struc® Structure

Product Information

Description


This assembly encloses the end of a module when a terminal panel is not used.

The end closure for a low module has a 48"-high panel and 2 48"-high end caps. The full-height end closure for a pass-through module has an 80"-high panel and 2 80"-high end caps. The half-height end closure for a pass-through module has a 48"-high panel and 2 80"-high end caps.

Notes

Order low module (C1112.) or pass-through module (C1114.) separately.


Dimensions



END CLOSURE
FOR LOW MODULE

SIDE


FRONT



FULL-HEIGHT END
CLOSURE FOR PASS-THROUGH MODULE

SIDE

FRONT



HALF-HEIGHT END
CLOSURE FOR PASS-THROUGH MODULE

SIDE


FRONT

End Closure Assembly *continued*

Specification Information		
Step 1.		
C145		
Step 2. Type		
1.4811	end closure for low module	
1.8011	full-height end closure for pass-through module	
2.8011	half-height end closure for pass-through module	
Prices for Steps 1-2.		
C1451.4811		\$184
C1451.8011		\$222
C1452.8011		\$197
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

End Closure for Wall-Attached
Module


C1453.



Co/Struc® Structure

Product Information
Description
This panel attaches to the left or right end of a wall-attached module run and encloses the service chase.
Notes
Order wall-attached connector (C1293.) separately.
Dimensions

Specification Information		
Step 1.		
C1453.3809		\$346
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0





Center Closure Panel

C1455.



Product Information
Description
This panel attaches to the inside of a module run and encloses the end of a service chase. Attachment hardware is included.
Notes
To enclose service chase on either end of cup sink opening or 14"-high openings created by partial chases at panel connections, order 16"-high panel.
To enclose lower interior module chase, order 43"-high panel.
To enclose upper section of pass-through module, order 32"-high panel.
Dimensions




Specification Information
Step 1.
C1455.
Step 2. Height
1607 16" high
3207 32" high
4307 43" high
Prices for Steps 1-2.
C1455. 1607 \$108
3207 \$124
4307 \$145
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Structure

Wire Chase Rail

C1310.
C1311.
C1312.



Product Information

Description

This rail attaches to 1 side of an equal-width module or support panel and provides a cable management channel. It can be divided into 2 separate channels. The rail holds work tools. It cannot be used with an Action Office® Series 1 or 2 connector.

The rail with type 1 or 2 receptacle covers has the following knockouts per width:

- 24"- and 30"-wide rails: 1 knockout
- 48"- and 60"-wide rails: 2 knockouts
- 72"-wide rail: 3 knockouts
- 96"-wide rail: 4 knockouts


Notes

Order connection components separately:


- Wire chase rail bridge (C1325.) to connect across mid-run support panel or back side of 3-way connector
- 90° inside corner wire chase rail (C1320.) to connect around 90° inside corner

Receptacles must be customer supplied.


Dimensions



SIDE



TYPE 1
FRONT




TYPE 2
FRONT

Specification Information		
Step 1.		
C131		
Step 2. Cover		
0.	full-length cover	
1.	type 1 receptacle cover	
2.	type 2 receptacle cover	
Step 3. Width		
24	24" wide	
30	30" wide	
48	48" wide	
60	60" wide	
72	72" wide	
96	96" wide	
Prices for Steps 1-3.		
C1310. 24		\$138
30		\$151
48		\$214
60		\$272
72		\$372
96		\$452
C1311. 24		\$184
30		\$198
48		\$297
60		\$381
72		\$486
96		\$605
C1312. 24		\$184
30		\$198
48		\$297
60		\$381
72		\$486
96		\$605
Step 4. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Wire Chase Rail, 90° Inside Corner

C1320.



Co/Struc® Structure

Product Information


Description

This rail attaches to 2 modules or support panels at a 90° inside corner. It provides a cable management channel and can be divided into 2 separate channels. The rail holds work tools. It cannot be used with an Action Office® Series 1 or 2 connector.

Notes

When installing receptacles, order wire chase rail receptacle cover (C1340., C1341., C1342., or C1343.) separately.

Dimensions



Specification Information

Step 1.

C1320.

Step 2. Width

24	24" wide
30	30" wide
48	48" wide

Step 3. Depth

For 24" wide (24)

24	24" deep
30	30" deep
48	48" deep

For 30" wide (30)

30	30" deep
48	48" deep

For 48" wide (48)

48	48" deep
-----------	----------

Prices for Steps 1-3.


	24	30	48
C1320. 24	\$355	372	438
30	—	\$389	452
48	—	—	\$508

Step 4. Surface Finish

HF	inner tone light		+\$0
LT	light tone		+\$0

Wire Chase Rail Bridge


C1325.




Product Information
Description
This bridge joins wire chase rails across a mid-run support panel or on the back side of a 3-way 90° connector. It allows the continuation of a cable management channel. The bridge cannot be used with an Action Office® Series 1 or 2 connector.
Notes
Order 1 of the following products separately:
<ul style="list-style-type: none">• 3-way 90° connector (C1231.)• Mid-run support (C1272. or C1273.)
Dimensions

Specification Information
Step 1.
C1325.
Step 2. Surface Finish
HF inner tone light
LT light tone

Co/Struc® Structure




Wire Chase Rail End CapC1330.



Co/Struc® Structure


Product Information
Description
This end cap finishes the exposed end of a wire chase rail. Package contains 2.
Dimensions

Specification Information		
Step 1.		
C1330.		\$69
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Receptacle Cover, Wire Chase Rail

C1340.
C1341.
C1342.
C1343.




Product Information

Description


This plate covers receptacle outlets on a wire chase rail. Package contains a receptacle holder and cover.

The type 1 cover accepts 1 standard duplex outlet; the type 2 cover accepts 2 standard duplex outlets. The type 3 cover accepts 1 rectangular-faced device; the type 4 cover accepts 2 rectangular-faced devices.


Dimensions




TYPE 1
COVER



TYPE 2
COVER



TYPE 3
COVER



TYPE 4
COVER

Specification Information

Step 1.

C134

Step 2. Type

- 0. type 1
- 1. type 2
- 2. type 3
- 3. type 4

Prices for Steps 1-2.


C1340.	\$24
C1341.	\$24
C1342.	\$25
C1343.	\$25

Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Work Surface

C2310.



Product Information

Description

This laminate, rectangular work surface hangs from a panel, module, or wall strips and has a general purpose top-loading capacity of up to 200 pounds. The support arm finish matches the surface finish. Attachment hardware is included.

The 96"-wide surface includes 1 additional support bracket. The 24"-wide surface cannot accept components.


Notes

When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

To order optional storage products, see Storage.

Order optional keyboard tray (G7718.) separately.

Dimensions



Specification Information

Step 1.

C2310.

Step 2. Depth

24 24" deep

30 30" deep

Step 3. Width

24 24" wide

30 30" wide

48 48" wide

60 60" wide

72 72" wide

96 96" wide

Prices for Steps 1-3.


	24	30	48	60	72	96
C2310. 24	\$245	272	348	447	523	702
30	\$272	321	429	557	654	835

Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Heavy-Duty Work Surface

C2210.



Product Information

Description

This rectangular work surface hangs from a panel, module, or wall strips and has a top-loading capacity of up to 400 pounds. It can hold drawers and includes drawer bearers. The support arm finish matches the surface finish. The work surface cannot accept an option work surface C frame adapter kit. Attachment hardware is included. Shipped knocked down.

The work surface has a laminate or Chemsurf® top. A 51"-wide surface has a 3" overhang on the left or right side to fill the gap created by a mid-run support, C-series 3-way connector, or interface connector. A 96"-wide surface includes 1 additional support bracket. The 24"-wide surface cannot accept components.

Notes

Work surface supports following quantities of A-, B-, or C-size drawers or storage cases:

- 30" wide: 1 drawer or storage case
- 48", 51", and 60" wide: 2 drawers or storage cases
- 72" wide: 3 drawers or storage cases
- 96" wide: 4 drawers or storage cases


Order storage components separately:

- Drawer (C0207, C0208, or C0209)
- Storage case (SM400.) and adapter rail (C1610.)

When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard tray (G7718.) separately.

Dimensions



Specification Information

Step 1.

C2210. [A]

Step 2. Depth

24 24" deep [A]

30 30" deep [A]

Step 3. Width

24 24" wide [A]

30 30" wide [A]

48 48" wide [A]

51L 51" wide, left overhang [A]

51R 51" wide, right overhang [A]

60 60" wide [A]

72 72" wide [A]

96 96" wide [A]

Prices for Steps 1-3.

		24	30	48	51L	51R	60	72
C2210.	24	\$245	272	400	415	415	491	576
	30	\$272	321	495	511	511	614	719
								96
C2210.	24							\$773
	30							\$921

Step 4. Surface Finish

Laminate Top

		2424	2430	2448	2451L	2451R
HF	inner tone light [A]	+\$0	0	0	0	0
LT	light tone [A]	+\$0	0	0	0	0
		2460	2472	2496	3024	3030
HF	inner tone light [A]	+\$0	0	0	0	0
LT	light tone [A]	+\$0	0	0	0	0
		3048	3051L	3051R	3060	3072
HF	inner tone light [A]	+\$0	0	0	0	0
LT	light tone [A]	+\$0	0	0	0	0
						3096
HF	inner tone light [A]					+\$0
LT	light tone [A]					+\$0


Heavy-Duty Work Surface *continued*

Co/Struc® Work Surfaces

Chemsurf® Top						
		2424	2430	2448	2451L	2451R
KH	light tone <input type="checkbox"/> A	+\$75	84	101	106	106
KR	inner tone light <input type="checkbox"/> A	+\$75	84	101	106	106
		2460	2472	2496	3024	3030
KH	light tone <input type="checkbox"/> A	+\$125	145	195	84	98
KR	inner tone light <input type="checkbox"/> A	+\$125	145	195	84	98
		3048	3051L	3051R	3060	3072
KH	light tone <input type="checkbox"/> A	+\$125	130	130	156	183
KR	inner tone light <input type="checkbox"/> A	+\$125	130	130	156	183
		3096				
KH	light tone <input type="checkbox"/> A					+\$233
KR	inner tone light <input type="checkbox"/> A					+\$233

Heavy-Duty Storage Work Surface

C2112.
C2122.



Product Information

Description

This rectangular work surface hangs from a module or wall strips and has a top-loading capacity of up to 400 pounds. It has a C frame adapter to hold C frames and drawers. Attachment hardware is included. Shipped knocked down.

The work surface has a laminate, Chemsurf®, or black resin top. The laminate or Chemsurf support arm finish matches the surface finish. 51"-wide laminate and Chemsurf surfaces have a 3" overhang on the left or right side to fill the gap created by a mid-run support, C-series 3-way connector, or interface connector. 96"-wide laminate and Chemsurf surfaces include 1 additional support bracket.

Notes

Laminate or Chemsurf work surface supports following quantities of A-, B-, or C-size drawers; A-style pencil drawers; C frames; or storage units:

- 24" wide: 1 drawer, 1 C frame
- 30" wide: 1 drawer, 1 C frame or 1 storage unit
- 48", 51", and 60" wide: 2 drawers, 2 C frames or 2 storage units
- 72" wide: 3 drawers, 3 C frames or 3 storage units
- 96" wide: 4 drawers, 4 C frames or 4 storage units

Resin work surface supports following quantities of A-, B-, or C-size drawers; A-style pencil drawers; or C frames:

- 24" and 30" wide: 1 drawer, 1 C frame
- 48" wide: 2 drawers, 2 C frames
- 72" wide: 3 drawers, 3 C frames


Order storage components separately:

- A-style pencil drawer (AO480.)
- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (SM100., SM101., SM102., or SM103.)

When hanging work surface from wall strips, order recommended work surface backplash (C2910.) separately.

Order optional keyboard tray (G7718.) separately for laminate or Chemsurf work surface.

Dimensions



Heavy-Duty Storage Work Surface

continued

Co/Struc® Work Surfaces

Specification Information

Step 1.

C21

A

Step 2. Surface Material

12. laminate or Chemsurf® top

A

22. resin top

A

Step 3. Depth

24 24" deep

A

30 30" deep

A

Step 4. Width

For laminate or Chemsurf® top (12.)

24 24" wide

A

30 30" wide

A

48 48" wide

A

51L 51" wide, left overhang

A

51R 51" wide, right overhang

A

60 60" wide

A

72 72" wide

A

96 96" wide

A

For resin top (22.)

24 24" wide

A

30 30" wide

A

48 48" wide

A

72 72" wide

A

Prices for Steps 1-4.

	24	30	48	51L	51R	60	72
C2112. 24	\$400	448	585	599	599	684	786
30	\$500	563	712	725	725	831	943
C2122. 24	\$950	1021	1344	—	—	—	1863
30	\$1100	1200	1640	—	—	—	2280
							96
C2112. 24							\$1024
30							\$1215
C2122. 24							—
30							—

Step 5.

Surface Finish: Laminate Top

For laminate or Chemsurf® top (12.)

		2424	2430	2448	2451L	2451R
HF	inner tone light <input type="checkbox"/>	+\$0	0	0	0	0
LT	light tone <input type="checkbox"/>	+\$0	0	0	0	0
		2460	2472	2496	3024	3030
HF	inner tone light <input type="checkbox"/>	+\$0	0	0	0	0
LT	light tone <input type="checkbox"/>	+\$0	0	0	0	0
		3048	3051L	3051R	3060	3072
HF	inner tone light <input type="checkbox"/>	+\$0	0	0	0	0
LT	light tone <input type="checkbox"/>	+\$0	0	0	0	0
						3096
HF	inner tone light <input type="checkbox"/>					+\$0
LT	light tone <input type="checkbox"/>					+\$0

Surface Finish: Chemsurf® Top

For laminate or Chemsurf® top (12.)


		2424	2430	2448	2451L	2451R
KH	light tone <input type="checkbox"/>	+\$85	95	128	130	130
KR	inner tone light <input type="checkbox"/>	+\$85	95	128	130	130
		2460	2472	2496	3024	3030
KH	light tone <input type="checkbox"/>	+\$148	168	220	75	117
KR	inner tone light <input type="checkbox"/>	+\$148	168	220	75	117
		3048	3051L	3051R	3060	3072
KH	light tone <input type="checkbox"/>	+\$159	159	159	178	202
KR	inner tone light <input type="checkbox"/>	+\$159	159	159	178	202
						3096
KH	light tone <input type="checkbox"/>					+\$265
KR	inner tone light <input type="checkbox"/>					+\$265

Support Finish

For resin top (22.)

HF	inner tone light <input type="checkbox"/>	+\$0
LT	light tone <input type="checkbox"/>	+\$0

Heavy-Duty Option Work Surface C2113. C2123.



Product Information

Description

This rectangular work surface hangs from a module or wall strips and has a top-loading capacity of up to 400 pounds. It can hold drawers and an option work surface C frame adapter kit. Attachment hardware is included. Shipped knocked down.

The work surface has a laminate, Chemsurf®, or resin top. The laminate or Chemsurf support arm finish matches the surface finish. 51"-wide laminate or Chemsurf surfaces have a 3" overhang on the left or right side to fill the gap created by a mid-run support, C-series 3-way connector, or interface connector. 96"-wide laminate or Chemsurf surfaces include 1 additional support bracket.

Notes

To hang components below work surface, order option work surface C frame adapter kit (C2193.) separately.

Work surface supports following quantities of A-style pencil drawers or A-, B-, or C-size drawers:

- 24" and 30" wide: 1 drawer
- 48", 51", and 60" wide: 2 drawers
- 72" wide: 3 drawers
- 96" wide: 4 drawers


Order storage components separately:

- A-style pencil drawer (AO480.)
- Drawer (CO207, CO208, or CO209)

When hanging work surface from wall strips, order recommended work surface backplash (C2910.) separately.

Order optional keyboard tray (G7718.) separately for laminate or Chemsurf work surface.

Dimensions



continued

Co/Struc® Work Surfaces

Step 1.

C21 A

13. laminate or Chemsurf® top

23. resin top ☐ A

24 24" deep

30 30" deep ☐ A

For laminate or Chemsurf® top (13.)

24 24" wide A

30 30" wide ☐ A

48 48" wide ☐ A

51L 51" wide, left overhang A

51R 51" wide, right overhang A

60 60" wide

72 72" wide

96 96" wide

For resin top (23.)

24 24" wide ☐ A

30 30" wide ☐ A

48 48" wide ☐ A

72 72" wide

Prices for Steps 1-4.

	24	30	48	51L	51R	60	72
C2113. 24	\$335	359	463	477	477	549	633
30	\$400	426	570	589	589	674	792

C2123. 24	\$950	1021	1344	—	—	—	1863
30	\$1100	1200	1640	—	—	—	2280

96

C2113. 24	\$850
30	\$1014

C2123. 24	—
30	—

Surface Finish: Laminate Top

For laminate or Chemsurf® top (13.)

			2424	2430	2448	2451L	2451R
HF	inner tone light	A	+\$0	0	0	0	0
LT	light tone	A	+\$0	0	0	0	0
			2460	2472	2496	3024	3030
HF	inner tone light	A	+\$0	0	0	0	0
LT	light tone	A	+\$0	0	0	0	0
			3048	3051L	3051R	3060	3072
HF	inner tone light	A	+\$0	0	0	0	0
LT	light tone	A	+\$0	0	0	0	0
							3096
HF	inner tone light	A					+\$0
LT	light tone	A					+\$0

Surface Finish: Chemsurf® Top

For laminate or Chemsurf® top (13.)

			2424	2430	2448	2451L	2451R
KH	light tone	A	+\$75	91	117	121	121
KR	inner tone light	A	+\$75	91	117	121	121
			2460	2472	2496	3024	3030
KH	light tone	A	+\$135	160	214	80	108
KR	inner tone light	A	+\$135	160	214	80	108
			3048	3051L	3051R	3060	3072
KH	light tone	A	+\$144	148	148	170	199
KR	inner tone light	A	+\$144	148	148	170	199
							3096
KH	light tone	A					+\$253
KR	inner tone light	A					+\$253


Support Finish

For resin top (23.)

HF	inner tone light	<input type="text" value="A"/>	+\$0
LT	light tone	<input type="text" value="A"/>	+\$0

Rail-Hanging Work Surface

C2410.



Product Information


Description

This 24"-deep rectangular work surface hangs on a standard rail. It has predrilled holes for hanging drawers, suspended pedestals, or keyboard trays. The work surface can be installed, relocated, or removed without tools and without disturbing other rail components. The support arm finish matches the surface finish.

Notes

Order optional keyboard tray (G7718.) separately.

Dimensions



Specification Information

Step 1.

C2410.24

Step 2. Width

48	48" wide
72	72" wide

Prices for Steps 1-2.

C2410.24	48	\$523
	72	\$711


Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Work Surfaces

C Frame Support

C2192.



Product Information

Description

This assembly attaches to a support panel, a module, or wall strips and provides support for a C frame and components hanging on the C frame. Attachment hardware is included. Shipped knocked down.





The assembly holds the following quantities of C frames:

- 24" and 30" wide: 1 C frame
- 48" and 60" wide: 2 C frames
- 72" wide: 3 C frames

Notes

For 51"-wide work surface, order 48"-wide C frame support.

Dimensions



Specification Information

Step 1.

C2192.

Step 2. Work Surface Depth

24 for 24"-deep work surface

30 for 30"-deep work surface

Step 3. Width

24 24" wide

30 30" wide

48 48" wide

60 60" wide

72 72" wide


Prices for Steps 1-3.

	24	30	48	60	72
C2192. 24	\$200	235	279	334	387
30	\$285	294	337	393	448

Step 4. Surface Finish

HF	inner tone light <input type="text" value="A"/>				+\$0
LT	light tone <input type="text" value="A"/>				+\$0

C Frame Adapter Kit, Option Work Surface C2193.



Product Information

Description

This kit attaches to a heavy-duty option work surface and allows C frames to hang below the work surface. Attachment hardware is included.

The kit holds the following quantities of C frames or storage units:

- 24" wide: 1 C frame
- 30" wide: 1 C frame or storage unit
- 48" and 60" wide: 2 C frames or storage units
- 72" wide: 3 C frames or storage units
- 96" wide: 4 C frames or storage units


Notes

For 51"-wide work surface, order 48"-wide kit.

Order storage components separately:

- C frame (CO205)
- Storage unit (SM100., SM101., SM102., or SM103.)

Dimensions



Specification Information

Step 1.

C2193.

Step 2. Width

24	24" wide
30	30" wide
48	48" wide
60	60" wide
72	72" wide
96	96" wide


Prices for Steps 1-2.

C2193. 24	\$110
30	\$130
48	\$165
60	\$199
72	\$234
96	\$302

Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Work Surface BacksplashC2910.




Product Information

Description

This plastic barrier attaches to the back edge of a C-series rectangular or corner work surface or to a process table and creates a backsplash. Package contains 6.

Dimensions



Specification Information

Step 1.

C2910.

Step 2. Width	
24	24" wide
30	30" wide
48	48" wide
51	51" wide
60	60" wide
72	72" wide
96	96" wide

Prices for Steps 1-2.	
C2910. 24	\$88
30	\$110
48	\$172
51	\$182
60	\$204
72	\$236
96	\$301


Step 3. Surface Finish		
For 24" wide (24), 30" wide (30), 48" wide (48), or 72" wide (72)		
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
For 51" wide (51), 60" wide (60), or 96" wide (96)		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Work Surfaces




Backsplash Filler

C2920.



Product Information
Description
This component attaches to a Casework work surface backsplash. It fills the gap created when a work surface is placed against terminal panels or wall strips. Package contains 6. Attachment hardware is included.
Notes
Order work surface (CG210.) or corner work surface (CG217.) separately.
Dimensions



Specification Information
Step 1.
C2920.96
Step 2. Surface Finish
HF inner tone light
LT light tone


\$120

+\$0

+\$0

Co/Struc® Work Surfaces

Heavy-Duty Corner Work Surface C2215.



Product Information

Description

This 48"-wide work surface hangs from modules, C-series support panels, or wall strips and provides additional top-loading capacity. It includes 1 drawer bearer package. The work surface cannot attach to Action Office® Series 1 or Series 2 connectors. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top. The support arm finish matches the surface finish.

Notes


When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard tray (G7718.) separately.


Dimensions



TOP



TOP




FRONT

Heavy-Duty Corner Work Surface

continued

Specification Information			
Step 1.			
C2215. <input type="checkbox"/> A			
Step 2. Depth			
2448	24" deep	<input type="checkbox"/> A	
3048	30" deep	<input type="checkbox"/> A	
Prices for Steps 1-2.			
C2215. 2448			\$659
3048			\$765
Step 3. Surface Finish			
Laminate Top			
		2448	3048
HF	inner tone light	<input type="checkbox"/> A	
		+\$0	0
LT	light tone	<input type="checkbox"/> A	
		+\$0	0
Chemsurf® Top			
		2448	3048
KH	light tone	<input type="checkbox"/> A	
		+\$167	194
KR	inner tone light	<input type="checkbox"/> A	
		+\$167	194

Heavy-Duty Corner Work Surface, C2216.
Action Office® Series 1 or 2 90°
Corner Connector



Product Information

Description

This 48"-wide work surface hangs from panels or wall strips connected by an Action Office® Series 1 or 2 90° corner connector. It provides additional top-loading capacity and includes 1 drawer bearer package. The support arm finish matches the surface finish. The work surface cannot attach to a C-series 90° corner connector. Attachment hardware is included.


The work surface has a laminate or Chemsurf® top.

Notes


When hanging work surface from wall strips, order recommended work surface backsplash (C2910.) separately.

Order optional keyboard tray (G7718.) separately.


Dimensions



TOP



TOP



FRONT


Heavy-Duty Corner Work Surface,
Action Office® Series 1 or 2 90°
Corner Connector *continued*

Specification Information			
Step 1.			
C2216. <input type="checkbox"/> A			
Step 2. Depth			
2448	24" deep	<input type="checkbox"/> A	
3048	30" deep	<input type="checkbox"/> A	
Prices for Steps 1-2.			
C2216. 2448			\$659
3048			\$765
Step 3. Surface Finish			
Laminate Top			
		2448	3048
HF	inner tone light <input type="checkbox"/> A	+\$0	0
LT	light tone <input type="checkbox"/> A	+\$0	0
Chemsurf® Top			
		2448	3048
KH	light tone <input type="checkbox"/> A	+\$167	194
KR	inner tone light <input type="checkbox"/> A	+\$167	194

Process Table

C4115.

C4215.



Product Information

Description

This freestanding table adjusts from 29" to 38" high. It has stretchers to support C frames and storage units and drawer bearers to support drawer storage. The support finish matches the surface finish. Shipped knocked down.

The table has a laminate or Chemsurf® top. Legs have adjustable glides or 4" locking casters.

Notes

Table supports the following pairs of drawer bearers (each pair has a left and right bearer):

- 48"-wide table: 2 pair
- 60"-wide table: 2 pair
- 72"-wide, 24"-deep table: 3 pair
- 72"-wide, 30"-deep table: 3 pair
- 72"-wide, 36"-deep table: 4 pair (2 pair on each side)

Table dimensions determine number of C frames and storage units and their configuration:

Table Dimensions—Quantity and Configuration

24"d × 48"w—2 C frames or storage units side by side
 24"d × 60"w—2 C frames or storage units side by side
 24"d × 72"w—3 C frames or storage units side by side
 30"d × 48"w—2 C frames or storage units side by side
 30"d × 60"w—2 C frames or storage units side by side
 30"d × 72"w—3 C frames or storage units side by side
 36"d × 72"w—4 C frames or storage units, 2 on each side

Order storage components separately:





- C frame (CO205)
- Drawer (CO207, CO208, or CO209)
- Storage unit (SM100., SM101., SM102., or SM103.)

For table less than 32" high, C frames hang from top stretcher only.

For table less than 32" high, storage units cannot mount.

Order optional keyboard tray (G7718.) separately.

Dimensions



30" OR 36" DEEP
SIDE

24" DEEP
SIDE

FRONT

Process Table *continued*

Specification Information

Step 1.

C4 ☐

Step 2. Casters/Glides

215. casters ☐

115. glides ☐

Step 3. Depth

24 24" deep ☐

30 30" deep ☐

36 36" deep ☐

Step 4. Width

For 24" deep (24) or 30" deep (30)

48 48" wide ☐

60 60" wide ☐

72 72" wide ☐

For 36" deep (36)

72 72" wide ☐

Prices for Steps 1-4.

	48	60	72
C4215. 24	\$1273	1402	1531
30	\$1330	1467	1600
36	—	—	\$1714
C4115. 24	\$1145	1275	1403
30	\$1203	1340	1473
36	—	—	\$1587

Step 5. Surface Finish

Laminate Top


	2448	2460	2472	3048	3060
HF inner tone light <input type="checkbox"/>	+\$0	0	0	0	0
LT light tone <input type="checkbox"/>	+\$0	0	0	0	0
				3072	3672
HF inner tone light <input type="checkbox"/>				+\$0	0
LT light tone <input type="checkbox"/>				+\$0	0

Chemsurf® Top


	2448	2460	2472	3048	3060
KH light tone <input type="checkbox"/>	+\$172	193	210	183	200
KR inner tone light <input type="checkbox"/>	+\$172	193	210	183	200
				3072	3672
KH light tone <input type="checkbox"/>				+\$221	239
KR inner tone light <input type="checkbox"/>				+\$221	239

Sink-Adaptable Module

C2510.




Product Information	
Description	
This component hangs from a module, a support panel, or wall strips. It has a Chemsurf® top, extruded PVC backsplash, and holds a 12"-high sink. Attachment hardware is included. Shipped knocked down.	
The module accepts the following cutout dimensions:	
Depth:	
• 24": 18 ¹ / ₄ "-deep maximum cutout	
• 30": 24 ¹ / ₄ "-deep maximum cutout	
Width:	
• 24": 19 ⁵ / ₈ "-wide maximum cutout	
• 30": 25 ⁵ / ₈ "-wide maximum cutout	
• 48": 43 ⁵ / ₈ "-wide maximum cutout	
• 60": 55 ⁵ / ₈ "-wide maximum cutout	
Notes	
Sink must be customer supplied and field installed.	
Dimensions	



Specification Information					
Step 1.					
C2510. <input type="text" value="A"/>					
Step 2. Depth					
24	24" deep	<input type="text" value="A"/>			
30	30" deep	<input type="text" value="A"/>			
Step 3. Width					
24	24" wide	<input type="text" value="A"/>			
30	30" wide	<input type="text" value="A"/>			
48	48" wide	<input type="text" value="A"/>			
60	60" wide	<input type="text" value="A"/>			
Prices for Steps 1-3.					
		24	30	48	60
C2510. 24		\$610	758	804	901
30		\$733	909	965	1082
Step 4. Surface Finish					
KH	light tone	<input type="text" value="A"/>	+\$0		
KR	inner tone light	<input type="text" value="A"/>	+\$0		

Sink Unit for Storage


C2511.



Product Information
Description
This 24"-deep, 48"-wide surface hangs from 2 24"-wide panels or freestanding modules or from 3 wall strips. It has a Chemsurf® top and extruded PVC backsplash. The surface supports drawer storage on the left or right side and a sink unit on the other side. The maximum cutout for the sink is 18 ¹ / ₄ " deep × 19 ⁵ / ₈ " wide. The surface includes a front panel and side panels to conceal the bottom of the sink and left and right work surface supports. Attachment hardware is included. Shipped knocked down.

Notes
Order storage components for open space opposite sink unit separately:
<ul style="list-style-type: none">• Drawer bearers (CO228 or CO396)• Drawer caddy with casters (CO692. or CO693.)• Gripper package (CO266) and C frame (CO205) for Co/Struc® drawers• L cart (CO212) and single counter top (CO555) for C frames and drawers when surface is mounted 38" above floor• Storage case (SM400.2422) and 24"-wide Co/Struc adapter rail (C1610.24)
Sink must be customer supplied and field installed.
SM100. components cannot mount under sink unit.

Dimensions




Specification Information
Step 1.
C2511.2448 [A] \$800
Step 2. Surface Finish
KH light tone [A] +\$0
KR inner tone light [A] +\$0



Co/Struc® Work Surfaces

ADA Sink Support, Systems Mounted

C2520.




Product Information
Description This sink support hangs from a module or wall strips. It supports a 5½"-deep sink and a Casework work surface or alternative material top. The sink support has a 4"-high laminate front trim and a removable front panel to allow access to plumbing. Attachment hardware is included. Shipped knocked down.
Notes Order Casework work surface (CG210.) separately. Sink must be customer supplied and field installed. A Just SL-ADA-1921-A-GR sink, with drain located in rear center, fits in a 24" × 24" sink support.
Dimensions



Specification Information					
Step 1.					
C2520. <input type="text" value="A"/>					
Step 2. Depth					
24	24" deep	<input type="text" value="A"/>			
30	30" deep	<input type="text" value="A"/>			
Step 3. Width					
24	24" wide	<input type="text" value="A"/>			
30	30" wide	<input type="text" value="A"/>			
48	48" wide	<input type="text" value="A"/>			
60	60" wide	<input type="text" value="A"/>			
Prices for Steps 1-3.					
		24	30	48	60
C2520. 24		\$800	900	1000	1100
30		\$850	950	1100	1200
Step 4. Base Finish					
HF	inner tone light	<input type="text" value="A"/>	+\$0		
LT	light tone	<input type="text" value="A"/>	+\$0		
Step 5. Edge Finish					
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.					
Price Category 1					+\$0
Step 6. Front Finish					
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.					
Price Category 1					+\$0

ADA Sink Support for Storage, Systems Mounted C2521.



Product Information

Description

This sink support hangs from a module or wall strips. It supports a 5½"-deep sink and a Casework work surface or alternative material top. The sink support has a 4"-high laminate front trim and a removable front panel to allow access to plumbing. It includes a pair of work surface support arms to hold drawer storage on the left or right side and a sink unit on the other side. Attachment hardware is included. Shipped knocked down.

Notes

Order Casework work surface (CG210.) separately.


Order storage components for open space opposite sink separately:

- Drawer bearer package (CO228) or lockable drawer bearer (CO396)
- Drawer caddy with casters (CO692. or CO693.)
- Gripper package (CO266) and C frame (CO205) for Co/Struc® drawers
- Storage case (SM400.2422) and 24"-wide Co/Struc adapter rail (C1610.24)

Sink must be customer supplied and field installed. A Just SL-ADA-1921-A-GR sink, with drain located in rear center, fits in a 24" × 24" sink support.

SM100. components cannot mount under sink unit.

Dimensions



Specification Information

Step 1.

C2521.2448 \$1000

Step 2. Base Finish

HF inner tone light +\$0

LT light tone +\$0

Step 3. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 4. Front Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Storage Unit

SM100.
SM101.
SM102.
SM103.




Product Information

Description
This cabinet hangs below a heavy-duty storage work surface or process table. It includes a shelf that adjusts in 3" increments. The cabinet cannot hang below a 24"-wide work surface. Shipped knocked down.
The cabinet is available without a door, with a door, or with a door and drawer slot. Doors attach on the left or right side.


Notes
Order heavy-duty storage work surface (C2112.) or process table (C4115. or C4215.) separately.
Order additional shelves (SM104.) separately.
For door with 3" drawer slot, order A-size drawer (CO207) separately.
For door with 2 3" drawer slots, order 2 A-size drawers (CO207) or 1 B-size drawer (CO208) separately.

Dimensions




Co/Struc® Storage

Specification Information		
Step 1.		
SM10		
Step 2. Door		
0.1622	no door	
1.1622	full-height door	
2.1622	door with 3" drawer slot	
3.1622	door with 2 3" drawer slots	
Prices for Steps 1-2.		
SM100.1622		\$655
SM101.1622		\$750
SM102.1622		\$787
SM103.1622		\$811
Step 3. Cabinet Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
Step 4. Door Finish		
<i>For full-height door (1.1622), door with 3" drawer slot (2.1622), or door with 2 3" drawer slots (3.1622)</i>		
HF	inner tone light	+\$0
LT	light tone	+\$0
W6	brittany blue	+\$0
W8	victorian teal	+\$0



Shelf for Storage Unit

SM104.




Product Information

Description

This shelf fits into the slots of a storage unit that mounts below the work surface. It provides additional shelf space.

Dimensions



TOP

Specification Information

Step 1.

SM104.1520	\$150
-------------------	-------


Step 2. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage

Overhead Storage Unit

SM105.




Product Information

Description
This unit hangs on a module, support panel, or wall strips above a work surface. It has sliding glass doors. Shipped knocked down.

Notes
To divide interior space, order add-on shelf with lip (C3019.13) separately.
Order optional task light separately:
• Utility task light (G6132.)
• Energy-efficient task light (G6120.)

Dimensions



Specification Information

Step 1.

SM105.16

Step 2. Width

30 30" wide
48 48" wide

Prices for Steps 1-2.

SM105.16 30 \$792
48 \$928


Step 3. Surface Finish

HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Storage

Add-On Shelf with Lip

C3019.




Product Information

Description
This shelf fits into an overhead storage unit, A- or B-style shelf, or B-style rail-hanging shelf to add a second shelf. Attachment hardware is included.

Notes

- Order support components separately:
- A-style (C3010.) or B-style (C3011.) shelf
 - B-style (CO521.) rail-hanging shelf
 - Overhead storage unit (SM105.)

Dimensions



Specification Information

Step 1.
C3019.

Step 2. Depth

- 13** 13" deep
16 16" deep

Step 3. Width

- 24** 24" wide
30 30" wide
36 36" wide
42 42" wide
48 48" wide
60 60" wide

Prices for Steps 1-3.


	24	30	36	42	48	60
C3019. 13	\$71	75	84	84	97	109
16	\$105	112	123	129	150	189

Step 4. Surface Finish

- HF** inner tone light +\$0
LT light tone +\$0

Shelf Unit

SM303.
SM305.



Product Information

Description

This 3- or 5-shelf unit hangs from a module, panel, or wall strips. It has 9"- or 16"-deep shelves that adjust in 1/4" increments. 9"-deep shelves are flat or slanted at 2 angles; 16"-deep shelves are flat or slanted at 4 angles for gravity-feed dispensing. Shipped knocked down.

Notes

For additional shelf, order add-on shelf (SM900.) separately.

To cover shelf unit with dust cover, order add-on shelf (SM900.) and place in top position of unit.

To enclose back of unit, order back panel (SM940.) separately.

To hang unit on architectural wall and enclose back of unit, order overhead cabinet back panel (CG393.) separately.


To enclose front of unit, order overhead cabinet door (CG390., CG391., or CG392.) separately.

Order optional task light for 16"-deep unit separately:

- Utility task light (G6132.)
- Energy-efficient task light (G6120.)

To attach task light to slanted shelf, order light bracket (SM910.) separately.

Dimensions



Specification Information

Step 1.

SM30

Step 2. Shelves/Height

3. 3 shelves and 24" high**5.** 5 shelves and 34" high

Step 3. Depth

09 9" deep**16** 16" deep

Step 4. Width

24 24" wide**30** 30" wide**36** 36" wide**48** 48" wide

Prices for Steps 1-4.


	24	30	36	48
SM303.09	\$244	258	264	278
16	\$306	323	341	383
SM305.09	\$407	421	438	467
16	\$489	522	553	613

Step 5. Side Panel/Shelf Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Shelf Unit, Diagonal Corner

SM313.
SM315.



Product Information

Description

This unit hangs from 24"-wide panels connected by an Action Office® 2-way 90° connector or wall strips used in a corner position. Shelves adjust in 1 1/4" increments and are flat or slanted at 4 angles for gravity-feed dispensing. Shipped knocked down.

Notes

For additional shelf in center section, order 24"-wide add-on shelf (SM900.1624). For additional shelf in side section, order 14"-wide add-on shelf (SM900.1614).




To cover back of unit, order back panel (SM940.) separately.

Order optional 24"-wide task light for center section separately:

- Utility task light (G6132.24)
- Energy-efficient task light (G6120.24)

To attach task light to slanted center shelf, order light bracket (SM910.) separately.

Dimensions



Shelf Unit, Diagonal Corner *continued*

Specification Information

Step 1.

SM31

Step 2. Shelves/Height

3.1648 3 shelves and 24" high

5.1648 5 shelves and 34" high

Prices for Steps 1-2.

SM313.1648 \$915

SM315.1648 \$1460

Step 3. Side Panel/Shelf Finish


HF inner tone light +\$0

LT light tone +\$0

Co/Struc® Storage

Shelf Unit, 90° Corner

SM329.



Product Information


Description

This 72"-high shelf unit hangs from 24"- or 48"-wide panels connected by 1 or 2 Action Office® 2-way 90° connectors or hangs from wall strips used in 1 or 2 90° corners. Shelves adjust in 1 1/4" increments and are flat or slanted at 4 angles for gravity-feed dispensing. Shipped knocked down.

Notes

For additional shelf, order add-on shelf (SM900.1631) separately.
To cover shelf unit with dust cover, order add-on shelf (SM900.1631) and place in top position of unit.

Dimensions



Specification Information

Step 1.

SM329.1648 \$2073


Step 2. Side Panel/Shelf Finish

HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Storage

Add-On Shelf

SM900.




Product Information

Description
This shelf fits inside a shelf unit, diagonal corner shelf unit, or 90° corner shelf unit. It is used as an additional shelf or dust cover.

Dimensions

Co/Struc® Storage



Specification Information

Step 1.
SM900.

Step 2. Depth
09 9" deep
16 16" deep

Step 3. Width

For 9" deep (09)
24 22³/₄" wide
30 28³/₄" wide
36 34³/₄" wide
48 46³/₄" wide


For 16" deep (16)
14 13¹/₂" wide
24 22³/₄" wide
30 28³/₄" wide
31 30¹/₂" wide
36 34³/₄" wide
48 46³/₄" wide

Prices for Steps 1-3.

	14	24	30	31	36	48
SM900.09	—	\$58	60	—	63	78
16	\$57	69	75	78	80	92


Step 4. Shelf Finish

HF	inner tone light	+\$0
LT	light tone	+\$0



Shelf Label Clip


C0398



Product Information

Description
This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.

Dimensions




Specification Information

Step 1.
C039827 \$94

Back Panel

SM940.




Product Information

Description
This panel attaches to a 3- or 5-shelf unit or diagonal corner shelf unit and closes the back. The shelf unit's top and bottom shelves must be flat.

Dimensions

Co/Struc® Storage



Specification Information

Step 1.
SM940.

Step 2. Height
24 24" high
34 34" high


Step 3. Width
24 22³/₄" wide
30 28³/₄" wide
36 34³/₄" wide
48 46³/₄" wide
60 58³/₄" wide

Prices for Steps 1-3.					
	24	30	36	48	60
SM940.24	\$128	138	148	177	174
34	\$155	169	183	212	206

Step 4. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Light Bracket

SM910.




Product Information

Description
These brackets fasten a task light to the bottom of an overhead cabinet, a flat or slanted shelf in a 3- or 5-shelf unit, or a center shelf in a diagonal corner shelf unit. Package contains 1 pair.
The task light attaches to a 9"-deep shelf in the same position as the shelf, either flat or slanted. The task light attaches to a 16"-deep slanted shelf in a flat position. The 16"-deep flat shelf does not require brackets for attaching the task light.

Notes
Order task light separately:
• For 9"-deep shelf, order utility task light (G6132.)
• For 16"-deep shelf, order energy-efficient task light (G6120.)

Dimensions



Specification Information

Step 1.
SM910.

Step 2. Usage
09 for 9" deep flat or slanted shelf
16 for 16" deep slanted shelf


Prices for Steps 1-2.
SM910.09 \$37
16 \$37

Step 3. Surface Finish
For 16" deep slanted shelf (16)
HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Storage

Overhead Cabinet Door

CG390.
CG391.



Product Information

Description

This 24"-wide door encloses a shelf unit. Attachment hardware is included.


Notes

For 24"-high door, order 3-shelf unit (SM303.0924 or SM303.1624) separately.

For 34"-high door, order 5-shelf unit (SM305.0924 or SM305.1624) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG39

Step 2. Door Swing

- 0.** left-hand door swing
- 1.** right-hand door swing

Step 3. Height

- 2424** 24" high
- 3424** 34" high

Step 4. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

Prices for Steps 1-4.

	N	L	K
CG390. 2424	\$129	169	169
3424	\$178	218	218
CG391. 2424	\$129	169	169
3424	\$178	218	218

Step 5. Pull Option/Finish

NN	none		+\$0
BK	black		+\$0
HF	inner tone light		+\$0
LT	light tone		+\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.


Price Category 1 +\$0

Step 7. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0
Price Category 2 +\$15

Overhead Cabinet Double Doors CG392.



Product Information

Description

These doors enclose a shelf unit. Attachment hardware is included.


Notes

For 24"-high doors, order 3-shelf unit (SM303.) separately.

For 34"-high doors, order 5-shelf unit (SM305.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG392.

Step 2. Height

24 24" high

34 34" high

Step 3. Width

30 30" wide

36 36" wide

48 48" wide

Step 4. Lock

N no lock

L keyed differently

K keyed alike

Prices for Steps 1-4.

	N	L	K
CG392. 24 30	\$177	217	217
36	\$207	247	247
48	\$256	296	296
34 30	\$230	270	270
36	\$271	311	311
48	\$356	396	396

Step 5. Pull Option/Finish

NN none +\$0

BK black +\$0

HF inner tone light +\$0

LT light tone +\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 7. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Price Category 2 +\$25

Overhead Cabinet Back PanelCG393.




Product Information
Description
This panel encloses the back of a 3- or 5-shelf unit and allows the shelf unit to hang on an architectural wall.
Notes
Order 3-shelf unit (SM303.) or 5-shelf unit (SM305.) separately.
Dimensions

Specification Information
Step 1.
CG393.
Step 2. Height
2424" high
3434" high
Step 3. Width
2424" wide
3030" wide
3636" wide
4848" wide

Prices for Steps 1-3.	24	30	36	48
CG393. 24	\$110	135	180	210
34	\$120	140	185	220


Step 4. Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0



Co/Struc® Storage

Storage Case

SM400.



Product Information

Description

This case mounts under a heavy-duty rectangular or corner work surface. It has 2 side panels, a back panel, a top brace, and a bottom spanner. Interior components can be height adjusted in 1 1/4" increments. The case cannot mount under a 24"-wide work surface or resin work surface. A lockable Co/Struc® drawer cannot mount in the case. Shipped knocked down.

Notes

Order heavy-duty work surface (C2210.) or heavy-duty corner work surface (C2215. or C2216.) separately.

To support case mounted under rectangular work surface, order adapter rail (C1610.); to support case mounted under corner work surface, order corner storage case rail (C1611.) separately.

Order optional storage case door (SM447. or SM448.) separately.

Order optional storage case fixed shelf (SM411.) separately.

Order optional drawer separately:


- Storage case drawer (SM433. or SM434.)
- Storage case bottle drawer (SM436.)

For top slot of case, specify 16"-deep drawer, shelf, or Co/Struc drawer to clear support beam under heavy-duty work surface.

To enclose bottom of case, order 24"- or 30"-deep storage case fixed shelf (SM411.2422 or SM411.3022); shelf not needed when drawer is in bottom position.

When ordering Co/Struc drawer for bottom position of case, order 24"- or 30"-deep storage case fixed shelf (SM411.2422 or SM411.3022) to enclose bottom of case.

Dimensions




Co/Struc® Storage


Specification Information		
Step 1.		
SM400.		
Step 2. Depth		
2422	23" deep	
3022	29" deep	
Prices for Steps 1-2.		
SM400.2422		\$456
3022		\$521
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Storage Case Fixed Shelf


SM411.



Product Information
Description
This shelf mounts inside a storage case.
The 16"-deep shelf is flat or slanted at 4 angles for gravity-feed dispensing.
Notes
To clear support beam of heavy-duty work surface, specify 16"-deep shelf for top slot of case.
Dimensions



SIDE




FRONT

Specification Information
Step 1.
SM411.
Step 2. Depth
1622 16" deep
2422 23" deep
3022 29" deep
Prices for Steps 1-2.
SM411.1622 \$69
2422 \$82
3022 \$93
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Storage

Storage Case Pullout Shelf, FlatSM421.



Product Information

Description

This flat, pullout shelf is 21³/₄" deep and mounts inside a storage case.

Notes

Position shelf to allow clearance for 2¹/₂"-high support beam of heavy-duty work surface.

Dimensions

Specification Information

Step 1.


SM421.2422\$128

Step 2. Surface Finish

HFinner tone light+\$0


LTlight tone+\$0

Co/Struc® Storage



Storage Case Drawer

SM433.
SM434.



Product Information

Description


This drawer mounts inside a storage case and is available with or without a lock. The 24"-deep drawer cannot be used in the top position of a storage case.

Notes

Specify fixed shelf (SM411.) to enclose top of lockable drawer in storage case.

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

SM43

Step 2. Depth

- 3. 16" deep
- 4. 22" deep

Step 3. Height

- 0422 3³/₄" high
- 0822 7¹/₂" high
- 1222 11¹/₄" high

Step 4. Lock

For no lock (*), skip this step.

- L keyed differently
- K keyed alike

Prices for Steps 1-4.

	*	L	K
SM433.0422	\$172	289	289
0822	\$201	306	306
1222	\$229	330	330
SM434.0422	\$187	297	297
0822	\$215	319	319
1222	\$244	343	343

Step 5. Front Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
W6	brittany blue	+\$0
W8	victorian teal	+\$0

Step 6. Key Number

For keyed alike (K)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0


Storage Case Drawer *continued*

Co/Struc® Storage

12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			

Storage Case Bottle Drawer

SM436.




Product Information


Description
This drawer mounts in a storage case and has a partial-height front for easy access to contents. It includes 3 adjustable dividers.

Notes
Order additional dividers (SM930.) separately.


Dimensions




SIDE



FRONT



SIDE



FRONT

Specification Information

Step 1.
SM436.

Step 2. Height
0822 6" high
1222 9 7/8" high


Prices for Steps 1-2.
SM436.0822 \$323
1222 \$375

Step 3. Front Finish
HF inner tone light +\$0
LT light tone +\$0
W6 brittany blue +\$0
W8 victorian teal +\$0

Co/Struc® Storage

Storage Case Door

SM447.
SM448.



Product Information

Description

This door attaches to the front of a storage case. Attachment hardware is included.

Hinges for the left-hand door swing are on the right side when facing the door; hinges for the right-hand door swing are on the left side. The door is available with or without a lock.


Notes

Co/Struc®, storage case, and storage case bottle drawers cannot be specified behind doors.

When specifying 20"- or 24"-high lockable door that has drawers above it, specify fixed shelf (SM411.) to prevent access to cabinet below.

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

SM44

Step 2. Door Swing

- 8.** left-hand door swing
- 7.** right-hand door swing

Step 3. Height

- 2022** 18⁷/₈" high
- 2422** 22⁵/₈" high
- 2822** 26¹/₂" high

Step 4. Lock

For no lock (), skip this step.*

- L** keyed differently
- K** keyed alike

Prices for Steps 1-4.

	*	L	K
SM448.2022	\$114	190	190
2422	\$143	202	202
2822	\$172	221	221
SM447.2022	\$114	190	190
2422	\$143	202	202
2822	\$172	221	221

Step 5. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0
W6	brittany blue	+\$0
W8	victorian teal	+\$0

Step 6. Key Number

For keyed alike (K)


01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0

Storage Case Door *continued*

12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			

Co/Struc® Storage


Lockable Drawer BearerSM491.



Product Information
Description
These drawer bearers attach Co/Struc® A-, B-, or C-size drawers in a storage case. They have integral stops to prevent drawers from falling when pulled out. Mounting hardware is included. Finish is black. Package contains 1 pair.
Notes
To lock drawer to lockable drawer bearer, order lockable drawer separately.
Dimensions


Specification Information
Step 1.
SM491.
\$42

Co/Struc® Storage



Bottle Drawer Divider


SM930.



Product Information

Description
This divider fits inside a storage case bottle drawer. Finish is soft white. Package contains 3.

Dimensions




Specification Information

Step 1.		
SM930.		
Step 2. Height		
0824	5 1/4" high	
1224	9 1/8" high	
Prices for Steps 1-2.		
SM930.0824		\$64
1224		\$79

Co/Struc® Storage


Storage Case Rail, CornerC1611.



Product Information
Description
This support rail attaches to wall strips or to 24"-wide panels or modules joined by an Action Office® Series 1 or 2 corner connector or a C-series corner connector. It supports a storage case hung from a 48"-wide, 24"-or 30"-deep heavy-duty corner work surface.
Dimensions


Specification Information
Step 1.
C1611.48
\$172
Step 2. Surface Finish
HF inner tone light
LT light tone

Co/Struc® Storage



Storage Case Leg

SM495.



Product Information


Description

This leg supports a storage case mounted under a heavy-duty work surface. It provides support when the work surface holds heavy loads. The leg adjusts in 1" increments to allow a work surface height of 34" to 38".

Notes

Order storage case (SM400.) separately.

Dimensions




Specification Information

Step 1.

SM495.

\$139

B-Style Rail-Hanging ShelfC0521.



Product Information

Description

Notes

Dimensions

This 15½"-high shelf with a lip hangs from a standard rail or adapter rail. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

For fully enclosed storage, order flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately.

Order optional add-on shelves of equal depth and width separately:

- Add-on shelf (AO521.1524 or AO521.1548)
- Add-on shelf with lip (C3019.)

Order optional task light separately:

- Utility task light (G6132.)
- Energy-efficient task light (G6120.)

Specification Information

Step 1.

Step 2. Depth

Step 3. Width

Prices for Steps 1-3.

Step 4. Surface Finish


C0521.

13	13" deep
16	16" deep

24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide


	24	30	36	42	48
C0521. 13	\$186	191	200	204	215
16	\$258	264	277	284	308

HF	inner tone light				+\$0
LT	light tone				+\$0



A-Style Shelf

C3010.



Product Information

Description

This 15½"-high shelf with a lip hangs from a panel, a freestanding module, or wall strips. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes

For enclosed storage, order A-style flipper door (AO550.) separately.

For fully enclosed storage of 24"- to 48"-wide shelf, order flipper door back panel (A3390.) and A-style flipper door (AO550.) separately.

Order optional add-on shelves of equal depth and width separately:


- Add-on shelf (AO521.1524 or AO521.1548)
- Add-on shelf with lip (C3019.)

For 60"-wide shelf, order 60"-wide task light.

Order optional task light separately:

- Utility task light (G6132.)
- Energy-efficient task light (G6120.)

Dimensions



Specification Information

Step 1.

C3010.

Step 2. Depth

13	13" deep
16	16" deep

Step 3. Width

24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide
60	60" wide

Prices for Steps 1-3.


	24	30	36	42	48	60
C3010. 13	\$123	127	135	139	148	162
16	\$172	178	189	195	215	258

Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

B-Style Shelf

C3011.



Product Information

Description

This 15½"-high shelf with a lip hangs from a panel, a freestanding module, or wall strips. It can hold up to 2 add-on shelves. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.

Notes

For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.

For fully enclosed storage of 24"- to 48"-wide shelf, order flipper door back panel (A3390.) and B-style flipper door (A3310., A3312., or A3313.) separately.

Order optional add-on shelves of equal depth and width separately:


- Add-on shelf (AO521.1524 or AO521.1548)
- Add-on shelf with lip (C3019.)

For 60"-wide shelf, order 60"-wide task light.

Order optional task light separately:


- Utility task light (G6132.)
- Energy-efficient task light (G6120.)

Dimensions




SIDE

TOP



SIDE

TOP




FRONT

B-Style Shelf *continued*



Specification Information						
Step 1.						
C3011.						
Step 2. Depth						
13	13" deep					
16	16" deep					
Step 3. Width						
24	24" wide					
30	30" wide					
36	36" wide					
42	42" wide					
48	48" wide					
60	60" wide					
Prices for Steps 1-3.						
		24	30	36	42	48 60
C3011. 13		\$121	126	133	138	147 162
16		\$172	178	189	195	215 259
Step 4. Surface Finish						
HF	inner tone light					+\$0
LT	light tone					+\$0

B-Style Shelf with Sloped TopC3091.




Product Information
Description
This 16"-deep shelf hangs from a panel, a freestanding module, or wall strips and has a sloped metal top. A filler piece is included for enclosed storage applications. The underside of the shelf accepts a task light. Attachment hardware is included. Shipped knocked down.
Notes
For enclosed storage, order B-style flipper door (A3310., A3312., or A3313.) separately.
To enclose back of shelf, order flipper door back panel (A3390.) separately.
Order optional task light separately:
<ul style="list-style-type: none">• Utility task light (G6132.)• Energy-efficient task light (G6120.)
Dimensions

Specification Information
Step 1.
C3091.16
Step 2. Width
2424" wide
3030" wide
4242" wide
4848" wide
Prices for Steps 1-2.
C3091.16 24\$345
30\$360
42\$420
48\$440
Step 3. Surface Finish
HFinner tone light+\$0
LTlight tone+\$0



Add-On Shelf with Lip

C3019.



Product Information

Description


This shelf fits into an overhead storage unit, A- or B-style shelf, or B-style rail-hanging shelf to add a second shelf. Attachment hardware is included.

Notes

Order support components separately:

- A-style (C3010.) or B-style (C3011.) shelf
- B-style (CO521.) rail-hanging shelf
- Overhead storage unit (SM105.)

Dimensions



Specification Information

Step 1.

C3019.

Step 2. Depth

13	13" deep
16	16" deep

Step 3. Width


24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide
60	60" wide

Prices for Steps 1-3.

	24	30	36	42	48	60
C3019. 13	\$71	75	84	84	97	109
16	\$105	112	123	129	150	189


Step 4. Surface Finish

HF	inner tone light					+\$0
LT	light tone					+\$0



Shelf Label Clip

C0398



Product Information


Description
This 3"-wide clip fastens to the edge of a metal shelf and holds content labels. Stick-on labels can be applied to the front of the clip or slide-in labels can be inserted from the back. Package contains 36 clear clips.


Specification Information

Step 1.
C039827 \$94

Dimensions

Co/Struc® Storage





Product Information

Description

This enclosed storage unit hangs on rails for stationary storage, or it is used with the TR3 cart for mobile storage. The locker has 20 drawer slots in 3" increments and holds interchangeable components.

Notes

Order rail components separately:

- Rail for roller rail (CO345)
- Roller for roller rail (CO34731)
- Standard rail assembly (CO282)

Order transport/storage components separately:


- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)
- TR3 cart (CO342)

Order external components separately:

- Chest tube rack (CO466)
- Fogarty catheter rack (CO467)

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Locker, with Tambour Door *continued*

Co/Struc® Storage

Specification Information

Step 1.

C0561

Step 2. Lock

FF	no lock
FL	keyed differently
KA	keyed alike

Prices for Steps 1-2.

C0561 FF	\$1358
FL	\$1452
KA	\$1452

Step 3. Locker Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 4. Door Finish

Inner tone light (HF) or light tone (LT) door finish must match locker finish.

DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
TC	teal medium	+\$0

Step 5. Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0


20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0
43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0
50	key number 50	+\$0
51	key number 51	+\$0
52	key number 52	+\$0
53	key number 53	+\$0
54	key number 54	+\$0
55	key number 55	+\$0
56	key number 56	+\$0
57	key number 57	+\$0
58	key number 58	+\$0
59	key number 59	+\$0
60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0

Locker, with Tambour Door *continued*

67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0
90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0

Locker

CO562



Product Information

Description

This storage unit hangs on rails for stationary storage, or it is used with the TR3 cart for mobile storage. The locker has 20 drawer slots in 3" increments and holds interchangeable components.

Notes

Order rail components separately:

- Rail for roller rail (CO345)
- Roller for roller rail (CO34731)
- Standard rail assembly (CO282)

Order transport/storage components separately:

- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Drawer (CO207, CO208, or CO209)
- Pullout catheter rack (CO469.)
- TR3 cart (CO342)


Order external components separately:

- Chest tube rack (CO466)
- Fogarty catheter rack (CO467)


Order optional locker tambour door (CO219) separately.

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



TOP



SIDE

FRONT

Specification Information

Step 1.

C0562

Step 2. Lock

FF	no lock
FL	keyed differently
KA	keyed alike

Prices for Steps 1-2.

C0562 FF	\$1236
FL	\$1329
KA	\$1329

Step 3. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 4. Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0


28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0
43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0
50	key number 50	+\$0
51	key number 51	+\$0
52	key number 52	+\$0
53	key number 53	+\$0
54	key number 54	+\$0
55	key number 55	+\$0
56	key number 56	+\$0
57	key number 57	+\$0
58	key number 58	+\$0
59	key number 59	+\$0
60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0
67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0

Locker *continued*

Co/Struc® Storage	75	key number 75	+\$0
	76	key number 76	+\$0
	77	key number 77	+\$0
	78	key number 78	+\$0
	79	key number 79	+\$0
	80	key number 80	+\$0
	81	key number 81	+\$0
	82	key number 82	+\$0
	83	key number 83	+\$0
	84	key number 84	+\$0
	85	key number 85	+\$0
	86	key number 86	+\$0
	87	key number 87	+\$0
	88	key number 88	+\$0
	89	key number 89	+\$0
	90	key number 90	+\$0
	91	key number 91	+\$0
	92	key number 92	+\$0
	93	key number 93	+\$0
	94	key number 94	+\$0
	95	key number 95	+\$0
	96	key number 96	+\$0
	97	key number 97	+\$0
	98	key number 98	+\$0
	99	key number 99	+\$0

Locker on Wheel Base, with Tambour Door

CO564



Product Information

Description

This enclosed, mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters. The casters include 3 swivel and 1 swivel with brake. The tambour door slides vertically for access to locker materials. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.

Notes

Order internal storage components separately:


- A-, B-, or C-size drawers (CO207, CO208, or CO209)
- C tray/shelf (CO203)
- C wire shelf (CO252)
- Cassette assembly (CO565)
- Pullout catheter rack (CO469.)

Order external components separately:

- Chest tube rack (CO466)
- Fogarty catheter rack (CO467)

For general information on locks, see Keyed-Alike Information in Appendices.

Dimensions



Locker on Wheel Base, with Tambour Door *continued*

Co/Struc® Storage

Specification Information

Step 1.

C0564

Step 2. Lock

FF	no lock
FL	keyed differently
KA	keyed alike

Prices for Steps 1-2.

C0564 FF	\$1950
FL	\$2039
KA	\$2039

Step 3. Locker Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 4. Door Finish

Inner tone light (HF) or light tone (LT) door finish must match locker finish.

DA	deep blue	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
TC	teal medium	+\$0

Step 5. Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0


20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0
43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0
50	key number 50	+\$0
51	key number 51	+\$0
52	key number 52	+\$0
53	key number 53	+\$0
54	key number 54	+\$0
55	key number 55	+\$0
56	key number 56	+\$0
57	key number 57	+\$0
58	key number 58	+\$0
59	key number 59	+\$0
60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0

Locker on Wheel Base, with
Tambour Door *continued*

67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0
90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0

Locker on Wheel Base


C0563



Product Information
Description
This mobile storage unit holds interchangeable components to store and transport materials. The locker has 20 drawer slots in 3" increments and a counterweighted base with 5"-diameter casters. The casters include 3 swivel and 1 swivel with brake. The locker does not include a handle; it is pushed or pulled by grasping the sides. It cannot be linked for multiple towing and cannot be transported in a truck.
Notes
Order internal storage components separately: <ul style="list-style-type: none">• A-, B-, or C-size drawers (CO207, CO208, or CO209)• C tray/shelf (CO203)• C wire shelf (CO252)• Cassette assembly (CO565)• Pullout catheter rack (CO469.) Order external components separately: <ul style="list-style-type: none">• Chest tube rack (CO466)• Fogarty catheter rack (CO467) Order optional locker tambour door (CO219) and locker lock kit (CO277) separately.
Dimensions


Specification Information		
Step 1.		
C0563FF		\$1834
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage



Locker Tambour Door

C0219




Product Information


Description

This vertical door provides closure on a locker. It opens by sliding into the locker bottom.

Dimensions



SIDE



FRONT

Specification Information

Step 1.

C0219FF \$130


Step 2. Surface Finish

DA	deep blue	+\$0
HF	inner tone light	+\$0
TC	teal medium	+\$0
LT	light tone	+\$0

Co/Struc® Storage


Locker Lock Kit

CO277



Product Information	
Description	This kit includes a lock and converts a locker with tambour door to a lockable storage unit.
Notes	For general information on locks, see Keyed-Alike Information in Appendices. Kit must be field installed.
Dimensions	

Specification Information		
Step 1.		
CO277		
Step 2. Lock		
FF	keyed differently	
KA	keyed alike	
Prices for Steps 1-2.		
CO277	FF	\$108
	KA	\$108
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
Step 4. Key Number		
For keyed alike (KA)		
01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0



Co/Struc® Storage


Locker Lock Kit *continued*

30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			
74	key number 74	+\$0			
75	key number 75	+\$0			
76	key number 76	+\$0			

Co/Struc® Storage

Cabinet

CV300.



Product Information


Description
This 92"-high, enclosed cabinet attaches to an architectural wall and must stand on the floor. It has a tambour door that slides vertically for access to cabinet materials. The cabinet is available with or without a lock and holds components that support endoscopes, CS tray/shelves, and CS wire shelves. Attachment hardware is included.

Notes


- Order interior components separately:
- Cabinet drip pan (CV313.)
 - Cabinet floor (CV314.)
 - Cabinet interior adapter (CV312.)
 - Cabinet pullout catheter rack (CV316.)
 - CS tray/shelf (CO204)
 - CS wire shelf (CO253)
 - Distal tube bracket (CV311.)
 - Endoscope bracket (CV310.)

For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions



TOP



SIDE



FRONT

Co/Struc® Storage

Specification Information

Step 1.

CV300.9214

Step 2. Lock

*For no lock, skip this step.***L** keyed differently**K** keyed alike

Prices for Steps 1-2.

CV300.9214 \$1718**CV300.9214 L** \$1771**CV300.9214 K** \$1771

Step 3. Surface/Door Finish

HF inner tone light +\$0**LT** light tone +\$0

Step 4. Key Number


*For keyed alike (K)***01** key number 01 +\$0**02** key number 02 +\$0**03** key number 03 +\$0**04** key number 04 +\$0**05** key number 05 +\$0**06** key number 06 +\$0**07** key number 07 +\$0**08** key number 08 +\$0**09** key number 09 +\$0**10** key number 10 +\$0**11** key number 11 +\$0**12** key number 12 +\$0**13** key number 13 +\$0**14** key number 14 +\$0**15** key number 15 +\$0**16** key number 16 +\$0**17** key number 17 +\$0**18** key number 18 +\$0**19** key number 19 +\$0**20** key number 20 +\$0**21** key number 21 +\$0**22** key number 22 +\$0**23** key number 23 +\$0**24** key number 24 +\$0**25** key number 25 +\$0**26** key number 26 +\$0**27** key number 27 +\$0**28** key number 28 +\$0**29** key number 29 +\$0**30** key number 30 +\$0**31** key number 31 +\$0**32** key number 32 +\$0**33** key number 33 +\$0**34** key number 34 +\$0**35** key number 35 +\$0**36** key number 36 +\$0**37** key number 37 +\$0**38** key number 38 +\$0**39** key number 39 +\$0**40** key number 40 +\$0**41** key number 41 +\$0**42** key number 42 +\$0**43** key number 43 +\$0**44** key number 44 +\$0**45** key number 45 +\$0**46** key number 46 +\$0**47** key number 47 +\$0**48** key number 48 +\$0**49** key number 49 +\$0**50** key number 50 +\$0**51** key number 51 +\$0**52** key number 52 +\$0**53** key number 53 +\$0**54** key number 54 +\$0**55** key number 55 +\$0**56** key number 56 +\$0**57** key number 57 +\$0**58** key number 58 +\$0**59** key number 59 +\$0**60** key number 60 +\$0**61** key number 61 +\$0**62** key number 62 +\$0**63** key number 63 +\$0**64** key number 64 +\$0**65** key number 65 +\$0**66** key number 66 +\$0**67** key number 67 +\$0**68** key number 68 +\$0**69** key number 69 +\$0**70** key number 70 +\$0**71** key number 71 +\$0**72** key number 72 +\$0**73** key number 73 +\$0**74** key number 74 +\$0

Cabinet *continued*


Co/Struc® Storage	75	key number 75	+\$0
	76	key number 76	+\$0
	77	key number 77	+\$0
	78	key number 78	+\$0
	79	key number 79	+\$0
	80	key number 80	+\$0
	81	key number 81	+\$0
	82	key number 82	+\$0
	83	key number 83	+\$0
	84	key number 84	+\$0
	85	key number 85	+\$0
	86	key number 86	+\$0
	87	key number 87	+\$0
	88	key number 88	+\$0
	89	key number 89	+\$0
	90	key number 90	+\$0
	91	key number 91	+\$0
	92	key number 92	+\$0
	93	key number 93	+\$0
	94	key number 94	+\$0
	95	key number 95	+\$0
	96	key number 96	+\$0
	97	key number 97	+\$0
	98	key number 98	+\$0
	99	key number 99	+\$0

Cabinet Lock Kit

CV315.



Product Information
Description
This kit converts a cabinet without lock to a lockable cabinet.
Notes
Order cabinet (CV300.) separately.
For general information on locks, see Keyed-Alike Information in Appendices.
Dimensions



Specification Information
Step 1.
CV315.
Step 2. Lock
L keyed differently
K keyed alike
Prices for Steps 1-2.
CV315. L \$115
K \$115
Step 3. Surface Finish
HF inner tone light +\$0
LT light tone +\$0
Step 4. Key Number
For keyed alike (K)
01 key number 01 +\$0
02 key number 02 +\$0
03 key number 03 +\$0
04 key number 04 +\$0
05 key number 05 +\$0
06 key number 06 +\$0
07 key number 07 +\$0
08 key number 08 +\$0
09 key number 09 +\$0
10 key number 10 +\$0
11 key number 11 +\$0
12 key number 12 +\$0
13 key number 13 +\$0
14 key number 14 +\$0
15 key number 15 +\$0
16 key number 16 +\$0
17 key number 17 +\$0
18 key number 18 +\$0
19 key number 19 +\$0
20 key number 20 +\$0
21 key number 21 +\$0
22 key number 22 +\$0
23 key number 23 +\$0
24 key number 24 +\$0
25 key number 25 +\$0
26 key number 26 +\$0
27 key number 27 +\$0
28 key number 28 +\$0
29 key number 29 +\$0

Co/Struc® Storage


Cabinet Lock Kit *continued*

Co/Struc® Storage

30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			
74	key number 74	+\$0			
75	key number 75	+\$0			
76	key number 76	+\$0			

Cabinet Pullout Catheter Rack


CV316.




Product Information
Description
This 22"-wide, powder-coated rack mounts inside a cabinet to organize catheters. It has 3 independent, pullout slides. Each pullout slide has 2 prongs, a side label clip for each prong, and 1 front label holder. Each prong holds up to 14 catheters. Finish is soft white. Mounting hardware is included.
Notes
Order cabinet (CV300.) separately.
Dimensions

Specification Information
Step 1.
CV316.
\$342

Co/Struc® Storage




Endoscope BracketCV310.



Product Information
Description
This bracket mounts at any height in a cabinet. It has 5 slots and adjustable cradles for vertical storage of endoscopes. Each slot is 3¾" deep. Mounting hardware is included.
Notes
Order following products separately: <ul style="list-style-type: none">• Cabinet (CV300.)• Cabinet drip pan (CV313.)• Distal tube bracket (CV311.)
Recommended installation height is 5¼" below top of cabinet back panel.
Dimensions


Specification Information		
Step 1.		
CV310.		\$289
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage



Distal Tube Bracket

CV311.



Product Information

Description
This bracket mounts at any height in a cabinet. It has 5 slots for vertical storage of an endoscope distal tube. Each slot is 7/8" deep. Mounting hardware is included.


Notes

Order following products separately:

- Cabinet (CV300.)
- Cabinet drip pan (CV313.)
- Endoscope bracket (CV310.)

Recommended installation height is 44 1/4" below top of cabinet back panel.

Dimensions



Specification Information

Step 1.

CV311. \$134

Step 2. Surface Finish


HF inner tone light +\$0

LT light tone +\$0

Co/Struc® Storage


Cabinet Interior Adapter

CV312.



Product Information	
Description	
This adapter mounts at any height in a cabinet. It has 5 slots at 3" intervals to support a CS tray/shelf or CS wire shelf. The adapter cannot hold drawers. Package contains 2 non-handed adapters. Mounting hardware is included.	
Notes	
Order cabinet (CV300.) and cabinet floor (CV314.) separately.	
Order storage components separately:	
• CS tray/shelf (CO204)	
• CS wire shelf (CO253)	
Dimensions	


Specification Information		
Step 1.		
CV312.05		\$89
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Co/Struc® Storage

Cabinet Drip Pan


CV313.




Product Information
Description
This plastic pan mounts at any height in a cabinet. It has a concave bottom to collect moisture from washed endoscopes. Attachment hardware is included.
Notes
Order cabinet (CV300.) separately.
Order endoscope support components separately:
<ul style="list-style-type: none">• Distal tube bracket (CV311.)• Endoscope bracket (CV310.)
Dimensions

Specification Information
Step 1.
CV313.
\$89
Step 2. Surface Finish
HF inner tone light
LT light tone
+\$0
+\$0

Co/Struc® Storage







Product Information
Description
This metal floor encloses the bottom of a cabinet when a cabinet interior adapter is installed. Attachment hardware is included.
Notes
Order cabinet (CV300.) and cabinet interior adapter (CV312.) separately.
Dimensions

Specification Information		
Step 1.		
CV314.		\$46
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0


Co/Struc® Storage



TOP




SIDE



FRONT

CST Frame

CO206




Product Information


Description
This frame hangs on a standard rail or an adapter rail to support and enclose CS components and accessories. It has 6 slots at 3" intervals. The frame cannot accept drawers or mount under work surfaces.

Notes
Order CS tray/shelf (CO2o4) separately.
For mobile storage, order L cart (CO212) separately.
To enclose frame, order CS lid (CO2o2) and CST flipper (CO214) separately.


Dimensions



TOP



SIDE




FRONT

Specification Information

Step 1.		
CO206FF		\$196
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage



Product Information

Description

This lid provides a top surface for a CST frame. It has raised edges to keep items from sliding off.

Dimensions

Specification Information

Step 1.


CO202FF \$105

Step 2. Surface Finish


HF inner tone light +\$0

LT light tone +\$0

Co/Struc® Storage




TOP



FRONT

CST Flipper


CO214



Product Information

Description
This component provides front closure for a CST frame. The flipper swings outward, then recedes into the frame.


Dimensions



Specification Information

Step 1.		
CO214FF		\$65
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0


Co/Struc® Storage




Product Information	
Description	
This component provides storage within a CST frame or cabinet interior adapter. It has 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.	
Notes	
Order CST frame (CO2o6) or cabinet interior adapter (CV312.) separately.	
Tray/shelf surface with raised edges can be configured with subcontainers; order subcontainers separately.	
Dimensions	

Specification Information		
Step 1.		
CO204FF		\$65
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0


Co/Struc® Storage



TOP




FRONT



CS Wire Shelf

C0253




Product Information


Description
This wire shelf fits into the slots of a CST frame or cabinet interior adapter. It stores items requiring air circulation. The shelf cannot be used with a subcontainer, subdivider vane, or subdivider.

Notes
Order CST frame (C0206) or cabinet interior adapter (CV312.) separately.

Dimensions



TOP



FRONT

Specification Information


Step 1.
C0253FF \$103

Step 2. Surface Finish
LU soft white +\$0

Co/Struc® Storage

C Frame

CO205






Product Information

Description
This frame hangs on a standard rail or adapter rail, or it is supported by an L cart, process table, counter top, or heavy-duty storage work surface. It has 4 slots at 3" intervals to enclose and support storage components and accessories. The C frame cannot mount under a 24"-wide work surface and it is not recommended for overhead storage.

Notes

- Order support components separately:
- Adapter rail (C1610.)
 - Counter top (CO217, CO555, or CO556)
 - L cart (CO212)
 - Process table (C4115. or C4215.)
 - Standard rail (CO282)
- Order storage components separately:
- C tray/shelf (CO203)
 - Drawer (CO207, CO208, or CO209)

Dimensions




Specification Information

Step 1.		
CO205FF		\$173
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage

Counter Top

C0217
C0555
C0556



Product Information


Description

This surface mounts on a C frame to provide a work surface or storage surface. It includes a gripper package for attaching the top to a C frame. The single or extended counter top fits on 1 C frame; the double counter top fits on 2 C frames. The extended counter top extends 4" past each C frame side.


Notes

Order C frame (C0205) separately.
Double counter top is not recommended for use on L cart (C0212).


Dimensions




TOP



SINGLE
FRONT




EXTENDED
FRONT



DOUBLE
FRONT

Co/Struc® Storage


Specification Information		
Step 1.		
CO		
Step 2. Configuration/Width		
555FF	single	
556FF	extended	
217FF	double	
Prices for Steps 1-2.		
CO555FF		\$122
CO556FF		\$151
CO217FF		\$262
Step 3. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Product Information
Description
This lid provides a top for a C frame. It has raised side edges to keep items from sliding off.
Notes
Use lid when storing L cart below 38"-high process table or cantilevered work surface; existing drawer bearers must be removed.
Dimensions


Specification Information			
Step 1.			
CO201FF			\$122
Step 2. Surface Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0

Co/Struc® Storage



C Flipper

CO213



Product Information

Description
This component provides front closure for a C frame. The flipper swings outward, then recedes into the frame.

Dimensions


Specification Information


Step 1.
CO213FF \$50

Step 2. Surface Finish

HF inner tone light +\$0
LT light tone +\$0

Co/Struc® Storage





Product Information

Description
This component provides storage within a C frame, locker, or interior cabinet adapter. It is designed with 2 usable surfaces: 1 surface has raised edges to prevent items from rolling off; the reverse surface is flat for applications where raised edges would interfere with material storage.


Notes

Order support component separately:


- C frame (CO205)
- Interior cabinet adapter (CG590.)
- Locker or locker on wheel base (CO561, CO562, CO563, or CO564)

Tray/shelf surface with raised edges can be configured with subcontainers or subdivider vane. Order subdivider vane (CO242) or subcontainers separately.

Dimensions



TOP



FRONT


Specification Information

Step 1.		
CO203FF		\$82
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage

Subdivider Vane

C0242



Product Information

Description
This component organizes IV bags, forms, tape rolls, and other small items within a drawer or C tray/shelf.


Dimensions


Specification Information

Step 1.
C0242FF \$28

Step 2. Surface Finish
LU soft white +\$0

Co/Struc® Storage







Product Information
Description
This wire shelf fits into the slots of a C frame, locker, or interior cabinet adapter and stores items requiring air circulation. It cannot be used with a subcontainer, subdivider vane, or subdivider.
Notes
Order support component separately:
<ul style="list-style-type: none">• C frame (CO205)• Interior cabinet adapter (CG590.)• Locker or locker on wheel base (CO561, CO562, CO563, or CO564)
Dimensions

Specification Information		
Step 1.		
CO252FF		\$103
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage







Product Information
Description
This hardware package mounts a C frame under a work surface. 1 gripper package is required to mount each C frame to a work surface. Attachment screws are included.
Notes
C frame cannot mount with grippers under heavy-duty work surface.
Dimensions

Specification Information		
Step 1.		
CO266FF		\$47
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Co/Struc® Storage




TOP



FRONT

Drawer Caddy with Casters

C0692.
C0693.



Product Information

Description

This mobile caddy fits below a hanging or freestanding work surface. It has full-width pulls, casters, 2 3"-high drawers, 1 6"-high drawer, and 1 9"-high drawer. The drawers can be removed and reconfigured within the caddy or used in lockers, C frames, or Series 20 carts.


The 9"-high drawer is available with a lock.

Notes


For caddy without drawers, order drawers separately.

For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions



TOP



SIDE



FRONT

Specification Information

Step 1.

C069

Step 2. Configuration

2.20 no drawers

3.20 drawers

Step 3. Lock

For no lock, skip this step.

For no drawers (2.20), skip this step.

For drawers (3.20)

L 9"-high drawer keyed differently

K 9"-high drawer keyed alike

Prices for Steps 1-3.

C0692.20 \$740

C0693.20 \$862

C0693.20 L \$950

C0693.20 K \$950

Step 4. Surface Finish

HF inner tone light +\$0

LT light tone +\$0

Step 5. Drawer Finish

For drawers (3.20)

B8 green tone +\$0

B9 green tone dark +\$0

DA deep blue +\$0

HF inner tone light +\$0

L3 lavender light +\$0

L4 lavender medium +\$0

L5 lavender dark +\$0

LT light tone +\$0

LU soft white +\$0

TC teal medium +\$0

VM deep magenta +\$0

VR vivid red +\$0

W5 bordeaux +\$0

YT yellow +\$0

Drawer Caddy with Casters *continued*

Co/Struc® Storage

Step 6. Key Number

For 9"-high drawer keyed alike (K)


01	key number 01	+\$0	46	key number 46	+\$0
02	key number 02	+\$0	47	key number 47	+\$0
03	key number 03	+\$0	48	key number 48	+\$0
04	key number 04	+\$0	49	key number 49	+\$0
05	key number 05	+\$0	50	key number 50	+\$0
06	key number 06	+\$0	51	key number 51	+\$0
07	key number 07	+\$0	52	key number 52	+\$0
08	key number 08	+\$0	53	key number 53	+\$0
09	key number 09	+\$0	54	key number 54	+\$0
10	key number 10	+\$0	55	key number 55	+\$0
11	key number 11	+\$0	56	key number 56	+\$0
12	key number 12	+\$0	57	key number 57	+\$0
13	key number 13	+\$0	58	key number 58	+\$0
14	key number 14	+\$0	59	key number 59	+\$0
15	key number 15	+\$0	60	key number 60	+\$0
16	key number 16	+\$0	61	key number 61	+\$0
17	key number 17	+\$0	62	key number 62	+\$0
18	key number 18	+\$0	63	key number 63	+\$0
19	key number 19	+\$0	64	key number 64	+\$0
20	key number 20	+\$0	65	key number 65	+\$0
21	key number 21	+\$0	66	key number 66	+\$0
22	key number 22	+\$0	67	key number 67	+\$0
23	key number 23	+\$0	68	key number 68	+\$0
24	key number 24	+\$0	69	key number 69	+\$0
25	key number 25	+\$0	70	key number 70	+\$0
26	key number 26	+\$0	71	key number 71	+\$0
27	key number 27	+\$0	72	key number 72	+\$0
28	key number 28	+\$0	73	key number 73	+\$0
29	key number 29	+\$0	74	key number 74	+\$0
30	key number 30	+\$0	75	key number 75	+\$0
31	key number 31	+\$0	76	key number 76	+\$0
32	key number 32	+\$0	77	key number 77	+\$0
33	key number 33	+\$0	78	key number 78	+\$0
34	key number 34	+\$0	79	key number 79	+\$0
35	key number 35	+\$0	80	key number 80	+\$0
36	key number 36	+\$0	81	key number 81	+\$0
37	key number 37	+\$0	82	key number 82	+\$0
38	key number 38	+\$0	83	key number 83	+\$0
39	key number 39	+\$0	84	key number 84	+\$0
40	key number 40	+\$0	85	key number 85	+\$0
41	key number 41	+\$0	86	key number 86	+\$0
42	key number 42	+\$0	87	key number 87	+\$0
43	key number 43	+\$0	88	key number 88	+\$0
44	key number 44	+\$0	89	key number 89	+\$0
45	key number 45	+\$0	90	key number 90	+\$0
			91	key number 91	+\$0
			92	key number 92	+\$0

Drawer Caddy with Casters *continued*

93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0

Drawer

CO207
CO208
CO209




Product Information

Description
This drawer is used in a locker or C frame; it does not need any other support components. The drawer can also mount under a work surface. It is available in 3 heights and is available with or without a lock.


Notes

To mount drawer under table or work surface, order drawer bearer package (CO228) or lockable drawer bearer (CO396) separately.
Order optional A-size subcontainer (CO231, CO232, CO233, CO234, CO235, or CO369) or B-size subcontainer (CO236, CO238, or CO240) separately.
For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions




TOP



A SIZE
FRONT



B SIZE
FRONT



C SIZE
FRONT

Co/Struc® Storage

Specification Information		
Step 1.		
CO20		
Step 2. Size/Height		
7	a size	
8	b size	
9	c size	
Step 3. Lock		
FF	no lock	
FL	keyed differently	
KA	keyed alike	
Prices for Steps 1-3.		
CO207	FF	\$76
	FL	\$154
	KA	\$154
CO208	FF	\$98
	FL	\$175
	KA	\$175
CO209	FF	\$125
	FL	\$207
	KA	\$207
Step 4. Surface Finish		
B8	green tone	+\$0
B9	green tone dark	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
L3	lavender light	+\$0
L4	lavender medium	+\$0
L5	lavender dark	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TC	teal medium	+\$0
VM	deep magenta	+\$0
VR	vivid red	+\$0
W5	bordeaux	+\$0
YT	yellow	+\$0
Step 5. Key Number		
<i>For keyed alike (KA)</i>		
01	key number 01	+\$0
02	key number 02	+\$0

03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0
43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0

Drawer *continued*

Co/Struc® Storage


50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			
54	key number 54	+\$0			
55	key number 55	+\$0			
56	key number 56	+\$0			
57	key number 57	+\$0			
58	key number 58	+\$0			
59	key number 59	+\$0			
60	key number 60	+\$0			
61	key number 61	+\$0			
62	key number 62	+\$0			
63	key number 63	+\$0			
64	key number 64	+\$0			
65	key number 65	+\$0			
66	key number 66	+\$0			
67	key number 67	+\$0			
68	key number 68	+\$0			
69	key number 69	+\$0			
70	key number 70	+\$0			
71	key number 71	+\$0			
72	key number 72	+\$0			
73	key number 73	+\$0			
74	key number 74	+\$0			
75	key number 75	+\$0			
76	key number 76	+\$0			
77	key number 77	+\$0			
78	key number 78	+\$0			
79	key number 79	+\$0			
80	key number 80	+\$0			
81	key number 81	+\$0			
82	key number 82	+\$0			
83	key number 83	+\$0			
84	key number 84	+\$0			
85	key number 85	+\$0			
86	key number 86	+\$0			
87	key number 87	+\$0			
88	key number 88	+\$0			
89	key number 89	+\$0			
90	key number 90	+\$0			
91	key number 91	+\$0			
92	key number 92	+\$0			
93	key number 93	+\$0			
94	key number 94	+\$0			
95	key number 95	+\$0			
96	key number 96	+\$0			

Drawer with Lockable Lid

C0370

C0372

C0374



Product Information

Description


This drawer has a hinged steel lid that can be locked. The lid can open when the drawer is pulled out to a stop position; the lid cannot be removed.

The drawer is 3", 6", or 9" high.


Notes

For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions




TOP



A SIZE
FRONT



B SIZE
FRONT



C SIZE
FRONT

Specification Information

Step 1.

C037

Step 2. Size/Height

0	a size
2	b size
4	c size

Step 3. Lid Lock

FF	keyed differently
KA	keyed alike

Prices for Steps 1-3.

C0370	FF	\$227
	KA	\$227
C0372	FF	\$235
	KA	\$235
C0374	FF	\$248
	KA	\$248

Step 4. Drawer Finish

B8	green tone	+\$0
B9	green tone dark	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
L3	lavender light	+\$0
L4	lavender medium	+\$0
L5	lavender dark	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TC	teal medium	+\$0
VM	deep magenta	+\$0
VR	vivid red	+\$0
W5	bordeaux	+\$0
YT	yellow	+\$0

Step 5. Lid Finish

LU	soft white	+\$0
----	------------	------

Step 6. Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0

Drawer with Lockable Lid *continued*

Co/Struc® Storage

04	key number 04	+\$0	51	key number 51	+\$0
05	key number 05	+\$0	52	key number 52	+\$0
06	key number 06	+\$0	53	key number 53	+\$0
07	key number 07	+\$0	54	key number 54	+\$0
08	key number 08	+\$0	55	key number 55	+\$0
09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0

Drawer with Lockable Lid *continued*


98	key number 98	+\$0
99	key number 99	+\$0

Locked Drawer with Lid

C0549

C0550

C0551



Product Information


Description

This 3"-, 6"-, or 9"-high drawer has a lid and a drawer lock. The lid can open when the drawer is pulled out to a stop position; the lid cannot be removed.


Notes

For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions




TOP



A SIZE
FRONT



B SIZE
FRONT



C SIZE
FRONT

Specification Information

Step 1.

C05

Step 2. Size/Height

49 a size

50 b size

51 c size

Step 3. Drawer Lock

FL keyed differently

KA keyed alike

Prices for Steps 1-3.

C0549 FL	\$237
KA	\$237

C0550 FL	\$245
KA	\$245

C0551 FL	\$259
KA	\$259

Step 4. Drawer Finish

B8	green tone	+\$0
B9	green tone dark	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
L3	lavender light	+\$0
L4	lavender medium	+\$0
L5	lavender dark	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TC	teal medium	+\$0
VM	deep magenta	+\$0
VR	vivid red	+\$0
W5	bordeaux	+\$0
YT	yellow	+\$0

Step 5. Drawer Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0

Locked Drawer with Lid *continued*

07	key number 07	+\$0	54	key number 54	+\$0
08	key number 08	+\$0	55	key number 55	+\$0
09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			

Co/Struc® Storage

Locked Drawer with Lid *continued*

Step 6. Lid Finish		
LU	soft white	+\$0


Co/Struc® Storage

Locked Drawer with Lockable Lid

C0371

C0373

C0375



Product Information

Description


This drawer has a lockable lid and a drawer lock. The lid can open when the drawer is pulled out to a stop position; the lid cannot be removed.

The drawer is 3", 6", or 9" high.


Notes

For general information on locks, see Keyed-Alike Information in Appendices.


Dimensions




TOP



A SIZE
FRONT



B SIZE
FRONT



C SIZE
FRONT

Specification Information

Step 1.

C037

Step 2. Size/Height

- | | |
|---|--------|
| 1 | a size |
| 3 | b size |
| 5 | c size |

Step 3. Drawer/Lid Lock

- | | |
|----|-------------------|
| FF | keyed differently |
| KA | keyed alike |

Prices for Steps 1-3.

C0371	FF	\$257
	KA	\$257
C0373	FF	\$266
	KA	\$266
C0375	FF	\$281
	KA	\$281

Step 4. Drawer Finish

B8	green tone	+\$0
B9	green tone dark	+\$0
DA	deep blue	+\$0
HF	inner tone light	+\$0
L3	lavender light	+\$0
L4	lavender medium	+\$0
L5	lavender dark	+\$0
LT	light tone	+\$0
LU	soft white	+\$0
TC	teal medium	+\$0
VM	deep magenta	+\$0
VR	vivid red	+\$0
W5	bordeaux	+\$0
YT	yellow	+\$0

Step 5. Drawer Key Number

For keyed alike (KA)

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0

Locked Drawer with Lockable Lid

continued

Co/Struc® Storage

07	key number 07	+\$0	54	key number 54	+\$0
08	key number 08	+\$0	55	key number 55	+\$0
09	key number 09	+\$0	56	key number 56	+\$0
10	key number 10	+\$0	57	key number 57	+\$0
11	key number 11	+\$0	58	key number 58	+\$0
12	key number 12	+\$0	59	key number 59	+\$0
13	key number 13	+\$0	60	key number 60	+\$0
14	key number 14	+\$0	61	key number 61	+\$0
15	key number 15	+\$0	62	key number 62	+\$0
16	key number 16	+\$0	63	key number 63	+\$0
17	key number 17	+\$0	64	key number 64	+\$0
18	key number 18	+\$0	65	key number 65	+\$0
19	key number 19	+\$0	66	key number 66	+\$0
20	key number 20	+\$0	67	key number 67	+\$0
21	key number 21	+\$0	68	key number 68	+\$0
22	key number 22	+\$0	69	key number 69	+\$0
23	key number 23	+\$0	70	key number 70	+\$0
24	key number 24	+\$0	71	key number 71	+\$0
25	key number 25	+\$0	72	key number 72	+\$0
26	key number 26	+\$0	73	key number 73	+\$0
27	key number 27	+\$0	74	key number 74	+\$0
28	key number 28	+\$0	75	key number 75	+\$0
29	key number 29	+\$0	76	key number 76	+\$0
30	key number 30	+\$0	77	key number 77	+\$0
31	key number 31	+\$0	78	key number 78	+\$0
32	key number 32	+\$0	79	key number 79	+\$0
33	key number 33	+\$0	80	key number 80	+\$0
34	key number 34	+\$0	81	key number 81	+\$0
35	key number 35	+\$0	82	key number 82	+\$0
36	key number 36	+\$0	83	key number 83	+\$0
37	key number 37	+\$0	84	key number 84	+\$0
38	key number 38	+\$0	85	key number 85	+\$0
39	key number 39	+\$0	86	key number 86	+\$0
40	key number 40	+\$0	87	key number 87	+\$0
41	key number 41	+\$0	88	key number 88	+\$0
42	key number 42	+\$0	89	key number 89	+\$0
43	key number 43	+\$0	90	key number 90	+\$0
44	key number 44	+\$0	91	key number 91	+\$0
45	key number 45	+\$0	92	key number 92	+\$0
46	key number 46	+\$0	93	key number 93	+\$0
47	key number 47	+\$0	94	key number 94	+\$0
48	key number 48	+\$0	95	key number 95	+\$0
49	key number 49	+\$0	96	key number 96	+\$0
50	key number 50	+\$0	97	key number 97	+\$0
51	key number 51	+\$0	98	key number 98	+\$0
52	key number 52	+\$0	99	key number 99	+\$0
53	key number 53	+\$0			

Locked Drawer with Lockable Lid

continued

Step 6. Lid Finish		
LU	soft white	+\$0

Step 7. Lid Key Number		
For keyed alike (KA)		

01	key number 01	+\$0
02	key number 02	+\$0
03	key number 03	+\$0
04	key number 04	+\$0
05	key number 05	+\$0
06	key number 06	+\$0
07	key number 07	+\$0
08	key number 08	+\$0
09	key number 09	+\$0
10	key number 10	+\$0
11	key number 11	+\$0
12	key number 12	+\$0
13	key number 13	+\$0
14	key number 14	+\$0
15	key number 15	+\$0
16	key number 16	+\$0
17	key number 17	+\$0
18	key number 18	+\$0
19	key number 19	+\$0
20	key number 20	+\$0
21	key number 21	+\$0
22	key number 22	+\$0
23	key number 23	+\$0
24	key number 24	+\$0
25	key number 25	+\$0
26	key number 26	+\$0
27	key number 27	+\$0
28	key number 28	+\$0
29	key number 29	+\$0
30	key number 30	+\$0
31	key number 31	+\$0
32	key number 32	+\$0
33	key number 33	+\$0
34	key number 34	+\$0
35	key number 35	+\$0
36	key number 36	+\$0
37	key number 37	+\$0
38	key number 38	+\$0
39	key number 39	+\$0
40	key number 40	+\$0
41	key number 41	+\$0
42	key number 42	+\$0

43	key number 43	+\$0
44	key number 44	+\$0
45	key number 45	+\$0
46	key number 46	+\$0
47	key number 47	+\$0
48	key number 48	+\$0
49	key number 49	+\$0
50	key number 50	+\$0
51	key number 51	+\$0
52	key number 52	+\$0
53	key number 53	+\$0
54	key number 54	+\$0
55	key number 55	+\$0
56	key number 56	+\$0
57	key number 57	+\$0
58	key number 58	+\$0
59	key number 59	+\$0
60	key number 60	+\$0
61	key number 61	+\$0
62	key number 62	+\$0
63	key number 63	+\$0
64	key number 64	+\$0
65	key number 65	+\$0
66	key number 66	+\$0
67	key number 67	+\$0
68	key number 68	+\$0
69	key number 69	+\$0
70	key number 70	+\$0
71	key number 71	+\$0
72	key number 72	+\$0
73	key number 73	+\$0
74	key number 74	+\$0
75	key number 75	+\$0
76	key number 76	+\$0
77	key number 77	+\$0
78	key number 78	+\$0
79	key number 79	+\$0
80	key number 80	+\$0
81	key number 81	+\$0
82	key number 82	+\$0
83	key number 83	+\$0
84	key number 84	+\$0
85	key number 85	+\$0
86	key number 86	+\$0
87	key number 87	+\$0
88	key number 88	+\$0
89	key number 89	+\$0

Co/Struc® Storage

Locked Drawer with Lockable Lid


continued

Co/Struc® Storage

90	key number 90	+\$0
91	key number 91	+\$0
92	key number 92	+\$0
93	key number 93	+\$0
94	key number 94	+\$0
95	key number 95	+\$0
96	key number 96	+\$0
97	key number 97	+\$0
98	key number 98	+\$0
99	key number 99	+\$0


Clear Lid


C0508.



Product Information
Description
This clear lid fits on a Co/Struc drawer and provides top closure.
Notes
Order drawer (C0207, C0208, or C0209) separately.
Dimensions

Specification Information
Step 1.
C0508.20
\$66





Product Information


Description
This package is used to mount a nonlocking drawer under a work surface. The drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.

Dimensions


Specification Information

Step 1.		
C0228FF		\$25
Step 2. Surface Finish		
BK	black	+\$0


Co/Struc® Storage



TOP




FRONT



Lockable Drawer Bearer


C0396



Product Information

Description
This component is used to mount a drawer with lock under a work surface. The lockable drawer bearer has integral stops to prevent drawers from falling when pulled out. Mounting hardware is included.

Dimensions




Specification Information

Step 1.	
C0396FF	\$59
Step 2. Surface Finish	
BK black	+\$0

Cash Drawer Insert

C0337



Product Information


Description
This insert fits inside an A-size drawer with lock for use as a cash storage drawer.

Dimensions


Specification Information

Step 1.
C033738 \$73


Co/Struc® Storage



TOP



FRONT




Product Information
Description
This assembly holds drawers in place on an L cart and indicates that the drawers are unaltered. It must be used with 2 C frames.
Notes
Order lock tabs separately.

Specification Information		
Step 1.		
C0473FF		\$260
Step 2. Surface Finish		
LU	soft white	+\$0

Label Clip


C0274
C0462



Product Information

Description
This clip snaps onto a drawer, tray/shelf, or flipper to hold content labels. The clip accepts pressure-sensitive labels from the front or slide-in labels from the back. Finish is clear. Package contains 12 clips; labels not included.

Dimensions



Specification Information

Step 1.

C0

Step 2. Width

27427 6" wide

462 19" wide

Prices for Steps 1-2.


C027427 \$20

C0462 \$39

Co/Struc® Storage

A-Size Drawer Subdivider


C0458




Product Information

Description
This subdivider has 35 separate compartments to organize small items within an A-size drawer. Finish is clear. Package contains 12.

Dimensions




TOP



FRONT

Specification Information

Step 1.
C0458 \$339



Product Information



Description

This assembly fits inside a locker or C frame or fits below a work surface with a drawer bearer. It has a wire cassette, cassette bins, and subdividers. The assembly fills 2 slots in a locker or C frame.

Notes

Order optional cassette cover (C0566) separately.

Dimensions




Co/Struc® Storage

Cassette Assembly *continued*

Specification Information		
Step 1.		
CO565FF		\$331
Step 2. Surface Finish		
LU	soft white	+\$0

Cassette Cover


C0566




Product Information
Description
This cover fits on a cassette assembly and provides top closure.
Notes
Order cassette assembly (C0565) separately.
Dimensions

Specification Information		
Step 1.		
C0566FF		\$73
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage




TOP



FRONT

Cassette Bin


C0567




Product Information
Description
This bin fits inside a cassette assembly. Package contains 10 bins, 20 dividers, and 10 label clips.
Notes
Order cassette assembly (C0565) separately.
Dimensions

Specification Information		
Step 1.		
C0567FF		\$116
Step 2. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage



BIN TOP



BIN FRONT


Dispensing Rail

C0365

C0366

C0367

C0368



Product Information

Description

This rail attaches to a panel, a freestanding module, or wall strips to support subcontainers at a 15° angle for gravity-feed dispensing.

Package contains 2.

The A-size rail holds A-size subcontainers; the B-size rail holds B-size subcontainers.

Notes

Dispensing rails hold the following number of subcontainers:

A-Size Rail Width—Number of Subcontainers

24" — 5 A-1, A-2, A-3 or A-4

24" — 2 A-5

24" — 10 A-8

48" — 10 A-1, A-2, A-3, or A-4

48" — 5 A-5

48" — 20 A-8

B-Size Rail Width—Number of Subcontainers

24" — 5 B-1 or B-3


24" — 2 B-5

48" — 10 B-1 or B-3


48" — 5 B-5

Order optional A-size subcontainer (C0231, C0232, C0233, C0234, C0235, or C0369) or B-size subcontainer (C0236, C0238, or C0240) separately.

Dimensions




SIDE




A-SIZE

FRONT



SIDE



B-SIZE

FRONT

Dispensing Rail *continued*

Specification Information

Step 1.

C036

Step 2. Size/Width

5FF a size and 24" wide

6FF a size and 48" wide

7FF b size and 24" wide

8FF b size and 48" wide

Prices for Steps 1-2.

C0365FF \$179

C0366FF \$205


C0367FF \$200

C0368FF \$226

Step 3. Surface Finish

LU soft white +\$0

Subcontainer



C0231
C0232
C0233
C0234
C0235
C0236
C0238
C0240
C0369

Product Information

Description

This storage container fits into a drawer or onto a dispensing rail to organize and store small items. Package contains 12.

The A-size subcontainer is used with an A-size drawer (without lid) or A-size dispensing rail. It is $2\frac{3}{16}$ " high; depth and width dimensions are:

Type—Depth—Width

A-1— $4\frac{9}{16}$ "— $4\frac{9}{16}$ "

A-2— $9\frac{3}{16}$ "— $4\frac{9}{16}$ "

A-3— $13\frac{13}{16}$ "— $4\frac{9}{16}$ "

A-4— $18\frac{7}{16}$ "— $4\frac{9}{16}$ "

A-5— $9\frac{3}{16}$ "— $9\frac{3}{16}$ "

A-8— $13\frac{13}{16}$ "— $2\frac{1}{4}$ "

The B-size subcontainer is used with a B-size drawer (without lid), B-size dispensing rail, or C-size drawer. It is $5\frac{3}{16}$ " high; depth and width dimensions are:

Type—Depth—Width

B-1— $4\frac{9}{16}$ "— $4\frac{9}{16}$ "

B-3— $13\frac{13}{16}$ "— $4\frac{9}{16}$ "


B-5— $9\frac{3}{16}$ "— $9\frac{3}{16}$ "

Notes


Order following products separately:

- Drawer (C0207, C0208, or C0209)
- Dispensing rail (C0365, C0366, C0367, or C0368)

Dimensions



TOP



A-SIZE
FRONT


B-SIZE
FRONT

Subcontainer *continued*

Specification Information		
Step 1.		
CO		
Step 2. Size		
231FF	a-1	
232FF	a-2	
233FF	a-3	
234FF	a-4	
235FF	a-5	
369FF	a-8	
236FF	b-1	
238FF	b-3	
240FF	b-5	
Prices for Steps 1-2.		
CO231FF		\$61
CO232FF		\$61
CO233FF		\$66
CO234FF		\$66
CO235FF		\$66
CO369FF		\$62
CO236FF		\$66
CO238FF		\$93
CO240FF		\$93
Step 3. Surface Finish		
LU	soft white	+\$0

A-Size Drawer Organizer Kit

CO540
CO541
CO542
CO543
CO544
CO545



Product Information

Description

This kit fits inside an A-size drawer (without lid) and has 6 combinations of subcontainers and subcontainer dividers.

A-size subcontainer is $2\frac{3}{16}$ " high. Width and depth dimensions are:

Type—Width—Depth

A-1— $4\frac{9}{16}$ "— $4\frac{9}{16}$ "

A-2— $4\frac{9}{16}$ "— $9\frac{3}{16}$ "



A-3— $4\frac{9}{16}$ "— $13\frac{13}{16}$ "

A-4— $4\frac{9}{16}$ "— $18\frac{7}{16}$ "

A-5— $9\frac{3}{16}$ "— $9\frac{3}{16}$ "


A-8— $2\frac{1}{4}$ "— $13\frac{13}{16}$ "

Dimensions




CO540FF


TOP




CO541FF




CO542FF



CO543FF



CO544FF



CO545FF

Specification Information

Step 1.

CO54

Step 2. Components

0FF 4 a-3 subcontainers, 6 a-3 dividers

1FF 1 a-3 subcontainer, 6 a-8 subcontainers, 2 a-3 dividers, 7 a-8 dividers

2FF 2 a-1 subcontainers, 2 a-2 subcontainers, 1 a-3 subcontainer, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers

3FF 3 a-3 subcontainers, 2 a-8 subcontainers, 4 a-3 dividers, 3 a-8 dividers

4FF 4 a-2 subcontainers, 1 a-4 subcontainer, 8 a-3 dividers

5FF 2 a-1 subcontainers, 2 a-3 subcontainers, 1 a-5 subcontainer, 4 a-3 dividers

Prices for Steps 1-2.


CO540FF	\$41
CO541FF	\$59
CO542FF	\$59
CO543FF	\$51
CO544FF	\$51
CO545FF	\$41

Step 3. Surface Finish

LU soft white	+\$0
----------------------	------

B-Size Drawer Organizer Kit

CO546
CO547



Product Information

Description

This kit fits inside a B-size drawer (without lid) or C-size drawer and has 2 combinations of subcontainers and subcontainer dividers.

The B-size subcontainer is 5 3/16" high. Width and depth dimensions are:


Type—Width—Depth

B-1—4 9/16"—4 9/16"

B-3—4 9/16"—13 13/16"

B-5—9 3/16"—9 3/16"

Dimensions



Specification Information

Step 1.

CO54

Step 2. Components

6FF 4 b-3 subcontainers, 4 b-3 dividers

7FF 2 b-1 subcontainers, 2 b-3 subcontainers, 1 b-5 subcontainer, 4 b-3 dividers

Prices for Steps 1-2.

CO546FF \$59


CO547FF \$59

Step 3. Surface Finish

LU soft white +\$0

Subcontainer Divider

C0378
C0379
C0380




Product Information

Description
This divider fits into subcontainers to separate and organize stored items. A divider can be adjusted by sliding it to the desired position. The divider is not available for A-5 and B-5 subcontainers. Package contains 12.

The divider for an A-8 subcontainer is 1⁵/₈" high and 2¹/₈" wide; the divider for an A-1, A-2, A-3, or A-4 subcontainer is 1⁵/₈" high and 4¹/₄" wide; and the divider for a B-1 or B-3 subcontainer is 4" high and 4¹/₄" wide.

Dimensions



SIDE

FRONT

Specification Information

Step 1.

C03

Step 2. Usage

- 78FF** for a-8 subcontainer
- 79FF** for a-1, a-2, a-3, or a-4 subcontainers
- 80FF** for b-1 or b-3 subcontainers

Prices for Steps 1-2.


C0378FF	\$56
C0379FF	\$56
C0380FF	\$58

Step 3. Surface Finish

LU	soft white	+\$0
-----------	------------	------

Subcontainer Handle

C0376
C0377



Product Information


Description
This snap-on handle has a surface for identification labels in a picking station. It also can convert a small or large subcontainer into a drawer when used on the flat surface of a C tray. The handle cannot be used with a subcontainer in an A-size drawer. Package contains 12.

Notes
Handle accepts following label:

- 2 1/8"-wide handle: 3/4" x 1" gummed label
- 4 1/2"-wide handle: 3/4" x 3" drop-in label

For 4 1/2"-wide handle, order label (C0459) separately.

Dimensions




Specification Information

Step 1.		
C037		
Step 2. Width		
6FF	2 1/8" wide	
7FF	4 1/2" wide	
Prices for Steps 1-2.		
C0376FF		\$28
C0377FF		\$34
Step 3. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage

Subcontainer Handle Label

C0459




Product Information
Description
This label slides into a 4½"-wide subcontainer handle. Package contains 100 cards with 10 perforated labels per card (1000 labels total).

Specification Information
Step 1.
C0459
\$24

Co/Struc® Storage

Subcontainer Label Flag


C0381
C0382



Product Information

Description
This flag snaps onto subcontainers that are hung from a dispensing rail. It provides a visible surface for identification labels. Package contains 12.
The flag for an A-8 subcontainer is 1 1/4" wide; the flag for an A-3 or B-3 subcontainer is 3 1/2" wide.

Dimensions




Specification Information

Step 1.		
C038		
Step 2. Usage		
1FF	for a-8 subcontainer	
2FF	for a-3 or b-3 subcontainer	
Prices for Steps 1-2.		
C0381FF		\$30
C0382FF		\$41
Step 3. Surface Finish		
LU	soft white	+\$0

Co/Struc® Storage


Subcontainer HangerC0389




Product Information
Description
This component hangs on pegs on either side of a locker. It holds 2 A-3 subcontainers or 1 A-3 and 1 B-3 subcontainer.
Dimensions

Specification Information		
Step 1.		
C0389FF		\$74
Step 2. Surface Finish		
LU	soft white	+\$0


Co/Struc® Storage



TOP




SIDE



FRONT

Pullout Catheter Rack


CO469.




Product Information
Description
This 20"-wide, powder-coated wire rack fits into the slots of a locker or C frame to organize catheters. It is available with 3 or 4 independent, pullout slides. Each pullout slide has 4 prongs, 4 side label clips, and 1 front label holder. Each prong holds up to 10 catheters. Finish is soft white.
Dimensions

Specification Information
Step 1.
CO469.20
Step 2. Slides
03 3 slides
04 4 slides
Prices for Steps 1-2.
CO469.20 03 \$478
04 \$515

Co/Struc® Storage






Product Information
Description
This rack mounts on the outside of a locker or cart that is at least 42" high. It holds 10 Fogarty catheters. Mounting hardware is included.
Dimensions


Specification Information
Step 1.
CO467FF
\$114
Step 2. Surface Finish
LU soft white
+\$0

Co/Struc® Storage



Chest Tube Rack

C0466





Product Information


Description
This rack mounts on the outside of a locker or cart. It holds 14 chest tubes. Mounting hardware is included.
Dimensions

Specification Information

Step 1.	
C0466FF	\$114
Step 2. Surface Finish	
LU soft white	+\$0

Co/Struc® Storage







Product Information

Description
This cart is used to maneuver and transport lockers. It has a foot-activated lever that lifts a locker from a rail to a cart or from a cart to a rail. The cart also has a towing handle and a train connector for linking multiple carts. A maximum of 3 carts can be linked when pulled manually; a maximum of 7 carts can be linked when pulled with a tug. Shipped knocked down.

Dimensions



TOP



SIDE


Specification Information

Step 1.

C0342FF \$1349

Step 2. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0



Product Information

Description
This mobile cart has a platform and 4 rail profiles for hanging frames. It has 4"-diameter, hard-rubber, swivel zinc casters.


Notes
To hang drawer or shelf storage components, order appropriate frame separately:

- C frame (CO205)
- CST frame (CO206)

Order optional components separately:

- Breakaway lock assembly (CO473)
- IV pole (CO478)
- IV pole clamp (CO465)

Dimensions




Specification Information

Step 1.
CO212FF \$582

Step 2. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

IV Pole


C0478



Product Information
Description
This 3-ring pole mounts on an L cart and holds IV bags.
Notes
Order IV pole clamp (C0465) separately.
Dimensions


Specification Information
Step 1.
C0478
\$260

Co/Struc® Transport




IV Pole Clamp

C0465



Product Information
Description
This clamp mounts on the back left or right side of an L cart. It holds a 3/4"-diameter IV pole. Mounting hardware is included.
Dimensions

Specification Information		
Step 1.		
CO465FF		\$99
Step 2. Surface Finish		
LU	soft white	+\$0




Storage

Work Surfaces



Base Cabinet with Drawers

CG601.
CG603.



Product Information

Description


This 35"-high, preconfigured cabinet has 4 drawers, a full back, and adjustable glides.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.

Lock option is for top drawer only. For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG60

Step 2. Configuration

- 4 8"-high drawers
- 1 4"-high drawer, 2 8"-high drawers, 1 12"-high drawer

Step 3. Depth

- 18** 18" deep
24 24" deep
30 30" deep

Step 4. Width

- 18** 18" wide
24 24" wide
30 30" wide

Step 5. Lock

- N** no lock
L keyed differently
K keyed alike

Prices for Steps 1-5.

		N	L	K
CG601.	18 18	\$890	920	920
	24	\$1002	1032	1032
	30	\$1114	1144	1144
	24 18	\$938	968	968
	24	\$1056	1086	1086
	30	\$1175	1205	1205
	30 18	\$989	1019	1019
	24	\$1108	1138	1138
	30	\$1226	1256	1256
CG603.	18 18	\$890	920	920
	24	\$1002	1032	1032
	30	\$1114	1144	1144
	24 18	\$938	968	968
	24	\$1056	1086	1086
	30	\$1175	1205	1205
	30 18	\$989	1019	1019
	24	\$1108	1138	1138
	30	\$1226	1256	1256

Base Cabinet with Drawers *continued*


Casework Storage

Step 6. Cabinet Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0
Step 7. Pull Option/Finish		
NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
Step 8. Drawer Edge Finish		
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>		
Price Category 1		+\$0
Step 9. Drawer Surface Finish		
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>		
Price Category 1		+\$0
Price Category 2		+\$64

Base Cabinet with Door

CG610.

CG611.



Product Information

Description


This 35"-high, preconfigured cabinet has an adjustable shelf, a door, a base cabinet bottom, a full back, and adjustable glides. The shelf is height adjustable in 1 1/4" increments.

Notes


Order work surface (CG210.) or corner work surface (CG217.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



TOP



SIDE

LEFT-HAND
DOOR SWING
FRONT

RIGHT-HAND
DOOR SWING
FRONT

Specification Information

Step 1.

CG61

Step 2. Door Swing

- 0. left-hand door swing
- 1. right-hand door swing

Step 3. Depth

- 18 18" deep
- 24 24" deep
- 30 30" deep

Step 4. Width

- 18 18" wide
- 24 24" wide

Step 5. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-5.

	N	L	K
CG610. 18 18	\$564	604	604
24	\$635	675	675
24 18	\$613	653	653
24	\$690	730	730
30 18	\$678	718	718
24	\$759	799	799
CG611. 18 18	\$564	604	604
24	\$635	675	675
24 18	\$613	653	653
24	\$690	730	730
30 18	\$678	718	718
24	\$759	799	799

Step 6. Cabinet Surface Finish

- HF inner tone light +\$0
- LT light tone +\$0

Step 7. Pull Option/Finish

- NN none +\$0
- BK black +\$0
- HF inner tone light +\$0
- LT light tone +\$0


Casework Storage

Step 8. Door Edge Finish	
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>	
Price Category 1	+\$0

Step 9. Door Surface Finish	
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>	
Price Category 1	+\$0
Price Category 2	+\$15

Base Cabinet with Double Doors

CG612.
CG613.



Product Information


Description

This 35"-high, preconfigured cabinet has 2 doors, a base cabinet bottom, a full back, and adjustable glides. It can be specified with 1 full-width shelf or with a fixed center divider and 2 half-width shelves. The shelves are height adjustable in 1 $\frac{1}{4}$ " increments.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.
For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Base Cabinet with Double Doors

continued


Casework Storage

Specification Information				
Step 1.				
CG61				
Step 2. Configuration				
2.	1 shelf			
3.	center divider and 2 shelves			
Step 3. Depth				
18	18" deep			
24	24" deep			
30	30" deep			
Step 4. Width				
For 1 shelf (2.)				
30	30" wide			
48	48" wide			
For center divider and 2 shelves (3.)				
48	48" wide			
60	60" wide			
Step 5. Lock				
N	no lock			
L	keyed differently			
K	keyed alike			
Prices for Steps 1-5.				
		N	L	K
CG612. 18 30		\$712	752	752
	48	\$907	947	947
	24 30	\$774	814	814
	48	\$986	1026	1026
	30 30	\$846	886	886
	48	\$1071	1111	1111
CG613. 18 48		\$1110	1150	1150
	60	\$1383	1423	1423
	24 48	\$1206	1246	1246
	60	\$1503	1543	1543
	30 48	\$1320	1360	1360
	60	\$1631	1671	1671
Step 6. Cabinet Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0

Step 7. Pull Option/Finish		
NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0
Step 8. Door Edge Finish		
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>		
Price Category 1		+\$0
Step 9. Door Surface Finish		
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>		
Price Category 1		+\$0
Price Category 2		+\$30

Base Cabinet with Drawer and Door

CG620.
CG621.



Product Information

Description

This 35"-high, preconfigured cabinet has an 8"-high drawer, an adjustable shelf, a door, a base cabinet bottom, a full back, and adjustable glides. The shelf is height adjustable in 1¹/₄" increments.


Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.

Base cabinet with lock includes lock on drawer and door.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG62

Step 2. Door Swing

- 0. left-hand door swing
- 1. right-hand door swing

Step 3. Depth

- 24 24" deep
- 30 30" deep

Step 4. Width

- 18 18" wide
- 24 24" wide

Step 5. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-5.

	N	L	K
CG620. 24 18	\$738	842	842
24	\$830	924	924
30 18	\$803	902	902
24	\$899	1016	1016
CG621. 24 18	\$738	842	842
24	\$830	924	924
30 18	\$803	902	902
24	\$899	1016	1016

Step 6. Cabinet Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 7. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 8. Door/Drawer Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------


Base Cabinet with Drawer and Door *continued*

Casework Storage

Step 9. Door/Drawer Surface Finish	
<i>See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.</i>	
Price Category 1	+\$0
Price Category 2	+\$31

Base Cabinet with Drawer and Double Doors

CG622.



Product Information


Description

This 35"-high, 30"-wide, preconfigured cabinet has an 8"-high drawer, an adjustable shelf, 2 doors, a base cabinet bottom, a full back, and adjustable glides. The shelf is height adjustable in 1¹/₄" increments.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.
For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG622.

Step 2. Depth

1830 18" deep
2430 24" deep
3030 30" deep

Step 3. Lock

N no lock
L keyed differently
K keyed alike

Prices for Steps 1-3.

	N	L	K
CG622. 1830	\$855	1065	1065
2430	\$929	1034	1034
3030	\$1002	1140	1140

Step 4. Cabinet Surface Finish

HF inner tone light +\$0
LT light tone +\$0

Step 5. Pull Option/Finish

NN none +\$0
BK black +\$0
HF inner tone light +\$0
LT light tone +\$0

Step 6. Door/Drawer Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Step 7. Door/Drawer Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Price Category 2 +\$46

Sink Cabinet with Foot Pedal Cutout and Door

CG630.
CG631.



Casework Storage

Product Information

Description


This 35"-high, 24"-wide preconfigured cabinet has a door, a base cabinet bottom with a foot pedal cutout, a full back, and adjustable glides.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.
Sink must be customer supplied and field installed.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG63

Step 2. Door Swing

- 0. left-hand door swing
- 1. right-hand door swing

Step 3. Depth

- 1824 18" deep
- 2424 24" deep
- 3024 30" deep

Step 4. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-4.

	N	L	K
CG630. 1824	\$574	611	611
2424	\$608	642	642
3024	\$664	698	698
CG631. 1824	\$574	611	611
2424	\$608	642	642
3024	\$664	698	698

Step 5. Cabinet Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 6. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 7. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Sink Cabinet with Foot Pedal Cutout and Door *continued*


Step 8. Door Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$15

Sink Cabinet with Foot Pedal Cutout and Double Doors

CG632.



Casework Storage

Product Information

Description


This 35"-high, preconfigured cabinet has 2 doors, a base cabinet bottom with a foot pedal cutout, a full back, and adjustable glides.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.
Sink must be customer supplied and field installed.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG632.

Step 2. Depth

24 24" deep
30 30" deep

Step 3. Width

30 30" wide
48 48" wide

Step 4. Lock

N no lock
L keyed differently
K keyed alike

Prices for Steps 1-4.

	N	L	K
CG632. 24 30	\$677	711	711
48	\$841	873	873
30 30	\$732	766	766
48	\$904	936	936

Step 5. Cabinet Surface Finish

HF inner tone light +\$0
LT light tone +\$0

Step 6. Pull Option/Finish

NN none +\$0
BK black +\$0
HF inner tone light +\$0
LT light tone +\$0

Step 7. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 8. Door Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.


Price Category 1 +\$0

Price Category 2 +\$30

ADA Sink Support, Wall Mounted CG508.



Product Information
Description
This 24"-deep sink support hangs from an architectural wall. It supports a 5½"-deep sink and a Casework work surface or alternative material top. The sink support has a 4"-high laminate front trim and a removable front panel to allow access to plumbing. Shipped knocked down.
Notes
Order Casework work surface (CG210.) separately. Sink must be customer supplied and field installed. A Just SL-ADA-1921-A-GR-3 sink, with drain located in rear center, fits in a 24" × 24" sink support.
Dimensions




Specification Information
Step 1.
CG508.24
Step 2. Width
24 24" wide
30 30" wide
48 48" wide
60 60" wide
Prices for Steps 1-2.
CG508.24 24 \$700
30 \$750
48 \$900
60 \$1000
Step 3. Base Finish
HF inner tone light +\$0
LT light tone +\$0
Step 4. Edge Finish
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.
Price Category 1 +\$0
Step 5. Front Finish
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.
Price Category 1 +\$0
Price Category 2 +\$50

Casework Storage

Storage Cabinet with Shelves and Door

CG710.
CG711.



Product Information

Description


This 80"-high, 24"-wide × 24"-deep preconfigured cabinet has 4 adjustable shelves, a 76"-high door, a base cabinet bottom, a full back, a top, and adjustable glides. The shelves are height adjustable in 1¼" increments.

Notes

Order additional base cabinet shelves (CG515.) separately.

Order optional sloped overhead cabinet top (CG397.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.



Storage Cabinet with Shelves and Door *continued*

Specification Information

Step 1.

CG71

Step 2. Door Swing

0.2424 left-hand door swing

1.2424 right-hand door swing

Step 3. Lock

N no lock

L keyed differently

K keyed alike

Prices for Steps 1-3.

CG710.2424	N	\$1507
	L	\$1546
	K	\$1546

CG711.2424	N	\$1507
	L	\$1546
	K	\$1546

Step 4. Cabinet Surface Finish

HF inner tone light +\$0

LT light tone +\$0

Step 5. Pull Option/Finish

NN none +\$0

BK black +\$0

HF inner tone light +\$0

LT light tone +\$0

Step 6. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 7. Door Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Price Category 2 +\$15

Storage Cabinet with Shelves and Double Doors CG712.



Casework Storage

Product Information

Description

This 80"-high, preconfigured cabinet has 4 adjustable shelves, 2 76"-high doors, a base cabinet bottom, a full back, a top, and adjustable glides. The shelves are height adjustable in 1 $\frac{1}{4}$ " increments.



Notes

Order additional base cabinet shelves (CG515.) separately.

Order optional sloped overhead cabinet top (CG397.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG712.

Step 2. Depth

18	18" deep
24	24" deep

Step 3. Width

30	30" wide
48	48" wide

Step 4. Lock

N	no lock
L	keyed differently
K	keyed alike

Prices for Steps 1-4.

	N	L	K
CG712. 18 30	\$1490	1530	1530
48	\$1978	2018	2018
24 30	\$1705	1744	1744
48	\$2243	2283	2283

Step 5. Cabinet Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 6. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 7. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------


Step 8. Door Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$30

Storage Cabinet

CG506.



Product Information

Description

This 80"-high cabinet has an open front, a full back, a top, and adjustable glides. It can be specified with doors, drawers, and shelves.

Notes

Order interior components separately:

- Base cabinet drawer (CG530.) for 18"-deep cabinet or base cabinet drawer (CG531.) for 24"-deep cabinet
- Base cabinet shelf (CG515.)

To enclose bottom of cabinet, order base cabinet bottom (CG511.) separately.

To enclose front of cabinet, order storage cabinet door (CG545. or CG546.) or storage cabinet double doors (CG547.) separately. Door cannot be used in front of cabinet drawer.


For cabinet with door and drawers, order 44"-high storage cabinet door (CG545.4424 or CG546.4424) or storage cabinet double doors (CG547.44), order at least 1 base cabinet shelf (CG515.) for behind the door, and order 1 of the following base cabinet drawer combinations:

- 4 8"-high drawers (CG530.08 or CG531.08)
- 1 4"-high drawer (CG530.04 or CG531.04), 2 8"-high drawers (CG530.08 or CG531.08), and 1 12"-high drawer (CG530.12 or CG531.12)
- 1 8"-high drawer (CG530.08 or CG531.08) and 2 12"-high drawers (CG530.12 or CG531.12)
- 2 4"-high drawers (CG530.04 or CG531.04) and 2 12"-high drawers (CG530.12 or CG531.12)


To enclose drawer and allow locking capability, order base cabinet shelf (CG515.) for each drawer separately.

Order optional sloped overhead cabinet top (CG397.) separately.

Dimensions



TOP



SIDE


FRONT

Casework Storage

Specification Information			
Step 1.			
CG506.			
Step 2. Depth			
18	18" deep		
24	24" deep		
Step 3. Width			
24	24" wide		
30	30" wide		
48	48" wide		
Prices for Steps 1-3.			
		24	3048
CG506. 18		\$620	654734
24		\$729	769863
Step 4. Surface Finish			
HF	inner tone light		+\$0
LT	light tone		+\$0

Base Cabinet

CG501.
CG503.



Product Information

Description

This 35"-high cabinet has an open front, a full back, and adjustable glides. It can be specified with doors, drawers, and shelves.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.

Order interior components separately:

- Base cabinet bottle drawer (CG536.)
- Base cabinet bottom (CG511. or CG512.)
- Base cabinet drawer (CG530.) for 18"-deep cabinet or base cabinet drawer (CG531.) for 24"- or 30"-deep cabinet
- Base cabinet pullout shelf (CG521.)
- Base cabinet shelf (CG515.)
- Interior cabinet adapter (CG590.)

To enclose 18"- to 24"-wide cabinet, order base cabinet door (CG540. or CG541.) separately.

To enclose 30"-wide cabinet, order base cabinet double doors (CG542.) separately.


To enclose 48"-wide cabinet, order base cabinet doors (CG540. and CG541.) or base cabinet double doors (CG542.) separately.

To enclose 24"-wide cabinet containing interior cabinet adapter, order base cabinet door, interior cabinet adapter (CG580. or CG581.) separately.

To enclose 60"-wide cabinet, order 2 sets of base cabinet double doors (CG542.3230) separately.

Doors cannot be used in front of cabinet drawers.

Dimensions



Specification Information

Step 1.

CG50

Step 2. Configuration

1. open
3. center divider

Step 3. Depth

- 18** 18" deep
24 24" deep
30 30" deep

Step 4. Width

For open (1.)

- 18** 18" wide
24 24" wide
30 30" wide
48 48" wide

For center divider (3.)

- 48** 48" wide
60 60" wide

Prices for Steps 1-4.


	18	24	30	48	60
CG501. 18	\$345	371	397	458	—
24	\$406	436	467	539	—
30	\$460	491	520	593	—
CG503. 18	—	—	—	\$662	702
24	—	—	—	\$779	826
30	—	—	—	\$865	912

Step 5. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Base Cabinet Bottom

CG511.
CG512.



Product Information


Description

This component encloses the bottom of a base cabinet or storage cabinet. The bottom with a cutout allows installation of foot pedals for a cabinet used as a sink unit. Attachment hardware is included.

Notes

Order base cabinet (CG501. or CG503.) or storage cabinet (CG506.) separately.

Dimensions




Base Cabinet Bottom *continued*

Casework Storage

Specification Information					
Step 1.					
CG51					
Step 2. Cutout					
1.	no cutout				
2.	foot pedal cutout				
Step 3. Depth					
18	18" deep				
24	24" deep				
30	30" deep				
Step 4. Width					
For no cutout (1.)					
18	18" wide				
24	24" wide				
30	30" wide				
48	48" wide				
For foot pedal cutout (2.)					
24	24" wide				
30	30" wide				
48	48" wide				
Prices for Steps 1-4.					
		18	24	30	48
CG511. 18		\$82	94	105	140
24		\$102	117	131	175
30		\$113	131	149	204
CG512. 18		—	\$98	110	144
24		—	\$122	137	180
30		—	\$137	155	210
Step 5. Surface Finish					
HF	inner tone light				+\$0
LT	light tone				+\$0

Base Cabinet Shelf

CG515.



Product Information


Description

This metal shelf mounts inside a base cabinet or storage cabinet. It is height adjustable in 1 1/4" increments. Attachment hardware is included.

Notes

Order base cabinet (CG501. or CG503.) or storage cabinet (CG506.) separately.

Dimensions



Specification Information

Step 1.

CG515.

Step 2. Depth

18	18" deep
24	24" deep
30	30" deep

Step 3. Width

18	18" wide
24	24" wide
30	30" wide
48	48" wide


Prices for Steps 1-3.

	18	24	30	48
CG515. 18	\$59	71	82	114
24	\$74	88	102	143
30	\$85	102	119	170

Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0


Base Cabinet Pullout ShelfCG521.



Casework Storage


Product Information
Description
This metal, pullout shelf mounts inside a 24"- or 30"-deep base cabinet or 24"-deep storage cabinet. It is height adjustable in 1¼" increments. The shelf cannot attach at a door hinge location. Attachment hardware is included.
Notes
Order base cabinet (CG501.24, CG501.30, CG503.24, or CG503.30) or storage cabinet (CG506.24) separately.
Dimensions

Specification Information
Step 1.
CG521.24
Step 2. Width
1818" wide
2424" wide
3030" wide
Prices for Steps 1-2.
CG521.24 18\$112
24\$132
30\$150
Step 3. Surface Finish
HFinner tone light+\$0
LTlight tone+\$0



Base Cabinet Drawer

CG530.
CG531.



Product Information

Description

This drawer mounts inside a base cabinet or storage cabinet.
Attachment hardware is included.


Notes

Order base cabinet (CG501. or CG503.) or storage cabinet (CG506.) separately.

When ordering locking drawer, order base cabinet shelf (CG515.) separately (except for drawer in top position).

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG53

Step 2. Depth

- 0. 15 3/4" deep
- 1. 21 3/4" deep

Step 3. Height

- 04 4" high
- 08 8" high
- 12 12" high

Step 4. Width

- 18 18" wide
- 24 24" wide
- 30 30" wide

Step 5. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-5.

	N	L	K
CG530. 04 18	\$120	149	149
24	\$140	170	170
30	\$160	189	189
08 18	\$133	163	163
24	\$155	184	184
30	\$177	207	207
12 18	\$148	177	177
24	\$171	200	200
30	\$194	223	223
CG531. 04 18	\$122	152	152
24	\$143	173	173
30	\$163	193	193
08 18	\$136	166	166
24	\$158	188	188
30	\$181	211	211
12 18	\$151	181	181
24	\$174	204	204
30	\$198	228	228

Base Cabinet Drawer *continued*

Casework Storage

Step 6. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 7. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------


Step 8. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$16

Base Cabinet Bottle Drawer

CG536.




Product Information

Description
This 6"-high drawer mounts in a base cabinet or storage cabinet. It has a partial-height front for easy access to contents. The drawer is available in 3 widths. The 15⁵/₈"-wide drawer includes 2 adjustable dividers; the 21⁵/₈"-wide drawer includes 3 adjustable dividers; and the 27⁵/₈"-wide drawer includes 4 adjustable dividers. Attachment hardware is included.

Notes
Order base cabinet (CG501.24, CG501.30, CG503.24, or CG503.30) or storage cabinet (CG506.24) separately.

Dimensions




Specification Information

Step 1.		
CG536.08		
Step 2. Width		
18	15 ⁵ / ₈ " wide	
24	21 ⁵ / ₈ " wide	
30	27 ⁵ / ₈ " wide	
Prices for Steps 1-2.		
CG536.08	18	\$170
	24	\$198
	30	\$226
Step 3. Edge Finish		
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.		
Price Category 1		+\$0
Step 4. Surface Finish		
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.		
Price Category 1		+\$0
Price Category 2		+\$16

Casework Storage

Base Cabinet Door

CG540.
CG541.



Casework Storage

Product Information

Description

This door encloses a base cabinet. Attachment hardware is included.


Notes

Order base cabinet (CG501. or CG503.) separately.

To enclose base cabinet without drawers, specify 32"-high door. To enclose base cabinet with 4"-high drawer, specify 28"-high door. To enclose base cabinet with 8"-high drawer or 2 4"-high drawers, specify 24"-high door.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG54

Step 2. Door Swing

- 0.** left-hand door swing
- 1.** right-hand door swing

Step 3. Height

- 24** 24" high
- 28** 28" high
- 32** 32" high

Step 4. Width

- 18** 18" wide
- 24** 24" wide

Step 5. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

Prices for Steps 1-5.

	N	L	K
CG540. 24 18	\$98	142	142
24	\$127	171	171
28 18	\$112	156	156
24	\$147	190	190
32 18	\$127	171	171
24	\$166	210	210
CG541. 24 18	\$98	142	142
24	\$127	171	171
28 18	\$112	156	156
24	\$147	190	190
32 18	\$127	171	171
24	\$166	210	210

Step 6. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Base Cabinet Door *continued*

Step 7. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Step 8. Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Price Category 2	+\$15
------------------	-------

Base Cabinet Door, Interior Cabinet Adapter

CG580.
CG581.



Product Information

Description

This 24"-wide door encloses a base cabinet containing an interior cabinet adapter. Attachment hardware is included.


Notes

Order base cabinet (CG501.2424 or CG501.3024) or base cabinet with center divider (CG503.2448 or CG503.3048) and base cabinet bottom (CG511.) separately.

Order interior cabinet adapter (CG590.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG58

Step 2. Door Swing

0.3224 left-hand door swing

1.3224 right-hand door swing

Step 3. Lock

N no lock

L keyed differently

K keyed alike

Prices for Steps 1-3.

CG580.3224	N	\$324
	L	\$434
	K	\$434

CG581.3224	N	\$324
	L	\$434
	K	\$434

Step 4. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 5. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------


Step 6. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$15

Base Cabinet Double Doors

CG542.



Product Information

Description

These doors enclose a base cabinet. Attachment hardware is included.

Notes


Order base cabinet (CG501. or CG503.) separately.

To enclose base cabinet without drawers, specify 32"-high doors. To enclose base cabinet with 4"-high drawer, specify 28"-high doors. To enclose base cabinet with 8"-high drawer or 2 4"-high drawers, specify 24"-high doors.

Double doors with lock option are controlled by 1 lock.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG542.

Step 2. Height

24	24" high
28	28" high
32	32" high

Step 3. Width

30	30" wide
48	48" wide

Step 4. Lock

N	no lock
L	keyed differently
K	keyed alike

Prices for Steps 1-4.

	N	L	K
CG542. 24 30	\$176	216	216
48	\$233	273	273
28 30	\$184	224	224
48	\$233	273	273
32 30	\$197	237	237
48	\$306	346	346

Step 5. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------


Step 7. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$30




Base Cabinet Filler

CG598.
CG599.



Casework Storage


Product Information
Description
This 6"-wide component fills space in a run of base cabinets. The straight filler fills the space between a run of base cabinets and the wall. The corner filler fills the corner space where 2 runs of base cabinets meet. Attachment hardware is included.
Dimensions



Specification Information
Step 1.
CG59
Step 2. Type
8. straight
9. corner
Step 3. Height
For straight (8.)
36 36" high
80 80" high
For corner (9.)
36 36" high
Prices for Steps 1-3.
CG598. 36 \$110
80 \$198
CG599. 36 \$88
Step 4. Surface Finish
HF inner tone light +\$0
LT light tone +\$0

Interior Cabinet Adapter

CG590.




Product Information

Description
This adapter mounts inside a 24"-wide base cabinet. It has 9 slots at 3" intervals to support Co/Struc® storage components and accessories. The adapter can hold lockable drawers. Attachment hardware is included.

Notes
Order base cabinet (CG501.2424 or CG501.3024) or base cabinet with center divider (CG503.2448 or CG503.3048) and base cabinet bottom (CG511.) separately.
Order optional base cabinet door, interior cabinet adapter (CG580. or CG581.) separately.
Order storage components separately:
• A-, B-, or C-size drawer (CO207, CO208, or CO209)
• C tray/shelf (CO203)
• C wire shelf (CO252)

Dimensions




Specification Information

Step 1.		
CG590.3218		\$302
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0

Casework Storage

Storage Cabinet Door

CG545.
CG546.



Casework Storage

Product Information

Description

This 24"-wide door encloses a storage cabinet. Attachment hardware is included.


Notes

Order storage cabinet (CG506.) separately.

For storage cabinet with drawers installed, specify 44"-high door. See Storage Cabinet for available drawer combinations.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG54

Step 2. Door Swing

- 5. left-hand door swing
- 6. right-hand door swing

Step 3. Height

- 4424** 44" high
- 7624** 76" high

Step 4. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

Prices for Steps 1-4.

	N	L	K
CG545. 4424	\$207	247	247
7624	\$352	392	392
CG546. 4424	\$207	247	247
7624	\$352	392	392

Step 5. Pull Option/Finish

NN	none		+\$0
BK	black		+\$0
HF	inner tone light		+\$0
LT	light tone		+\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 7. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0
Price Category 2 +\$15

Storage Cabinet Double Doors

CG547.



Product Information

Description

These doors enclose a storage cabinet. Attachment hardware is included.


Notes

Order storage cabinet (CG506.) separately.

For storage cabinet with drawers installed, specify 44"-high door. See Storage Cabinet for available drawer combinations.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG547.

Step 2. Height

44 44" high

76 76" high

Step 3. Width

For 44" high (44)

30 30" wide

For 76" high (76)

30 30" wide

48 48" wide

Step 4. Lock

N no lock

L keyed differently

K keyed alike

Prices for Steps 1-4.

	N	L	K
CG547. 44 30	\$381	425	425
76 30	\$490	534	534
48	\$773	817	817

Step 5. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 7. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Price Category 2 +\$30

CG321.



Dimensions

Overhead Cabinet with Door

continued

Casework Storage

Specification Information

Step 1.

CG3

Step 2. Configuration/Height

- 1** 1 shelf and 24" high
- 2** 2 shelves and 34" high

Step 3. Door Swing

- 0.** left-hand door swing
- 1.** right-hand door swing

Step 4. Depth

- 0924** 10¹/₂" deep
- 1624** 17" deep

Step 5. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

Prices for Steps 1-5.

	N	L	K
CG310. 0924	\$425	465	465
1624	\$501	541	541
CG311. 0924	\$425	465	465
1624	\$501	541	541
CG320. 0924	\$572	616	616
1624	\$694	738	738
CG321. 0924	\$572	616	616
1624	\$694	738	738

Step 6. Cabinet Surface Finish

- HF** inner tone light +\$0
- LT** light tone +\$0

Step 7. Pull Option/Finish

- NN** none +\$0
- BK** black +\$0
- HF** inner tone light +\$0
- LT** light tone +\$0

Step 8. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Step 9. Door Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0

Price Category 2 +\$15

Overhead Cabinet with Double Doors

CG312.
CG322.



Casework Storage

Product Information

Description

This preconfigured cabinet attaches to an architectural wall. It has 1 or 2 adjustable shelves, 2 doors, a full back, and a flat top. The shelves attach in a flat or slanted position and are height adjustable in 1/4" increments.




Notes

Order optional add-on shelf (SM900.) separately.

Order optional sloped overhead cabinet top (CG397.) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG3

Step 2. Configuration/Height

- 12. 1 shelf and 24" high
- 22. 2 shelves and 34" high

Step 3. Depth

- 09 10 1/2" deep
- 16 17" deep

Step 4. Width

- 30 30" wide
- 36 36" wide
- 48 48" wide

Step 5. Lock

- N no lock
- L keyed differently
- K keyed alike

Prices for Steps 1-5.

	N	L	K
CG312. 09 30	\$508	548	548
36	\$546	586	586
48	\$611	651	651
16 30	\$586	626	626
36	\$625	665	665
48	\$706	746	746
CG322. 09 30	\$657	697	697
36	\$698	738	738
48	\$779	819	819
16 30	\$784	824	824
36	\$834	874	874
48	\$935	975	975

Step 6. Cabinet Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0

Step 7. Pull Option/Finish

NN	none	+\$0
BK	black	+\$0
HF	inner tone light	+\$0
LT	light tone	+\$0

Overhead Cabinet with Double Doors *continued*

Step 8. Door Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Step 9. Door Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Price Category 2	+\$30
------------------	-------

Overhead Cabinet Door

CG390.
CG391.



Casework Storage

Product Information

Description

This 24"-wide door encloses a shelf unit. Attachment hardware is included.


Notes

For 24"-high door, order 3-shelf unit (SM303.0924 or SM303.1624) separately.

For 34"-high door, order 5-shelf unit (SM305.0924 or SM305.1624) separately.

For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.

Dimensions



Specification Information

Step 1.

CG39

Step 2. Door Swing

- 0.** left-hand door swing
- 1.** right-hand door swing

Step 3. Height

- 2424** 24" high
- 3424** 34" high

Step 4. Lock

- N** no lock
- L** keyed differently
- K** keyed alike

Prices for Steps 1-4.

	N	L	K
CG390. 2424	\$129	169	169
3424	\$178	218	218
CG391. 2424	\$129	169	169
3424	\$178	218	218

Step 5. Pull Option/Finish

NN	none		+\$0
BK	black		+\$0
HF	inner tone light		+\$0
LT	light tone		+\$0

Step 6. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.


Price Category 1 +\$0

Step 7. Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0
Price Category 2 +\$15

Overhead Cabinet Double Doors CG392.




Product Information
Description
These doors enclose a shelf unit. Attachment hardware is included.
Notes
For 24"-high doors, order 3-shelf unit (SM303.) separately.
For 34"-high doors, order 5-shelf unit (SM305.) separately.
For keyed-alike locks, order lock plugs separately. See Keyed-Alike Information in Appendices.
Dimensions



Specification Information			
Step 1.			
CG392.			
Step 2. Height			
24	24" high		
34	34" high		
Step 3. Width			
30	30" wide		
36	36" wide		
48	48" wide		
Step 4. Lock			
N	no lock		
L	keyed differently		
K	keyed alike		
Prices for Steps 1-4.			
	N	L	K
CG392. 24 30	\$177	217	217
36	\$207	247	247
48	\$256	296	296
34 30	\$230	270	270
36	\$271	311	311
48	\$356	396	396
Step 5. Pull Option/Finish			
NN	none		+\$0
BK	black		+\$0
HF	inner tone light		+\$0
LT	light tone		+\$0
Step 6. Edge Finish			
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.			
Price Category 1			+\$0
Step 7. Surface Finish			
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.			
Price Category 1			+\$0
Price Category 2			+\$25

CG393.



Product Information

This panel encloses the back of a 3- or 5-shelf unit and allows the shelf unit to hang on an architectural wall.

Order 3-shelf unit (SM303.) or 5-shelf unit (SM305.) separately.

Specification Information

CG393.

24 24" high

34 34" high

24 24" wide

30 30" wide

36 36" wide

48 48" wide

Prices for Steps 1-3.

	24	30	36	48
CG393. 24	\$110	135	180	210
34	\$120	140	185	220


Step 4. Surface Finish

HF	inner tone light	+\$0
LT	light tone	+\$0



Sloped Overhead Cabinet Top

CG397.



Product Information

Description


This top attaches to an overhead cabinet, storage cabinet, or shelf unit and is sloped to ease cleaning and prevent the stacking of materials on top of the cabinet. Attachment hardware is included.

Notes

Order 1 of the following support components separately:

- Overhead cabinet with door (CG310., CG311., CG320., or CG321.)
- Overhead cabinet with double doors (CG312. or CG322.)
- Shelf unit (SM303. or SM305.)
- Storage cabinet (CG506.)

Dimensions



Sloped Overhead Cabinet Top


continued

Casework Storage

Specification Information				
Step 1.				
CG397.				
Step 2. Depth				
09	10½" deep			
16	17" deep			
18	18" deep			
24	24" deep			
Step 3. Width				
For 10½" deep (09) or 17" deep (16)				
24	24" wide			
30	30" wide			
36	36" wide			
48	48" wide			
For 18" deep (18) or 24" deep (24)				
24	24" wide			
30	30" wide			
48	48" wide			
Prices for Steps 1-3.				
	24	30	36	48
CG397. 09	\$120	140	160	200
16	\$160	180	200	240
18	\$170	190	—	250
24	\$210	230	—	290
Step 4. Surface Finish				
HF	inner tone light			+\$0
LT	light tone			+\$0

Light Bracket

SM910.



Product Information

Description

These brackets fasten a task light to the bottom of an overhead cabinet, a flat or slanted shelf in a 3- or 5-shelf unit, or a center shelf in a diagonal corner shelf unit. Package contains 1 pair.

The task light attaches to a 9"-deep shelf in the same position as the shelf, either flat or slanted. The task light attaches to a 16"-deep slanted shelf in a flat position. The 16"-deep flat shelf does not require brackets for attaching the task light.

Notes

Order task light separately:

- For 9"-deep shelf, order utility task light (G6132.)
- For 16"-deep shelf, order energy-efficient task light (G6120.)

Dimensions

Specification Information

Step 1.

SM910.

Step 2. Usage

- 09** for 9" deep flat or slanted shelf
16 for 16" deep slanted shelf



Prices for Steps 1-2.

SM910.09	\$37
16	\$37


Step 3. Surface Finish

For 16" deep slanted shelf (16)

HF	inner tone light	+\$0
LT	light tone	+\$0




CG399.



Description

Notes

Dimensions



Step 2. Type

- Step 3. Height

- Step 4. Depth


- Prices for Steps 1-4.

Step 5. Surface Finish

- Herman Miller for Healthcare*

Work Surface

CG210.



Product Information

Description


This laminate work surface with a backsplash attaches to the top of 1 or more base cabinets or hangs from an architectural wall and abuts the side of a base cabinet.

Notes

Order base cabinet (CG501., CG503., CG601., CG603., CG610., CG611., CG612., CG613., CG620., CG621., or CG622.) or sink cabinet (CG630., CG631., or CG632.) separately.

For work surface attached to wall, order work surface support (CG509.) separately.

Dimensions



Specification Information

Step 1.

CG210.

Step 2. Depth

18	18" deep
24	24" deep
30	30" deep

Step 3. Width

24	24" wide
30	30" wide
36	36" wide
42	42" wide
48	48" wide
60	60" wide
72	72" wide
96	96" wide

Prices for Steps 1-3.

	24	30	36	42	48	60	72
CG210. 18	\$152	167	181	195	210	238	266
24	\$190	209	226	244	262	297	333
30	\$209	230	253	275	297	342	387
							96
CG210. 18							\$324
24							\$405
30							\$476

Step 4. Edge Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0


Step 5. Surface Finish

See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1 +\$0
 Price Category 2 +\$30
 Price Category 3 +\$60


Corner Work Surface

CG217.




Product Information
Description
This laminate corner work surface with a backsplash attaches to the top of 1 or more base cabinets or hangs from an architectural wall and abuts the side of a base cabinet.
Notes
Order base cabinet (CG501., CG503., CG601., CG603., CG610., CG611., CG612., CG613., CG620., CG621., or CG622.) or sink cabinet (CG630., CG631., or CG632.) separately.
For work surface attached to wall, order work surface support (CG509.) separately.
Dimensions

Specification Information
Step 1.
CG217.
Step 2. Depth/Width
1848 18" deep x 48" wide
2448 24" deep x 48" wide
3060 30" deep x 60" wide
Prices for Steps 1-2.
CG217. 1848 \$406
2448 \$507
3060 \$663
Step 3. Edge Finish
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.
Price Category 1 +\$0
Step 4. Surface Finish
See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.
Price Category 1 +\$0
Price Category 2 +\$45
Price Category 3 +\$90



Sidesplash

CG295.



Product Information


Description

This component completes a work surface backsplash at the end of a run of base cabinets. Attachment hardware is included.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.

Dimensions



Specification Information

Step 1.

CG295.

Step 2. Depth

18	18" deep
24	24" deep
30	30" deep

Prices for Steps 1-2.

CG295. 18	\$20
24	\$26
30	\$35

Step 3. Edge Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
------------------	------

Step 4. Surface Finish


See laminates list for laminate numbers. First 2 digits of number indicate laminate line; remaining digits indicate laminate color.

Price Category 1	+\$0
Price Category 2	+\$5
Price Category 3	+\$10




Backsplash Filler

C2920.



Product Information
Description
This component attaches to a Casework work surface backsplash. It fills the gap created when a work surface is placed against terminal panels or wall strips. Package contains 6. Attachment hardware is included.
Notes
Order work surface (CG210.) or corner work surface (CG217.) separately.
Dimensions


Specification Information		
Step 1.		
C2920.96		\$120
Step 2. Surface Finish		
HF	inner tone light	+\$0
LT	light tone	+\$0



Casework Work Surfaces

Work Surface Support

CG509.



Product Information


Description

This support bracket holds a work surface or corner work surface that hangs from an architectural wall and abuts the side of a base cabinet.

Notes

Order work surface (CG210.) or corner work surface (CG217.) separately.

Dimensions



Specification Information

Step 1.

CG509.

Step 2. Depth

18 18" deep

24 24" deep

Step 3. Width

24 24" wide

30 30" wide

36 36" wide

48 48" wide

Prices for Steps 1-3.

	24	30	36	48
CG509. 18	\$86	91	95	105
24	\$91	95	100	109

Step 4. Surface Finish

HF	inner tone light			+\$0
LT	light tone			+\$0

Packaging Information

Keyed-Alike Information

Packaging

Packaging varies according to individual product. A surcharge may be added for special packaging requests. Contact Customer Care for specific information.

Keyed-Alike Locks

All components with locks may be specified keyed alike.

Keyed-alike Co/Struc® components with “CO”, “CV”, or “SM” prefix will be shipped with lock plugs installed.

Keyed-alike Casework components with “CG” prefix will be shipped with lock chassis installed. Lock plugs will be shipped separately, tagged for product and location as specified. Installation of lock plugs will be on-site according to instructions provided. Herman Miller, Inc., will not be responsible for lock plug installation.

To order keyed-alike locks for Casework components, specify keyed alike. Then place a separate order for plugs and locks:

1. List plug part number and key number (232092-XXX); select a key number (or numbers) between 226 and 427, excluding numbers 408 and 412.
2. List quantity of locks per key number.
3. Indicate this is a no charge order.

Keyed-alike locks are field installed.

A Medeco master key can be used to open any lock on Co/Struc components with a “CO” or “CV” prefix. Customers interested in ordering master keys should contact Customer Care with a written request. To remove lock plugs from existing components, loosen the hex nut inside the product and remove the plug. These locks are available in key numbers 01 through 99.

A Hudson master key can be used to open any lock on Co/Struc components with an “SM” prefix or on Casework components with a “CG” prefix. To remove lock plugs from existing components, a removal key (#232228) or lock change tool (#234559) is required. Customers interested in ordering master keys, removal keys, or lock change tools should contact Customer Care with a written request. Locks for “CG” prefix components are available in key numbers 226 through 427, excluding 408 and 412. Locks for “SM” prefix components are available in key numbers 226 through 324 but are ordered as key numbers 01 through 99. Use the following chart to cross-reference the number selection.

SM lock order numbers	Actual SM lock plug numbers
1	226
2	227
3	228
4	229
5	230
6	231
7	232
8	233
9	234
10	235
11	236
12	237
13	238
14	239
15	240
16	241
17	242
18	243
19	244
20	245
21	246
22	247
23	248
24	249
25	250
26	251
27	252
28	253
29	254
30	255
31	256
32	257
33	258
34	259
35	260
36	261
37	262
38	263
39	264
40	265
41	266
42	267
43	268
44	269
45	270
46	271
47	272
48	273
49	274
50	275

SM lock order numbers	Actual SM lock plug numbers
51	276
52	277
53	278
54	279
55	280
56	281
57	282
58	283
59	284
60	285
61	286
62	287
63	288
64	289
65	290
66	291
67	292
68	293
69	294
70	295
71	296
72	297
73	298
74	299
75	300
76	301
77	302
78	303
79	304
80	305
81	306
82	307
83	308
84	309
85	310
86	311
87	312
88	313
89	314
90	315
91	316
92	317
93	318
94	319
95	320
96	321
97	322
98	323
99	324



Laminate Open-Line Program

Customer's Own Laminate Information

Laminate Open-Line Program

Program Overview

The Laminate Open-Line program expands the offering of laminates approved for application to Herman Miller Casework products ("CG" prefix). This program includes selected laminates manufactured by three of the industry's leading suppliers: Formica, Pionite, and Wilsonart.

Laminates included in this program have been selected to complement Herman Miller's proprietary laminate offering, providing customers with additional color, pattern, and finish choices. Herman Miller's selected group of colors and finishes within each laminate line is available for order within the 20-day or less lead-time program.

Laminate Open-Line Warranty

All laminates included in this open-line program are supported by the Herman Miller for Healthcare Casework warranty.

Laminate Open-Line Order Information

All laminates included in the open-line program are numbered and priced for direct specification through normal Herman Miller, Inc., order entry channels.

To order Herman Miller Casework products with any of the open-line laminates included in this price book, enter the 4-digit laminate number on your product order.

Approved Open-Line Laminates

Formica®

General-Purpose Laminate - Grade #10

- #58 Matte finish
- #90 Polished finish

Lab-Grade Laminate - Grade #12

- #58 Matte finish
- #90 Polished finish

Pionite®

Standard Laminate

- S Suede
- J Crystal
- H Hi-Brite

Chemguard Laminate

- S Suede

Wilsonart®

General-Purpose Laminate

- #60 Matte finish
- #90 Crystal finish

Chemsurf® Laminate

- #60 Matte finish
- #90 Crystal finish

Customer's Own Laminate

If you require a laminate not approved by Herman Miller, please contact your Herman Miller for Healthcare dealer or Customer Care representative. Customer's Own Laminate will be handled on an individual basis to determine requirements for performance, testing, pricing, and lead times.



Casework Laminates

Price Category 1

Herman Miller Laminates

Order Number	Description
9VHF	Inner Tone Light
9VLT	Light Tone

Price Category 2

Wilsonart® Laminates

Order Number	Reference Number	Description
9W01	10734-60-107	Limber Maple
9W03	10767-60-107	Maple Sugar
9W04	10768-60-107	Golden Anagre
9W05	1500N-60-107	Grey
9W06	1503-60-107	Putty
9W07	1511-60-107	Mandarin Red
9W08	1530-60-107	Beige
9W09	1531-60-107	Light Beige
9W10	1570-60-107	White
9W11	1572-60-107	Antique White
9W12	1573-60-107	Frosty White
9W13	1595-60-107	Black
9W14	1723-1-107	Bleu Windrift
9W15	1730A-1-107	Light Tan Marble
9W16	1741-1-107	Verdi Pompeii
9W17	1742-1-107	White Pompeii
9W18	1744-1-107	Tearose Illusion
9W19	1745-60-107	Classic Black
9W21	1753-1-107	Canyon Blu
9W22	1754-1-107	Canyon
9W23	1755-1-107	Canyon Black
9W24	1762-60-107	Aegean
9W25	1763-60-107	Brune Slate
9W26	1786-60-107	Verdi Oxide
9W27	1787-60-107	Oxide
9W28	1788-60-107	Sierra Marble
9W29	2902-60-107	Chargrey Suede
9W30	2904-60-107	White Suede
9W31	2932-60-107	Almond Leather
9W32	2970-60-107	White Leather
9W33	4105-60-107	Blue Hirado
9W34	4141-60-107	Mocha Glace
9W35	4142-60-107	Grey Glace

9W36	4143-60-107	Neutral Glace
9W37	4144-60-107	Green Glace
9W38	4166-60-107	Pampas
9W39	4168-60-107	Grey Pampas
9W40	4170-60-107	Beige Pampas
9W41	4202-60-107	Favorite Mix
9W42	4204-60-107	Symmetry
9W43	4205-60-107	Confetti
9W44	4206-60-107	Verdigris
9W45	4207-60-107	Fresco
9W46	4208-60-107	Blue Cashmere
9W47	4209-60-107	Beige Cashmere
9W48	4226-60-107	White Hieroglyph
9W49	4305-90-107	Erin Glen
9W50	4306-90-107	Erin
9W51	4406-50-107	Diamondhead
9W53	4509-60-107	Stratos
9W54	4550-1-107	Granite
9W56	4552-50-107	Ebony Star
9W57	4556-60-107	Scopia
9W58	4557-60-107	Dakota Ridge
9W59	4558-60-107	Artesia
9W60	4590-60-107	Grey Millstone
9W61	4604-50-107	Pearl Sand
9W62	4607-60-107	Caldera Grey
9W63	4608-60-107	Caldera Beige
9W64	4609-60-107	Caldera Rose
9W65	4612-90-107	Spektrum
9W67	4621-60-107	White Nebula
9W68	4622-60-107	Grey Nebula
9W69	4623-60-107	Graphite Nebula
9W70	4624-60-107	Beige Nebula
9W71	4625-60-107	Mauve Nebula
9W72	4626-60-107	Blue Nebula
9W73	4627-60-107	Hunter Nebula
9W74	4628-60-107	Sand Nebula
9W75	4629-60-107	Blush Nebula
9W76	4630-60-107	Cloud Nebula
9W77	4632-60-107	Teal Nebula
9W78	4633-60-107	Natural Nebula
9W79	4634-60-107	Storm Nebula
9W80	4640-60-107	Dove Moraine
9W81	4644-60-107	Alberta Tan
9W82	4645-60-107	Caldera Jade
9W83	4646-60-107	Blue Moraine
9W84	4647-60-107	Navajo Dusk

Casework Laminates

continued

9W85	4649-60-107	Adriatic Sands
9W86	4651-60-107	Navy Legacy
9W87	4654-60-107	Olive Legacy
9W88	4655-60-107	Natural Legacy
9W89	4656-60-107	Bronze Legacy
9W90	4657-60-107	Green Legacy
9W91	4661-60-107	Iron Legacy
9W92	4662-60-107	Heather Legacy
9W93	4663-60-107	Tawny Legacy
9W94	4667-60-107	Green Tigris
9W95	4669-60-107	Natural Tigris
9W96	4672-60-107	Curry Tigris
9W97	4673-60-107	Saffron Tigris
9W98	4674-60-107	Evening Tigris
9W99	4676-60-107	Lichen
9WA1	4678-60-107	Galvan
9WA5	4738-60-107	Ochre Roletta
9WA6	4739-60-107	Clay Roletta
9WA7	4740-60-107	Roletta
9WA8	7039-60-107	Windsor Mahogany
9WA9	7040A-60-107	Figured Mahogany
9WB1	7054-60-107	Wild Cherry
9WB2	7061-60-107	Natural Pear
9WB3	7062-60-107	Congo Spruce
9WB4	7110T-60-107	Montana Walnut
9WB5	7209-60-107	Nepal Teak
9WB6	7213-60-107	Ivorywood
9WB7	7806-60-107	Bannister Oak
9WB8	7816-60-107	Solar Oak
9WB9	7817-60-107	Oregon Oak
9WC1	7850-60-107	Beigewood
9WC2	7852-60-107	Woodland Rose
9WC4	7885T-60-107	English Oak
9WC5	7888-60-107	Golden Oak
9WC6	7907-1-107	Pearl White Birdseye
9WC7	7908-60-107	Sanibel Maple
9WC8	7909-60-107	Fusion Maple
9WC9	7911-60-107	Manitoba Maple
9WD2	D14-60-107	Port
9WD3	D15-60-107	Wineberry
9WD5	D25-60-107	Atlantis
9WD6	D28-60-107	Seaspray
9WD7	D29-60-107	Silverpine
9WD8	D301-60-107	Chiffon
9WD9	D304-60-107	Rose Buff
9WE1	D30-60-107	Natural Almond

9WE2	D307-60-107	Hollyberry
9WE3	D308-60-107	Ivory Rose
9WE6	D315-60-107	Platinum
9WE7	D318-60-107	Shadowblue
9WE8	D319-60-107	Blue Ice
9WE9	D321-60-107	Brittany Blue
9WF1	D322-60-107	Bluejay
9WF2	D327-60-107	Pepperdust
9WF3	D328-60-107	Larkspur
9WF4	D329-60-107	Bisque
9WF5	D331-60-107	Sand
9WF6	D337-60-107	Pebble
9WF7	D341-60-107	Marigold
9WG1	D351-60-107	Tyrol Green
9WG2	D354-60-107	Designer White
9WG4	D368-60-107	Taupetone
9WG5	D369-60-107	Burgundy
9WG9	D379-60-107	Indigo
9WH1	D381-60-107	Fashion Grey
9WH2	D391-60-107	Cloud
9WH7	D402-60-107	Oasis
9WH8	D403-60-107	White Sand
9WJ1	D412-60-107	Fortuneberry
9WJ2	D414-60-107	Shogun
9WJ4	D417-60-107	Lapis Blue
9WJ6	D421-60-107	Midori
9WJ7	D422-60-107	Ming Gold
9WJ8	D431-60-107	Alabaster
9WJ9	D432-60-107	Cashmere
9WK1	D433-60-107	Vicuna
9WK2	D435-60-107	Pongee
9WK3	D43-60-107	Doeskin
9WK4	D436-60-107	Fleece
9WK5	D437-60-107	Puma
9WK6	D439-60-107	Wallaby
9WK7	D442-60-107	Andes Clay
9WK8	D443-60-107	Tibetan Jade
9WL4	D448-60-107	Juarez Flower
9WL9	D455-60-107	Montpellier
9WM1	D456-60-107	Cafe' Creme
9WM2	D457-60-107	St. Croix Sky
9WN4	D47-60-107	Moccasin
9WN6	D477-60-107	Provence Pear
9WN7	D478-60-107	Amazon
9WN9	D481-60-107	Tuscan Olive
9WP1	D482-60-107	Black Forest

Casework Laminates

continued

9WP2	D48-60-107	Thistle
9WP3	D50-60-107	Khaki Brown
9WP4	D72-60-107	Featherstone
9WP5	D73-60-107	Pewter
9WP6	D75-60-107	Porcelain Rose
9WP7	D76-60-107	Mauve Blush
9WP8	D79-60-107	Hunter Green
9WP9	D90-60-107	North Sea
9WQ1	D91-60-107	Slate Grey
9WQ2	D92-60-107	Dove Grey
9WQ3	D96-60-107	Shadow
9WQ4	D97-60-107	Haze
9WQ5	D98-60-107	Wildrose
9WQ6	P361A-60-107	Valley Pecan
9WQ7	W313A-60-107	Gunstock Walnut

Formica® Laminates

Order Number	Reference Number	Description
9X10	453-10-58	Stone
9X12	459-10-58	Brite White
9X13	460-10-58	Ivory
9X14	462-10-58	Contract Gray
9X15	463-10-58	Sail White
9X16	464-10-58	Graystone
9X17	466-10-58	Porcelain
9X18	468-10-58	Light Mink
9X19	837-10-58	Graphite
9X20	839-10-58	Stop Red
9X22	845-10-58	Spectrum Red
9X24	851-10-58	Spectrum Blue
9X27	858-10-58	Pumice
9X30	866-10-58	Doeskin
9X31	875-10-58	Cordovan
9X37	898-10-58	Vanilla
9X38	899-10-58	Desert Beige
9X39	902-10-58	Platinum
9X41	909-10-58	Black
9X42	912-10-58	Storm
9X43	913-10-58	Indian Smoke
9X44	914-10-58	Marine Blue
9X45	918-10-58	Neutral White
9X46	920-10-58	Almond
9X47	921-10-58	Birch
9X48	923-10-58	Surf
9X50	925-10-58	Champagne

9X51	927-10-58	Folkstone
9X52	928-10-58	Mouse
9X53	929-10-58	Oyster Gray
9X55	931-10-58	Putty Gray
9X56	932-10-58	Antique White
9X57	933-10-58	Mission White
9X60	949-10-58	White
9X63	953-10-58	Ice White
9X65	956-10-58	Victorian Teal
9X67	958-10-58	Beige
9X68	961-10-58	Fog
9X72	966-10-58	Burgundy
9X74	969-10-58	Navy Blue
9X78	2483-10-58	Ultramarine
9X82	7025-10-58	Sunlight
9X87	7030-10-58	Silver Moss
9X89	7197-10-58	Dover White
9X97	126-10-90	Nile Dorian Marble
9X98	132-10-90	Classic Onyx
9X99	220-10-58	Beige Clear Sand
9XA3	224-10-58	Folkstone Clear Sand
9XA8	253-10-58	Multifleck
9XA9	254-10-58	Varifleck
9XB1	271-10-90	Blackstone
9XB2	658-10-90	Folkstone American Granite
9XB3	659-10-90	Nile American Granite
9XB7	663-10-58	Antique White Papercraft
9XB8	689-10-58	Stellar
9XB9	692-10-58	Folkstone Celesta
9XC1	715-10-58	Champagne Papyrus
9XC2	718-10-58	Moonstone Papyrus
9XC3	726-10-58	Almond Papyrus
9XC4	754-10-90	Black Faux Marble
9XC5	760-10-58	Ivory Brushing
9XC8	771-10-58	White Kid Leather
9XC9	827-10-90	White Onyx
9XD3	3177-10-58	Omega Granite
9XD4	7022-10-58	Natural Canvas
9XD5	7023-10-58	Gray Canvas
9XD7	7213-10-58	Earth Wash
9XD9	7215-10-58	Earthen Terra
9XE1	7216-10-58	Sienna Terra
9XE2	7217-10-58	Arctic Terra
9XE3	7218-10-58	Prairie Terra
9XE4	7219-10-58	Forest Terra
9XE5	7220-10-58	Tundra Terra

Casework Laminates

continued

9XE8	7240-10-58	Parchment
9XF1	7244-10-58	Sun Wash
9XF2	7245-10-58	Meadow Wash
9XF3	7246-10-58	Mountain Wash
9XF8	299-10-58	Ebony Oxide
9XF9	300-10-58	Sail White Oxide
9XG2	303-10-58	Antique White Oxide
9XG4	305-10-58	Bayou Oxide
9XH2	503-10-58	Stone Grafix
9XH4	505-10-58	Dusty Jade Grafix
9XH5	506-10-58	Beige Grafix
9XH6	507-10-58	Folkstone Grafix
9XH7	508-10-58	Mouse Grafix
9XH8	510-10-58	Bayou Grafix
9XJ1	513-10-58	Mink Grafix
9XJ2	515-10-58	Graphite Grafix
9XJ4	520-10-58	Shell Grafix
9XJ5	522-10-58	Blue Silk Grafix
9XJ6	680-10-58	Granito
9XK2	1782-10-58	Stardust
9XK3	1783-10-58	Firedust
9XK4	1816-10-58	Fogdust
9XK5	1817-10-58	Almond Dust
9XK8	2486-10-58	Bayou Dust
9XK9	2488-10-58	Mouse Dust
9XL1	2491-10-58	Berry Dust
9XL2	2493-10-58	Stone Dust
9XL3	2494-10-58	Dusty Jade Dust
9XL4	2496-10-58	Antique White Dust
9XL6	7013-10-58	Santa Fe Slate
9XL7	7014-10-58	Colorado Slate
9XL8	7015-10-58	Glacier Slate
9XM1	7017-10-58	Spring Foliage
9XM2	7018-10-58	Navy Grafix
9XM3	7019-10-58	Porcelain Grafix
9XM4	7032-10-58	Summer Foliage
9XM5	7109-10-58	Monterrey Slate
9XN3	7208-10-58	Mineral Metalcraft
9XN5	7226-10-58	Crayon
9XN9	7230-10-58	Desert Erosion
9XP1	7231-10-58	Fallen Leaves

Pionite® Laminates

Order Number	Reference Number	Description
9Yo3	SW826-S-G48	Angel White
9Yo4	SW836-S-G48	Winter White
9Yo5	SW811-S-G48	White
9Yo6	SW806-S-G48	Carnation White
9Yo7	SW813-S-G48	Ice White
9Yo8	SG254-S-G48	Ash Gray
9Yo9	SG241-S-G48	Folkstone Gray
9Y10	SG204-S-G48	Mercury
9Y11	SG246-S-G48	Dolphin Gray
9Y12	SG213-S-G48	Opti Gray
9Y13	SG214-S-G48	Bankers Gray
9Y14	SG211-S-G48	Ingot Gray
9Y15	SC228-S-G48	Slate
9Y16	SG208-S-G48	Limestone
9Y17	SG203-S-G48	Pumice
9Y18	SG251-S-G48	Feather Gray
9Y19	SG224-S-G48	Stonedust
9Y20	SW841-S-G48	Warm White
9Y21	SW803-S-G48	Eggshell White
9Y22	SW812-S-G48	Tawny White
9Y23	ST655-S-G48	Almond
9Y24	ST613-S-G48	Tapioca
9Y25	ST683-S-G48	Pebble
9Y26	ST606-S-G48	Taupe
9Y27	SG210-S-G48	Greige
9Y28	SG209-S-G48	Putty Gray
9Y29	ST653-S-G48	Bisque
9Y30	ST650-S-G48	Palomino
9Y31	ST617-S-G48	Beige
9Y32	ST610-S-G48	Squash
9Y33	ST622-S-G48	Caramel
9Y36	ST604-S-G48	Nubian Brown
9Y37	ST643-S-G48	Cameo
9Y40	SR509-S-G48	Rose Ash
9Y41	SR505-S-G48	Navajo Red
9Y42	SR521-S-G48	Burgundy
9Y43	SR516-S-G48	Cerise
9Y44	SR522-S-G48	Raspberry
9Y45	SP402-S-G48	Plum
9Y46	SP401-S-G48	Royal Burgundy
9Y47	SB005-S-G48	French Blue
9Y49	SP409-S-G48	Iris
9Y50	SB007-S-G48	Navy Blue

Casework Laminates

continued

9Y51	SB006-S-G48	Indigo Blue
9Y52	SB016-S-G48	Deep Teal
9Y53	SV720-S-G48	Sage
9Y54	SV718-S-G48	Pistachio
9Y55	SG240-S-G48	Moss Gray
9Y56	SV710-S-G48	Anthra Green
9Y57	SV712-S-G48	Blue Spruce
9Y58	SV711-S-G48	Neptune
9Y59	SB010-S-G48	Teal
9Y60	SV704-S-G48	Boxwood
9Y62	SY914-S-G48	Primary Yellow
9Y63	SY913-S-G48	Forsythia
9Y64	SO312-S-G48	Bittersweet
9Y66	SR520-S-G48	Primary Red
9Y67	SV709-S-G48	Springtime
9Y68	SB009-S-G48	Royal Blue
9Y70	SE101-S-G48	Black
9Y73	MW550-S-G48	White Fiesta
9Y74	AG321-S-G48	Gray Fiesta
9Y75	MG041-S-G48	Ash Fiesta
9Y76	MB060-S-G48	Vermont Granite
9Y77	MR130-S-G48	Rose Granite
9Y78	MR140-S-G48	Ruby Granite
9Y79	MO001-S-G48	Breccia Granite
9Y80	MG021-G-G48	Smoke Granite
9Y81	ME001-J-G48	Black Granite
9Y82	MB040-G-G48	Atlantis Pavia
9Y83	MV430-S-G48	Jade Pavia
9Y84	MB050-G-G48	Blue Pavia
9Y85	MT270-G-G48	Beige Pavia
9Y86	MW530-G-G48	Alabaster Marmor
9Y87	MG011-G-G48	Greystone
9Y88	ME011-G-G48	Pyrenees Marble
9Y89	MW540-S-G48	Bianco Marble
9Y90	AG281-S-G48	Storm Cirrus
9Y91	AT971-S-G48	Antique Cirrus
9Y92	AW841-S-G48	Ivory Cirrus
9Y93	AB251-S-G48	Venetian Fresco
9Y94	AV771-S-G48	Agean Fresco
9Y95	AB141-S-G48	Baltic Fresco
9Y96	AT121-S-G48	Moroccan Fresco
9Y97	AR301-S-G48	Melange
9Y98	AV731-S-G48	Celedon Impression
9Y99	AB121-S-G48	Opal Impression
9YA1	AP321-S-G48	Amethyst Impression
9YA2	AV741-S-G48	Slate Impression

9YA3	AW812-S-G48	White Tundra
9YA4	AT892-S-G48	Etched Gold
9YA5	AV645-S-G48	Etched Verdigris
9YA6	AG232-S-G48	Etched Platinum
9YA7	AT541-S-G48	Almond Woodsong
9YA8	AT521-S-G48	Winter Americana
9YA9	LG110-S-G48	Ivory Kid
9YB1	AT383-S-G48	Parchment
9YB2	MT300-S-G48	Natural Limestone
9YB3	AV781-S-G48	Sage Coral
9YB4	AW851-S-G48	White Talc
9YB5	AR331-S-G48	Blush Talc
9YB6	AG381-S-G48	Mineral Talc
9YB7	AG361-S-G48	Graphite Talc
9YB8	AV674-S-G48	Olive Organix
9YB9	AT931-S-G48	Terre Cotta Organix
9YC1	AB051-S-G48	Indigo Organix
9YC2	AW785-S-G48	White Spektrum
9YC3	AG251-S-G48	Gray Spektrum
9YC4	AE021-S-G48	Graphite Spektrum
9YC5	AV664-S-G48	Green Spektrum
9YC6	AT921-S-G48	Sand Spektrum
9YC7	AR281-S-G48	Blush Spektrum
9YC8	AB041-S-G48	Cloud Spektrum
9YC9	AW831-S-G48	White Chromatix
9YD1	AT641-S-G48	Sand Chromatix
9YD3	AT631-S-G48	Flaxseed Chromatix
9YD4	AG111-S-G48	Gray Chromatix
9YD5	AE011-J-G48	Charcoal Chromatix
9YD6	AT881-S-G48	Almond Jaguar
9YD7	AG241-S-G48	Opti Gray Jaguar
9YD8	AV654-S-G48	Verde Jaguar
9YD9	AR271-S-G48	Rose Jaguar
9YE2	AW775-S-G48	White Jaguar
9YE3	AW881-S-G48	Eggshell Mosaic
9YE4	AT271-S-G48	Pebble Mosaic
9YE9	AR291-S-G48	Rose Ethos
9YF2	AB071-S-G48	Horizon Ethos
9YF3	AB101-S-G48	Navy Frost
9YF4	AG271-S-G48	Gray Frost
9YF5	AV701-S-G48	Spruce Frost
9YF6	AB091-S-G48	Horizon Blue Frost
9YF7	AV711-S-G48	Emerald Santos
9YF8	AG291-S-G48	Gray Santos
9YG1	AT951-S-G48	Neutral Santos
9YG2	AG301-S-G48	Pebble Santos

Casework Laminates

continued

9YG3	AW871-S-G48	Vanilla Fiber
9YG4	AT991-S-G48	Wheat Fiber
9YG5	AT101-S-G48	Oatmeal Fiber
9YG6	AT291-S-G48	Chamomile Fiber
9YG7	AY102-S-G48	Curry Fiber
9YG8	AT981-S-G48	Coriander Fiber
9YG9	AT241-S-G48	Cinnamon Fiber
9YH1	AO101-S-G48	Chili Fiber
9YH2	AV861-S-G48	Kale Fiber
9YH3	AV721-S-G48	Thyme Fiber
9YH4	AV751-S-G48	Celery Fiber
9YH5	AG371-S-G48	Cardamom Fiber
9YH6	AT251-S-G48	Beige Crepe
9YH7	AT111-S-G48	Caramel Crepe
9YH8	AP311-S-G48	Burgundy Crepe
9YH9	AG351-S-G48	Slate Crepe
9YJ1	AB111-S-G48	Navy Crepe
9YJ2	AG341-S-G48	Opti Gray Crepe
9YJ3	AG331-S-G48	Stonedust Crepe
9YJ4	WM872-S-G48	White Decor Maple
9YJ5	WO991-S-G48	Scrubbed Oak
9YJ6	WO531-S-G48	Pêche Contempo
9YJ7	WO541-S-G48	Almond Contempo
9YJ8	WM861-S-G48	Seaside Maple
9YJ9	WE261-S-G48	White Elm
9YK1	WM791-S-G48	Hardrock Maple
9YK2	WM971-S-G48	Sugar Maple
9YK3	WA001-S-G48	Clear Ash
9YK4	WP101-S-G48	Knotty Pine
9YK5	WG101-S-G48	Acadia Beech
9YK6	WB201-S-G48	Lg Butcher Block
9YK7	WB071-S-G48	Select BB
9YK8	WO891-S-G48	Sherwood Oak
9YK9	WO361-S-G48	Amberwood
9YL1	WO951-S-G48	Fine Oak
9YL2	WA011-S-G48	Summit Ash
9YL3	WX041-S-G48	Anigre
9YL4	WM951-S-G48	Honey Maple
9YL6	WC431-H-G48	Millwork Cherry
9YL7	WX031-S-G48	Pearwood
9YL8	WC421-S-G48	Oiled Cherry
9YL9	WT741-S-G48	Taj Teak
9YM1	WT601-S-G48	Scandia Teak
9YM2	WO581-S-G48	Wisconsin Oak
9YM3	WO862-S-G48	Hearth Oak

9YM4	WW971-S-G48	Gunstock Savoy Walnut
9YM5	WC331-H-G48	Victorian Cherry
9YM6	WC411-S-G48	Yorkshire Cherry
9YM7	WY031-S-G48	Formal Mahogany

Price Category 3

Chemically Resistant Laminates

Order Number	Reference Number	Description
9ZKR	Herman Miller KR	Inner Tone Light
9ZKH	Herman Miller KH	Light Tone
9Z01	Pionite SW826-S-H48	Angel White
9Z02	Pionite SW836-S-H48	Winter White
9Z03	Pionite SW811-S-H48	White
9Z04	Pionite SW806-S-H48	Carnation White
9Z05	Pionite SW813-S-H48	Ice White
9Z06	Pionite SG254-S-H48	Ash Gray
9Z07	Pionite SG241-S-H48	Folkstone Gray
9Z08	Pionite SG204-S-H48	Mercury
9Z09	Pionite SG246-S-H48	Dolphin Gray
9Z10	Pionite SG213-S-H48	Opti Gray
9Z11	Pionite SG214-S-H48	Bankers Gray
9Z12	Pionite SG211-S-H48	Ingot Gray
9Z13	Pionite SC228-S-H48	Slate
9Z14	Pionite SG208-S-H48	Limestone
9Z15	Pionite SG203-S-H48	Pumice
9Z16	Pionite SG251-S-H48	Feather Gray
9Z17	Pionite SG224-S-H48	Stonedust
9Z18	Pionite SW841-S-H48	Warm White
9Z19	Pionite SW803-S-H48	Eggshell White
9Z20	Pionite SW812-S-H48	Tawny White
9Z21	Pionite ST655-S-H48	Almond
9Z22	Pionite ST613-S-H48	Tapioca
9Z23	Pionite ST683-S-H48	Pebble
9Z24	Pionite ST606-S-H48	Taupe
9Z25	Pionite SG210-S-H48	Greige
9Z26	Pionite SG209-S-H48	Putty Gray
9Z27	Pionite ST653-S-H48	Bisque
9Z28	Pionite ST650-S-H48	Palomino
9Z29	Pionite ST617-S-H48	Beige
9Z30	Pionite ST610-S-H48	Squash
9Z31	Pionite ST622-S-H48	Caramel
9Z34	Pionite ST604-S-H48	Nubian Brown
9Z35	Pionite ST643-S-H48	Cameo
9Z38	Pionite SR509-S-H48	Rose Ash

Casework Laminates

continued


9Z39	Pionite SR505-S-H48	Navajo Red
9Z40	Pionite SR521-S-H48	Burgundy
9Z41	Pionite SR516-S-H48	Cerise
9Z42	Pionite SR522-S-H48	Raspberry
9Z43	Pionite SP402-S-H48	Plum
9Z44	Pionite SP401-S-H48	Royal Burgundy
9Z45	Pionite SB005-S-H48	French Blue
9Z47	Pionite SP409-S-H48	Iris
9Z48	Pionite SB007-S-H48	Navy Blue
9Z49	Pionite SB006-S-H48	Indigo Blue
9Z50	Pionite SB016-S-H48	Deep Teal
9Z51	Pionite SV720-S-H48	Sage
9Z52	Pionite SV718-S-H48	Pistachio
9Z53	Pionite SG240-S-H48	Moss Gray
9Z54	Pionite SV710-S-H48	Anthra Green
9Z55	Pionite SV712-S-H48	Blue Spruce
9Z56	Pionite SV711-S-H48	Neptune
9Z57	Pionite SB010-S-H48	Teal
9Z58	Pionite SV704-S-H48	Boxwood
9Z60	Pionite SY914-S-H48	Primary Yellow
9Z61	Pionite SY913-S-H48	Forsythia
9Z62	Pionite SO312-S-H48	Bittersweet
9Z64	Pionite SR520-S-H48	Primary Red
9Z65	Pionite SV709-S-H48	Springtime
9Z66	Pionite SB009-S-H48	Royal Blue
9Z68	Pionite SE101-S-H48	Black
9ZAA	Wilsonart 1573-60-390	Frosty White
9ZAB	Wilsonart 1572-60-390	Antique White
9ZAC	Wilsonart D30-60-390	Natural Almond
9ZAD	Wilsonart 1531-60-390	Light Beige
9ZAF	Wilsonart D97-60-390	Haze
9ZAG	Wilsonart D318-60-390	Shadowblue
9ZAH	Wilsonart D381-60-390	Fashion Grey
9ZAJ	Wilsonart D92-60-390	Dove Grey
9ZAK	Wilsonart D96-60-390	Shadow
9ZAL	Wilsonart D91-60-390	Slate Grey
9ZAM	Wilsonart 4143-60-390	Neutral Glace
9ZAN	Wilsonart 4141-60-390	Mocha Glace
9ZAP	Wilsonart 4649-60-390	Adriatic Sands
9ZAQ	Wilsonart 4142-60-390	Grey Glace
9ZAR	Wilsonart 1595-60-390	Black
9ZAS	Wilsonart 4622-60-390	Grey Nebula
9ZAT	Wilsonart 4590-60-390	Grey Millstone
9ZAU	Wilsonart 4557-60-390	Dakota Ridge
9ZAV	Wilsonart 4617-60-390	Faience
9ZAW	Wilsonart 4623-60-390	Graphite Nebula
9ZFA	Formica 840-12-58	Black




20-Day or Less and Assigned Lead-Time Order Information

20-Day or Less and Assigned Lead-Time Programs

20-day or less and assigned lead-time products are designed for established Herman Miller dealers with pre-approved credit.

Items in this price book are distributed from multiple manufacturing facilities across the United States. All products and options in this price book not designated by an Assigned Lead-Time icon  will ship in 20 business days or less after being acknowledged by Herman Miller.

Products and options designated by an  are on the assigned lead-time program and will ship in 20 days, less than 20 days, or more than 20 days from order acknowledgement.

Ordering Procedure

Please place orders through Order Manager. If this is not available to you, place orders through mail or fax. Verbal purchase orders will not be accepted.

Order Entry fax number: (616) 654 3085.

For more information, contact your Customer Care representative at (877) 464 4681.

Shipments and Delivery

Per Herman Miller Terms and Conditions.

Changes and Cancellation

Per Herman Miller Terms and Conditions.



10-Day or Less Order Information

10-Day or Less Program

10-day or less products are designed for established Herman Miller dealers with pre-approved credit.

Items in the 10-day or less product list are distributed from Holland, Michigan. 10-day or less orders will be shipped within 10 business days or less from the date of Herman Miller acknowledgement.

Ordering Procedures

1. Please place orders through Order Manager. If this is not available to you, place orders through mail or fax. Verbal purchase orders will not be accepted. Order Entry fax number: (616) 654 3085.
2. Orders must be clearly marked "10-day or less program" and include only 10-day or less products. Orders not marked this way will be processed under standard lead times. Orders for standard lead times must be placed on a separate purchase order.
3. Only finishes listed in the 10-day or less product list are available through the 10-day or less program.
4. For any 10-day or less project over \$250,000 list price, please contact your Customer Care representative within 48 hours of order entry to reserve lead time.
5. For more information, contact your Customer Care representative at (877) 464 4681.

Shipments and Delivery

Per Herman Miller Terms and Conditions.

Changes and Cancellation

Per Herman Miller Terms and Conditions.

10-Day or Less Product List

Products listed below are available through the 10-day or less program.

CO201	C Lid	
CO203	C Tray/Shelf	
CO205	C Frame	
CO207FF	Drawer, No Lock	Finishes DA, HF, LT, LU, and TC only
CO208FF		
CO209FF		
CO212	L Cart	
CO219	Locker Tambour Door	Finishes HF and LT only
CO228	Drawer Bearer Package	
CO231	Subcontainer (package of 12)	
CO232		
CO233		
CO234		
CO235		
CO236		
CO238		
CO240		
CO242	Subdivider Vane	
CO252	C Wire Shelf	
CO266	Gripper Package	
CO277FF	Locker Lock Kit, Keyed Differently	
CO282	Standard Rail	
CO283	Rail Hardware (package of 300)	
CO284		
CO285		
CO342	TR3 Cart	
CO369	Subcontainer (package of 12)	
CO378	Subcontainer Divider (package of 12)	
CO379		
CO380		
CO396	Lockable Drawer Bearer	
CO479	Wall Rail, Seismic	
CO540	A-Size Drawer Organizer Kit	
CO541		
CO542		
CO543		
CO544		
CO545		
CO546	B-Size Drawer Organizer Kit	
CO547		
CO555	Counter Top	
CO556		
CO561FF	Locker, with Tambour Door, No Lock	Finishes HFHF and LTLT only
CO562FF	Locker, No Lock	

Index by Product Name

Casework	
ADA Sink Support, Wall Mounted	page(s) 255
Backsplash Filler	292
Base Cabinet	261
Base Cabinet Bottle Drawer	269
Base Cabinet Bottom	263
Base Cabinet Door	270
Base Cabinet Door, Interior Cabinet Adapter	272
Base Cabinet Double Doors	273
Base Cabinet Drawer	267
Base Cabinet Filler	274
Base Cabinet Pullout Shelf	266
Base Cabinet Shelf	265
Base Cabinet with Door	245
Base Cabinet with Double Doors	247
Base Cabinet with Drawer and Door	249
Base Cabinet with Drawer and Double Doors	251
Base Cabinet with Drawers	243
Corner Work Surface	290
Interior Cabinet Adapter	275
Light Bracket	287
Overhead Cabinet Back Panel	284
Overhead Cabinet Door	282
Overhead Cabinet Double Doors	283
Overhead Cabinet Filler	288
Overhead Cabinet with Door	278
Overhead Cabinet with Double Doors	280
Sidesplash	291
Sink Cabinet with Foot Pedal Cutout and Door	252
Sink Cabinet with Foot Pedal Cutout and Double Doors	254
Sloped Overhead Cabinet Top	285
Storage Cabinet	259
Storage Cabinet Door	276
Storage Cabinet Double Doors	277
Storage Cabinet with Shelves and Door	256
Storage Cabinet with Shelves and Double Doors	258
Work Surface	289
Work Surface Support	293
Co/Struc® System	
2-Way 90° Corner Connector	page(s) 40
2-Way Module Corner Connector	58
3-Way 90° Corner Connector	41
3-Way Module Corner Connector	59
4-Way Module Corner Connector	60
A-Size Drawer Organizer Kit	226
A-Size Drawer Subdivider	217

A-Style Shelf	151
Access Panel Assembly	81
ADA Sink Support for Storage, Systems Mounted	119
ADA Sink Support, Systems Mounted	118
Adapter Rail	14
Add-On Lab Module	52
Add-On Shelf	130
Add-On Shelf with Lip	124, 155
Anchor Kit	62
B-Size Drawer Organizer Kit	227
B-Style Rail-Hanging Shelf	150
B-Style Shelf	152
B-Style Shelf with Sloped Top	154
Back Panel	132
Backsplash Filler	109
Bottle Drawer Divider	147
Breakaway Lock Assembly	215
C Flipper	190
C Frame	186
C Frame Adapter Kit, Option Work Surface	107
C Frame Support	106
C Lid	189
C Tray/Shelf	191
C Wire Shelf	193
Cabinet	170
Cabinet Drip Pan	179
Cabinet Floor	180
Cabinet Interior Adapter	178
Cabinet Lock Kit	173
Cabinet Pullout Catheter Rack	175
Cash Drawer Insert	214
Cassette Assembly	218
Cassette Bin	221
Cassette Cover	220
Ceiling Entry Support Panel	35
Center Closure Panel	91
Chest Tube Rack	235
Clear Lid	211
Combination Access Panel Assembly	83
Counter Top	187
CS Lid	182
CS Tray/Shelf	184
CS Wire Shelf	185
CST Flipper	183
CST Frame	181
Diplomat Panel	55

Index by Product Name *continued*

Dispensing Rail	222
Distal Tube Bracket	177
Draw Rod	23
Drawer	198
Drawer Bearer Package	212
Drawer Caddy with Casters	195
Drawer with Lockable Lid	201
End Cap for Roller Rail	11
End Cap for Stand-Off Plate	13
End Cap for Standard or Seismic Rail	19
End Closure Assembly	88
End Closure for Wall-Attached Module	90
End Closure Panel	87
Endoscope Bracket	176
Equipment Rail, Systems Mounted	17
Equipment Rail, Wall Mounted	15
Filler Extrusion	76
Floor Anchor Adapter Kit	28
Floor Anchor Bracket, Action Office® Series 1	25
Floor Anchor Bracket, Action Office® Series 2	26
Floor Anchor Bracket, Co/Struc® Lab Module	24
Floor Anchor Setting Tool	63
Floor Anchor, Ethospace® Frame	27
Floor Entry Support Panel	33
Floor Track	77
Fogarty Catheter Rack	234
Frame Stiffener, Ethospace® Top Cap	30
Gripper Package	194
Heavy-Duty Corner Work Surface	110
Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2 90° Corner Connector	112
Heavy-Duty Option Work Surface	103
Heavy-Duty Storage Work Surface	101
Heavy-Duty Work Surface	99
Interface Connector	42
Interior Lab Bracket	64
Interior Shelf	67
Interior Shelf and Support Kit	69
IV Pole	238
IV Pole Clamp	239
L Cart	237
Label Clip	216
Light Bracket	133
Lockable Drawer Bearer	146, 213
Locked Drawer with Lid	204
Locked Drawer with Lockable Lid	207

Locker	160
Locker Lock Kit	168
Locker on Wheel Base	166
Locker on Wheel Base, with Tambour Door	163
Locker Tambour Door	167
Locker, with Tambour Door	157
Low Cup Sink Access Panel Assembly	85
Low Module	48
Low Module Terminal Panel	50
Mid-Run Support, Double Sided	45
Mid-Run Support, Single Sided	44
Module Support Leg	61
Module Wall Start	57
Open Support Panel	31
Overhead Cabinet Back Panel	136
Overhead Cabinet Door	134
Overhead Cabinet Double Doors	135
Overhead Storage Unit	123
Panel Connector	43
Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	29
Panel Wall Start	39
Pass-Through Module	49
Pass-Through Module Terminal Panel	53
Process Table	114
Pullout Catheter Rack	233
Rail-Hanging Work Surface	105
Rail for Roller Rail	9
Rail Hardware	8
Receptacle Cover, Wire Chase Rail	97
Roller for Roller Rail	10
Shelf Backsplash	70
Shelf End Cap	68
Shelf for Storage Unit	122
Shelf Label Clip	131, 156
Shelf Unit	125
Shelf Unit, 90° Corner	129
Shelf Unit, Diagonal Corner	127
Sink-Adaptable Module	116
Sink Unit for Storage	117
Stand-Off Plate for Roller Rail	12
Standard Rail	7
Storage Case	137
Storage Case Bottle Drawer	143
Storage Case Door	144
Storage Case Drawer	141
Storage Case Fixed Shelf	139

Index by Product Name *continued*

Storage Case Leg	149
Storage Case Pullout Shelf, Flat	140
Storage Case Rail, Corner	148
Storage Unit	120
Subcontainer	224
Subcontainer Divider	228
Subcontainer Handle	229
Subcontainer Handle Label	230
Subcontainer Hanger	232
Subcontainer Label Flag	231
Subdivider Vane	192
Support Panel	32
Support Panel Top Cap	46
Terminal Panel	37
Terminal Panel Bracket	47, 65
Toggle Wall Anchor	22
Top Shelf	66
Top Shelf for Wall-Attached Module	75
TR3 Cart	236
Vertical Chase Cover	56
Vinyl Base	78
Vinyl Base, 90° Outside Corner	79
Vinyl Base, End Cap	80
Wall-Attached Connector	73
Wall-Attached Module	71
Wall-Attached Wall Start Connector	74
Wall Attachment Channel Hardware	72
Wall Rail, Seismic	20
Wall Strip	21
Wire Chase Rail	92
Wire Chase Rail Bridge	95
Wire Chase Rail End Cap	96
Wire Chase Rail, 90° Inside Corner	94
Work Surface	98
Work Surface Backsplash	108

Index by Product Number

C1112.	Low Module	page(s) 48
C1114.	Pass-Through Module	49
C1119.	Add-On Lab Module	52
C1122.	Low Module Terminal Panel	50
C1122.	Pass-Through Module Terminal Panel	53
C1124.		
C1141.	Wall-Attached Module	71
C1146.	Diplomat Panel	55
C1150.	Open Support Panel	31
C1151.	Support Panel	32
C1155.	Floor Entry Support Panel	33
C1156.	Ceiling Entry Support Panel	35
C1161.	Terminal Panel	37
C1170.	Equipment Rail, Wall Mounted	15
C1172.	Equipment Rail, Systems Mounted	17
C1191.	Interior Lab Bracket	64
C1192.	Floor Track	77
C1193.	Vinyl Base	78
C1194.	Vinyl Base, 90° Outside Corner	79
C1195.	Vinyl Base, End Cap	80
C1211.	Panel Wall Start	39
C1212.	Module Wall Start	57
C1220.	2-Way Module Corner Connector	58
C1221.	2-Way 90° Corner Connector	40
C1230.	3-Way Module Corner Connector	59
C1231.	3-Way 90° Corner Connector	41
C1240.	4-Way Module Corner Connector	60
C1261.	Interface Connector	42
C1272.	Mid-Run Support, Single Sided	44
C1273.	Mid-Run Support, Double Sided	45
C1274.	Support Panel Top Cap	46
C1275.	Module Support Leg	61
C1276.	Anchor Kit	62
C1277.	Floor Anchor Setting Tool	63
C1281.	Panel Connector	43
C1292.	Wall Attachment Channel Hardware	72
C1293.	Wall-Attached Connector	73
C1295.	Wall-Attached Wall Start Connector	74
C1310.	Wire Chase Rail	92
C1311.		
C1312.		
C1320.	Wire Chase Rail, 90° Inside Corner	94
C1325.	Wire Chase Rail Bridge	95
C1330.	Wire Chase Rail End Cap	96

C1340.	Receptacle Cover, Wire Chase Rail	97
C1341.		
C1342.		
C1343.		
C1412.	Access Panel Assembly	81
C1422.	Combination Access Panel Assembly	83
C1432.	Low Cup Sink Access Panel Assembly	85
C1450.	End Closure Panel	87
C1451.	End Closure Assembly	88
C1452.		
C1453.	End Closure for Wall-Attached Module	90
C1455.	Center Closure Panel	91
C1460.	Vertical Chase Cover	56
C1510.	Shelf End Cap	68
C1510.	Top Shelf	66
C1511.	Interior Shelf	67
C1511.	Shelf End Cap	68
C1514.	Top Shelf for Wall-Attached Module	75
C1515.	Interior Shelf and Support Kit	69
C1520.	Shelf Backsplash	70
C1530.	Filler Extrusion	76
C1531.		
C1610.	Adapter Rail	14
C1611.	Storage Case Rail, Corner	148
C2112.	Heavy-Duty Storage Work Surface	101
C2113.	Heavy-Duty Option Work Surface	103
C2122.	Heavy-Duty Storage Work Surface	101
C2123.	Heavy-Duty Option Work Surface	103
C2192.	C Frame Support	106
C2193.	C Frame Adapter Kit, Option Work Surface	107
C2210.	Heavy-Duty Work Surface	99
C2215.	Heavy-Duty Corner Work Surface	110
C2216.	Heavy-Duty Corner Work Surface, Action Office® Series 1 or 2	
	90° Corner Connector	112
C2310.	Work Surface	98
C2410.	Rail-Hanging Work Surface	105
C2510.	Sink-Adaptable Module	116
C2511.	Sink Unit for Storage	117
C2520.	ADA Sink Support, Systems Mounted	118
C2521.	ADA Sink Support for Storage, Systems Mounted	119
C2910.	Work Surface Backsplash	108
C2920.	Backsplash Filler	109, 292
C3010.	A-Style Shelf	151
C3011.	B-Style Shelf	152
C3019.	Add-On Shelf with Lip	124, 155
C3091.	B-Style Shelf with Sloped Top	154

Index by Product Number *continued*

Index: Product Number

C3920. Terminal Panel Bracket	47, 65
C3921.	
C4115. Process Table	114
C4215.	
CG210. Work Surface	289
CG217. Corner Work Surface	290
CG295. Sidesplash	291
CG310. Overhead Cabinet with Door	278
CG311.	
CG312. Overhead Cabinet with Double Doors	280
CG320. Overhead Cabinet with Door	278
CG321.	
CG322. Overhead Cabinet with Double Doors	280
CG390. Overhead Cabinet Door	134, 282
CG391.	
CG392. Overhead Cabinet Double Doors	135, 283
CG393. Overhead Cabinet Back Panel	136, 284
CG397. Sloped Overhead Cabinet Top	285
CG398. Overhead Cabinet Filler	288
CG399.	
CG501. Base Cabinet	261
CG503.	
CG506. Storage Cabinet	259
CG508. ADA Sink Support, Wall Mounted	255
CG509. Work Surface Support	293
CG511. Base Cabinet Bottom	263
CG512.	
CG515. Base Cabinet Shelf	265
CG521. Base Cabinet Pullout Shelf	266
CG530. Base Cabinet Drawer	267
CG531.	
CG536. Base Cabinet Bottle Drawer	269
CG540. Base Cabinet Door	270
CG541.	
CG542. Base Cabinet Double Doors	273
CG545. Storage Cabinet Door	276
CG546.	
CG547. Storage Cabinet Double Doors	277
CG580. Base Cabinet Door, Interior Cabinet Adapter	272
CG581.	
CG590. Interior Cabinet Adapter	275
CG598. Base Cabinet Filler	274
CG599.	
CG601. Base Cabinet with Drawers	243
CG603.	

CG610. Base Cabinet with Door	245
CG611.	
CG612. Base Cabinet with Double Doors	247
CG613.	
CG620. Base Cabinet with Drawer and Door	249
CG621.	
CG622. Base Cabinet with Drawer and Double Doors	251
CG630. Sink Cabinet with Foot Pedal Cutout and Door	252
CG631.	
CG632. Sink Cabinet with Foot Pedal Cutout and Double Doors	254
CG710. Storage Cabinet with Shelves and Door	256
CG711.	
CG712. Storage Cabinet with Shelves and Double Doors	258
CO201 C Lid	189
CO202 CS Lid	182
CO203 C Tray/Shelf	191
CO204 CS Tray/Shelf	184
CO205 C Frame	186
CO206 CST Frame	181
CO207 Drawer	198
CO208	
CO209	
CO212 L Cart	237
CO213 C Flipper	190
CO214 CST Flipper	183
CO215. Draw Rod	23
CO217 Counter Top	187
CO219 Locker Tambour Door	167
CO228 Drawer Bearer Package	212
CO231 Subcontainer	224
CO232	
CO233	
CO234	
CO235	
CO236	
CO238	
CO240	
CO242 Subdivider Vane	192
CO252 C Wire Shelf	193
CO253 CS Wire Shelf	185
CO266 Gripper Package	194
CO274 Label Clip	216
CO277 Locker Lock Kit	168
CO282 Standard Rail	7

Index by Product Number *continued*

CO283	Rail Hardware	8
CO284		
CO285		
CO337	Cash Drawer Insert	214
CO338	End Cap for Standard or Seismic Rail	19
CO342	TR3 Cart	236
CO345	Rail for Roller Rail	9
CO346	Stand-Off Plate for Roller Rail	12
CO347	Roller for Roller Rail	10
CO348	End Cap for Roller Rail	11
CO365	Dispensing Rail	222
CO366		
CO367		
CO368		
CO369	Subcontainer	224
CO370	Drawer with Lockable Lid	201
CO371	Locked Drawer with Lockable Lid	207
CO372	Drawer with Lockable Lid	201
CO373	Locked Drawer with Lockable Lid	207
CO374	Drawer with Lockable Lid	201
CO375	Locked Drawer with Lockable Lid	207
CO376	Subcontainer Handle	229
CO377		
CO378	Subcontainer Divider	228
CO379		
CO380		
CO381	Subcontainer Label Flag	231
CO382		
CO389	Subcontainer Hanger	232
CO396	Lockable Drawer Bearer	213
CO398	Shelf Label Clip	131, 156
CO437	End Cap for Stand-Off Plate	13
CO458	A-Size Drawer Subdivider	217
CO459	Subcontainer Handle Label	230
CO462	Label Clip	216
CO465	IV Pole Clamp	239
CO466	Chest Tube Rack	235
CO467	Fogarty Catheter Rack	234
CO469	Pullout Catheter Rack	233
CO473	Breakaway Lock Assembly	215
CO478	IV Pole	238
CO479	Wall Rail, Seismic	20
CO481	Floor Anchor Bracket, Action Office® Series 1	25
CO482	Toggle Wall Anchor	22
CO483		
CO485	Floor Anchor Bracket, Action Office® Series 2	26

CO486	Floor Anchor, Ethospace® Frame	27
CO508	Clear Lid	211
CO521	B-Style Rail-Hanging Shelf	150
CO540	A-Size Drawer Organizer Kit	226
CO541		
CO542		
CO543		
CO544		
CO545		
CO546	B-Size Drawer Organizer Kit	227
CO547		
CO549	Locked Drawer with Lid	204
CO550		
CO551		
CO552	Floor Anchor Adapter Kit	28
CO553	Panel Stiffener, Action Office® Series 1 or Series 2 Top Cap	29
CO554	Frame Stiffener, Ethospace® Top Cap	30
CO555	Counter Top	187
CO556		
CO558	Wall Strip	21
CO559		
CO560		
CO561	Locker, with Tambour Door	157
CO562	Locker	160
CO563	Locker on Wheel Base	166
CO564	Locker on Wheel Base, with Tambour Door	163
CO565	Cassette Assembly	218
CO566	Cassette Cover	220
CO567	Cassette Bin	221
CO568	Floor Anchor Bracket, Co/Struc® Lab Module	24
CO692	Drawer Caddy with Casters	195
CO693		
CV300	Cabinet	170
CV310	Endoscope Bracket	176
CV311	Distal Tube Bracket	177
CV312	Cabinet Interior Adapter	178
CV313	Cabinet Drip Pan	179
CV314	Cabinet Floor	180
CV315	Cabinet Lock Kit	173
CV316	Cabinet Pullout Catheter Rack	175
SM100	Storage Unit	120
SM101		
SM102		
SM103		
SM104	Shelf for Storage Unit	122

Index: Product Number

Index by Product Number *continued*

Index: Product Number

SM105. Overhead Storage Unit	123
SM303. Shelf Unit	125
SM305.	
SM313. Shelf Unit, Diagonal Corner	127
SM315.	
SM329. Shelf Unit, 90° Corner	129
SM400. Storage Case	137
SM411. Storage Case Fixed Shelf	139
SM421. Storage Case Pullout Shelf, Flat	140
SM433. Storage Case Drawer	141
SM434.	
SM436. Storage Case Bottle Drawer	143
SM447. Storage Case Door	144
SM448.	
SM491. Lockable Drawer Bearer	146
SM495. Storage Case Leg	149
SM900. Add-On Shelf	130
SM910. Light Bracket	133, 287
SM930. Bottle Drawer Divider	147
SM940. Back Panel	132

For more information about our products and services or to see a list of dealers, please visit us at www.HermanMiller.com or call (800) 221 7158.

This book is recyclable. Printed in U.S.A. on recycled paper with soy-based inks. F.CO2006-D

© 2004 Herman Miller, Inc., Zeeland, Michigan

® HermanMiller, ●, Action Office, Aeron Chair Configuration, Aeron, Ambi, AO, Arrio, Aside, Caper, CLT, Co/Struc, DOT, Eames, Eames Lounge Chair Configuration, EnhancedAccess, Equa 2, Ergon 3, Ethospace, E-Wall, Formcoat, Ground Cloth, Herman Miller for Healthcare, Herman Miller Options, Kinemat, Kiva, Limerick, Meridian, Passage, Pedastool, Pellicle, Perspectives, PostureFit, Prospects, Reaction, Resolve, Scooter, V-Wall, Vary Easy, Wild Card, and Z-Axis are among the registered trademarks of Herman Miller, Inc., and its owned subsidiaries.

™ Avive, Boomerang, Bubbletack, Burdick Group, CoActive, Flex-Edge, FlexFront, Hopsak 2, Iota, Milafin, Mirra, Multiscrim, Multiscrim 2, Nelson, Q, Rapunzel, Recline Selector, Stackable, and TriFlex are among the trademarks of Herman Miller, Inc., and its owned subsidiaries.

SM TAP and Textile Alliance Program are among the service marks of Herman Miller, Inc.

Corian is a registered trademark of DuPont.

Formica brand laminate is a registered trademark of Formica Corporation.

GSA Advantage! is a registered trademark of General Services Administration.

Microsoft and Natural Keyboard are registered trademarks of Microsoft Corporation.

Pionite brand laminate is a registered trademark of Panolam Industries International, Inc.

Wilsonart brand laminate and Chemsurf are registered trademarks of Wilsonart International.